

# **GSENSE2020 channel merging description**

**REPORT**

XIMEA s.r.o.  
Lesna 52  
900 33 Marianka  
Slovakia

Email: [peter.schlosser@ximea.com](mailto:peter.schlosser@ximea.com)

Web: [www.ximea.com](http://www.ximea.com)

## HDR channel construction

The equation for construction of the HDR channel is as follows:

In case of GLOBAL SHUTTER (CHANNEL\_BITDEPTH=10):

$$\begin{aligned} val_{HG}(x, y) &= val_{HG\_SIG}(x, y) - val_{HG\_BLACK}(x, y) + BLACK\_OFFSET \\ val_{LG}(x, y) &= val_{LG\_SIG}(x, y) - val_{LG\_BLACK}(x, y) + BLACK\_OFFSET \end{aligned}$$

In case of ROLLING SHUTTER (CHANNEL\_BITDEPTH=12):

$$\begin{aligned} val_{HG}(x, y) &= val_{HG\_SIG}(x, y) - HG\_BLACK + BLACK\_OFFSET \\ val_{LG}(x, y) &= val_{LG\_SIG}(x, y) - LG\_BLACK + BLACK\_OFFSET \end{aligned}$$

For both cases:

$$\begin{aligned} \text{if } val_{HG}(x, y) < 0 &=> val_{HG}(x, y) = 0; \\ \text{if } val_{LG}(x, y) < 0 &=> val_{LG}(x, y) = 0; \end{aligned}$$

$$val_{LG\_LFSR}(x, y) = 2^{17-CHANNEL\_BITDEPTH} * val_{LG}(x, y) | (LFSR \& (2^{17-CHANNEL\_BITDEPTH} - 1))$$

$$\alpha = \frac{GAIN\_RATIO}{2^{17-CHANNEL\_BITDEPTH}}$$

$$\beta = (GAIN\_RATIO - 1) * BLACK\_OFFSET$$

$$val_{LG\_SCALED}(x, y) = \alpha * val_{LG\_LFSR}(x, y) - \beta$$

$$\gamma(x, y) = 1 - val_{HG}(x, y) / (HG\_RANGE + BLACK\_OFFSET)$$

$$\text{if } \gamma(x, y) < 0 &=> \gamma(x, y) = 0;$$

$$val_{hdr}(x, y) = \gamma(x, y) * val_{HG}(x, y) + (1 - \gamma(x, y)) * val_{LG\_SCALED}(x, y)$$

$$\delta = \frac{2^{16} - 1}{GAIN\_RATIO * (BLACK\_OFFSET + LG\_RANGE) - \beta}$$

$$val_{hdr\_norm}(x, y) = val_{hdr}(x, y) * \delta$$

$$\text{if } val_{hdr\_norm}(x, y) > 2^{16} - 1 &=> val_{hdr\_norm}(x, y) = 2^{16} - 1$$

### Calibrated parameters:

HG\_BLACK - black level of the high gain channel  
 LG\_BLACK - black level of the low gain channel  
 HG\_RANGE - valid range of the high gain channel (HG\_SATURATION\_VALUE - HG\_BLACK)  
 LG\_RANGE - valid range of the high gain channel (LG\_SATURATION\_VALUE - LG\_BLACK)  
 GAIN\_RATIO - ratio of the measured slopes of the response curves of high and low gain channels  
 BLACK\_OFFSET - black level offset of the merged image (arbitrary constant)