

XIMEA CamTool

Generated by Doxygen 1.8.10

Tue Sep 27 2016 15:36:28

Contents

1	Hierarchical Index	1
1.1	Class Hierarchy	1
2	Class Index	3
2.1	Class List	3
3	File Index	7
3.1	File List	7
4	Class Documentation	9
4.1	CxHsiCameraParams::CxBand Class Reference	9
4.2	CxBinaryLayer Class Reference	9
4.3	CxBinaryLayerContainer Class Reference	11
4.3.1	Member Function Documentation	12
4.3.1.1	deleteAllLayers()	12
4.4	CxCameraManager Class Reference	12
4.4.1	Detailed Description	14
4.4.2	Member Function Documentation	14
4.4.2.1	cameraClosed	14
4.4.2.2	cameraInfo(uint uiCameraIdx, SxCameraInfo &aInfo) const	14
4.4.3	Member Data Documentation	14
4.4.3.1	m_hLastCameraAcqError	14
4.4.3.2	m_hLastCameraLost	15
4.4.3.3	m_lock	15
4.4.3.4	m_mapCameras	15
4.4.3.5	m_setRunningCameras	15
4.5	CxCameraParam Class Reference	15
4.6	CxCameraParameters Class Reference	17
4.7	CxCameraSource Class Reference	20
4.7.1	Detailed Description	22
4.7.2	Member Function Documentation	22
4.7.2.1	currentOutputImageInfo(SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL)	22

4.7.2.2	run()	22
4.7.2.3	saveSettings(QDomDocument &aDoc, QDomElement &aElm) const	22
4.7.2.4	setupStart()	22
4.7.2.5	setupStop()	23
4.7.3	Member Data Documentation	23
4.7.3.1	m_aCurrentOutputFormat	23
4.7.3.2	m_bCurrentOutputFormatValid	23
4.7.3.3	m_bFakeTransportFormat	23
4.7.3.4	m_hCamera	23
4.7.3.5	m_iFrameNo	23
4.7.3.6	m_queLastFrameTimes	23
4.7.3.7	m_timerLastFrames	23
4.8	CxChain Class Reference	24
4.8.1	Detailed Description	26
4.8.2	Chains and chainable objects	26
4.8.2.1	Connecting non-chainble object to the chain	27
4.8.3	Member Function Documentation	27
4.8.3.1	branchAfter(CxChainable *pChainable, CxChainable *pAfter)	27
4.8.3.2	connectNonChainable(CxChainable *pChainable, QObject *pListener, const char *onDataAvailableSlot,	
4.8.3.3	disconnectNonChainable(CxChainable *pChainable, QObject *pListener, const char *onDataAvailableS	
4.8.3.4	duplicateBranch(CxChainable *pRoot, CxChainable *pLast, CxChainable *pNewRoot, bool bIncludeRo	
4.8.3.5	duplicateChainWithViews(CxChainable *pRoot, CxChainable *pNewRoot, bool bIncludeRoot)	28
4.8.3.6	insertAfter(CxChainable *pChainable, CxChainable *pAfter)	28
4.8.3.7	insertBefore(CxChainable *pChainable, CxChainable *pBefore)	29
4.8.3.8	loadSettings(const QDomElement &aDom, QString *pErrorMsg=NULL)	29
4.8.3.9	saveChainableSettings(const CxChainable *pChainable, QDomDocument &aXml, QDomElement &aDom	
4.8.3.10	start()	29
4.8.4	Member Data Documentation	29
4.8.4.1	m_eState	29
4.8.4.2	m_lock	29
4.8.4.3	m_pRoot	30
4.9	CxChainable Class Reference	30
4.9.1	Detailed Description	34
4.9.2	Implementation of chainable object	34
4.9.3	Member Enumeration Documentation	35
4.9.3.1	ExChainableFlag	35
4.9.3.2	ExDataSourceType	35
4.9.4	Member Function Documentation	36
4.9.4.1	chainDataSourceTypes() const	36
4.9.4.2	configurationWidget()	36

4.9.4.3	connect(QObject *pReceiver, const char *method)	36
4.9.4.4	data(QObject *pReceiver, qint32 iFrameNo)	37
4.9.4.5	dataAvailable	37
4.9.4.6	dataSources()	37
4.9.4.7	discardData(QObject *pReceiver, qint32 iFrameNo)	37
4.9.4.8	disconnect(QObject *pReceiver, const char *method)	37
4.9.4.9	enabled() const	38
4.9.4.10	errorMessage()	38
4.9.4.11	help(QString &slmagesPath)	38
4.9.4.12	init(QString *pErrMsg=NULL)=0	38
4.9.4.13	memoryPool()	38
4.9.4.14	objectBuffersInMemoryPool()	38
4.9.4.15	onFreePoolBuffers	39
4.9.4.16	processData(IxChainData *pReceivedData)	39
4.9.4.17	run	39
4.9.4.18	saveSettings(QDomDocument &aDoc, QDomElement &aElm) const	39
4.9.4.19	setChain(CxChain *pChain, bool bSetSuccessors=false)	40
4.9.4.20	setupStart()	40
4.9.4.21	setupStop()	40
4.9.4.22	start(CxChain::ExChainErrorState *pErrState=NULL, QString *pErrMsg=NULL)	40
4.9.4.23	startRun	40
4.9.4.24	stop()	41
4.9.4.25	waitForStarted()	41
4.9.4.26	waitForStopped()	41
4.9.5	Member Data Documentation	41
4.9.5.1	m_bDataConsumer	41
4.9.5.2	m_bEnabled	41
4.9.5.3	m_bLoosingFrames	41
4.9.5.4	m_dMaxFps	42
4.9.5.5	m_eFlags	42
4.9.5.6	m_eReqState	42
4.9.5.7	m_eState	42
4.9.5.8	m_hMemoryPool	42
4.9.5.9	m_iFirstLostFrame	42
4.9.5.10	m_iMaxOutputDataCount	42
4.9.5.11	m_lock	42
4.9.5.12	m_lockOutputData	42
4.9.5.13	m_mapObjectOutputData	42
4.9.5.14	m_mapOutputData	42
4.9.5.15	m_pMyChain	43

4.9.5.16	<code>m_pThread</code>	43
4.9.5.17	<code>m_queFpsCounter</code>	43
4.9.5.18	<code>m_semStartStop</code>	43
4.9.5.19	<code>m_sErrorMessage</code>	43
4.9.5.20	<code>m_setPrecedessors</code>	43
4.9.5.21	<code>m_setSuccessors</code>	43
4.9.5.22	<code>m_timerFps</code>	43
4.10	<code>CxChainManager</code> Class Reference	43
4.10.1	Detailed Description	45
4.10.2	Member Function Documentation	46
4.10.2.1	<code>branchAfter(CxChain *pChain, CxChainable *pChainable, CxChainable *pAfter)</code>	46
4.10.2.2	<code>cameraSource(XICORE_HANDLE hCamera)</code>	46
4.10.2.3	<code>closeChain(CxChain *pChain)</code>	46
4.10.2.4	<code>connectNonChainable(CxChain *pChain, CxChainable *pChainable, QObject *pListener, const char *on)</code>	46
4.10.2.5	<code>createCameraSource(XICORE_HANDLE hCamera)</code>	46
4.10.2.6	<code>createImageSource(IxImageDataLoader *pLoader)</code>	47
4.10.2.7	<code>createStaticImageSource(CxImageData *pImageData, const QString &sTitle)</code>	47
4.10.2.8	<code>duplicateChainBranch(CxChain *pChain, CxChainable *pRoot, CxChainable *pLast, CxChainable *pNewRoot, bool bIncludeViews)</code>	47
4.10.2.9	<code>duplicateChainWithViews(CxChain *pChain, CxChainable *pRoot, CxChainable *pNewRoot, bool bIncludeViews)</code>	47
4.10.2.10	<code>importToExistingChain(CxChain *pChain, const QDomDocument &aDocChain)</code>	47
4.10.2.11	<code>insertAfter(CxChain *pChain, CxChainable *pChainable, CxChainable *pAfter)</code>	47
4.10.2.12	<code>insertBefore(CxChain *pChain, CxChainable *pChainable, CxChainable *pBefore)</code>	48
4.10.2.13	<code>remove(CxChain *&pChain, CxChainable *pChainable, bool bDeleteBranchWhenNoListener=false)</code>	48
4.10.3	Member Data Documentation	48
4.10.3.1	<code>m_lock</code>	48
4.10.3.2	<code>m_mapCameraSources</code>	48
4.10.3.3	<code>m_setChainables</code>	48
4.10.3.4	<code>m_setChains</code>	48
4.11	<code>CxCircularBuffer</code> Class Reference	48
4.11.1	Detailed Description	50
4.11.2	Member Function Documentation	50
4.11.2.1	<code>init(QString *pErrMsg=NULL)</code>	50
4.11.2.2	<code>loadedData(qint32 iSeqIdx=0)</code>	50
4.12	<code>CxClickableLabel</code> Class Reference	51
4.13	<code>CxColorPickerLabel</code> Class Reference	51
4.14	<code>CxLut::CxCompContrast</code> Class Reference	52
4.15	<code>CxCorePlugin</code> Class Reference	53
4.15.1	Detailed Description	53
4.16	<code>CxHsiCameraParams::CxCorrectionMatrix</code> Class Reference	53
4.17	<code>CxDemosaic</code> Class Reference	54

4.17.1 Detailed Description	55
4.17.2 Member Function Documentation	55
4.17.2.1 help(QString &slImagesPath)	55
4.17.2.2 processData(IxChainData *pReceivedData)	55
4.17.2.3 setupStart()	56
4.17.2.4 setupStop()	56
4.18 CxExportToolButton Class Reference	56
4.19 CxFFT Class Reference	57
4.19.1 Detailed Description	58
4.19.2 Member Function Documentation	58
4.19.2.1 Calc1D(float *pInput, float *pReal, float *pImaginary, void *pCfg, void *pFFTin, void *pFFTout)	58
4.19.2.2 Calc1D(float *pInput, float *pReal, float *pImaginary, int iLength)	58
4.19.2.3 Calc2D(const SxPicBuf &picSrcMono, SxPicBuf &picRe, SxPicBuf &picIm, void *pCfg, void *pFFTin, void *pFFTout)	58
4.19.2.4 Calc2D(const SxPicBuf &picSrcMono, SxPicBuf &picRe, SxPicBuf &picIm)	58
4.20 CxFFTImageData Class Reference	59
4.20.1 Detailed Description	59
4.21 CxHsiCameraParams::CxFilterZone Class Reference	60
4.22 CxImageFormatManager::CxFormatItem Class Reference	60
4.23 CxGuiFunctions Class Reference	61
4.24 CxHistoData Class Reference	61
4.24.1 Detailed Description	62
4.24.2 Member Enumeration Documentation	62
4.24.2.1 ExHistoDataFlag	62
4.25 CxHsiCameraParams Class Reference	62
4.25.1 Member Enumeration Documentation	64
4.25.1.1 ExCorrectionMatrixAlgorithm	64
4.25.1.2 ExCorrectionMatrixType	64
4.26 CxHsiCorrectionImage Class Reference	64
4.27 CxImageData Class Reference	65
4.27.1 Detailed Description	66
4.27.2 Member Data Documentation	66
4.27.2.1 m_aPic	66
4.28 CxImageFormat Class Reference	66
4.28.1 Member Enumeration Documentation	67
4.28.1.1 EType	67
4.29 CxImageFormatManager Class Reference	67
4.29.1 Detailed Description	69
4.30 CxImageMetadata Class Reference	69
4.30.1 Detailed Description	70
4.31 CxImageProvider Class Reference	70

4.31.1 Detailed Description	71
4.31.2 Member Function Documentation	71
4.31.2.1 currentOutputImageInfo(SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL) 71	
4.31.2.2 init(QString *pErrMsg=NULL)	72
4.31.2.3 queryImagesCountInMemoryPool(const SxPicBufInfo &picInfo, TxMemPoolHandle hPool) 72	
4.32 CxImageSource Class Reference	72
4.32.1 Detailed Description	73
4.33 CxShadingPreset::CxItem Class Reference	73
4.34 CxLut Class Reference	73
4.34.1 Member Enumeration Documentation	75
4.34.1.1 ExLutDataType	75
4.34.1.2 ExMapMode	75
4.35 CxMdiView Class Reference	76
4.35.1 Detailed Description	77
4.35.2 Member Enumeration Documentation	77
4.35.2.1 ExViewCap	77
4.36 CxMemoryManager Class Reference	77
4.36.1 Detailed Description	79
4.36.2 Member Function Documentation	79
4.36.2.1 allocPicBuf(TxMemPoolHandle hPool, QObject *pClient, SxPicBuf &picBuf, bool bReset) 79	
4.36.2.2 bufferInfo(TxMemPoolHandle hPool, SxPicBufInfo &picInfo)	80
4.36.2.3 buffersCount(TxMemPoolHandle hPool)	80
4.36.2.4 checkLimits(ExAppOperation eOperation, XICORE_HANDLE hCam, QString &sMessage, const QString	
4.36.2.5 createPool()	80
4.36.2.6 deletePool(TxMemPoolHandle hPool)	80
4.36.2.7 finishAll	81
4.36.2.8 freeBuffersCount(TxMemPoolHandle hPool)	81
4.36.2.9 freePicBuf(SxPicBuf &picBuf)	81
4.36.2.10 freePool(TxMemPoolHandle hPool)	81
4.36.2.11 initPool(TxMemPoolHandle hPool, const SxPicBufInfo &picInfo, int iCount, int &rAllocatedCount) 81	
4.36.3 Member Data Documentation	82
4.36.3.1 m_lock	82
4.36.3.2 m_setMemPools	82
4.37 CxHsiCameraParams::CxOpticalComponent Class Reference	82
4.37.1 Detailed Description	82
4.38 CxPicBufAPI Class Reference	83
4.38.1 Member Function Documentation	85
4.38.1.1 AllocPicBufFromHeap(SxPicBuf &pic, quint32 uiWidth, quint32 uiHeight, quint32 uiChannels, ExImageD	
4.38.1.2 AllocPicBufFromPool(TxMemPoolHandle hPool, QObject *pMemPoolClient, SxPicBuf &aPic, const SxPi	
4.38.1.3 ConvertFFTTToRGB(const SxPicBuf &picRe, const SxPicBuf &picIm, SxPicBuf &picRGB, SxPicBuf *pTm	

4.38.1.4	ConvertFFTTToRGB24(const SxPicBuf &picRe, const SxPicBuf &picIm, SxPicBuf &picRGB24, SxPicBuf &picRGB32)	86
4.38.1.5	ConvertFFTTToRGB32(const SxPicBuf &picRe, const SxPicBuf &picIm, SxPicBuf &picRGB32, SxPicBuf &picRGB565)	86
4.38.1.6	ConvertFFTTToRGB565(const SxPicBuf &picRe, const SxPicBuf &picIm, SxPicBuf &picRGB565, SxPicBuf &picRGB24)	86
4.38.1.7	ConvertPicToDifBpc(const SxPicBuf &picSrc, SxPicBuf &picDst, bool bStretchValues=false)	86
4.38.1.8	DetectColorBitDepth(const SxPicBuf &picBuf, uint &uiBitDepth)	86
4.38.1.9	FilterAvgOnMosaicPicBuf(const SxPicBuf &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight)	86
4.38.1.10	FilterMaxOnMosaicPicBuf(const SxPicBuf &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight)	86
4.38.1.11	FilterMinOnMosaicPicBuf(const SxPicBuf &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight)	86
4.38.1.12	FilterNormalizeOnMosaicPicBuf(const SxPicBuf &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight)	86
4.38.1.13	FreePicBuf(SxPicBuf &pic, bool bReset=true)	87
4.38.1.14	Median(const SxPicBuf &picSrc, SxPicBuf &picDst, int iKernelRadius, void *pBuffer=NULL)	88
4.38.1.15	ReverseColOrder(SxPicBuf &pic)	88
4.38.1.16	ReverseRowOrder(SxPicBuf &pic)	88
4.39	CxPicBufMeas Class Reference	88
4.40	CxHistoData::CxPixelStats Class Reference	89
4.40.1	Detailed Description	89
4.41	CxPlugin Class Reference	89
4.41.1	Detailed Description	90
4.42	CxProfileData Class Reference	90
4.42.1	Detailed Description	91
4.43	CxLut::CxPseudoColorTable Class Reference	91
4.44	CxRgb32Convertor Class Reference	91
4.44.1	Detailed Description	93
4.44.2	Member Function Documentation	93
4.44.2.1	processData(IxChainData *pReceivedData)	93
4.44.2.2	setupStop()	93
4.44.3	Member Data Documentation	93
4.44.3.1	m_lockLuts	93
4.45	CxRTTI Class Reference	94
4.45.1	Detailed Description	95
4.45.2	Member Function Documentation	95
4.45.2.1	createInstanceAndCast(const QString &sClass)	95
4.46	CxSegmentedControl Class Reference	95
4.46.1	Detailed Description	96
4.47	CxShadingPreset Class Reference	96
4.48	CxStaticImageSource Class Reference	97
4.48.1	Detailed Description	98
4.48.2	Member Function Documentation	98
4.48.2.1	currentOutputImageInfo(SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL)	98
4.48.2.2	data(QObject *pReceiver, qint32 iFrameNo)	99

4.48.2.3	run()	99
4.48.2.4	saveSettings(QDomDocument &aDoc, QDomElement &aElm) const	99
4.49	CxTuningBar Class Reference	99
4.49.1	Detailed Description	101
4.50	CxUnits Class Reference	101
4.50.1	Detailed Description	102
4.50.2	Member Enumeration Documentation	102
4.50.2.1	ExRatio	102
4.51	CxValue Class Reference	102
4.51.1	Detailed Description	103
4.51.2	Member Enumeration Documentation	103
4.51.2.1	Type	103
4.52	CxValueRegInvalidator Class Reference	103
4.52.1	Detailed Description	104
4.53	CxValueRegister Class Reference	104
4.53.1	Detailed Description	104
4.54	CxVecObjContainer Class Reference	104
4.54.1	Detailed Description	105
4.55	CxVectorObject Class Reference	105
4.55.1	Detailed Description	107
4.55.2	Member Enumeration Documentation	107
4.55.2.1	ExCameraFormatChangeBehavior	107
4.56	CxVectorObjectEllipse Class Reference	107
4.56.1	Detailed Description	108
4.57	CxVectorObjectLine Class Reference	108
4.57.1	Detailed Description	109
4.57.2	Member Enumeration Documentation	109
4.57.2.1	ExLineStyle	109
4.58	CxVectorObjectRect Class Reference	110
4.58.1	Detailed Description	111
4.59	CxVideoCodecManager Class Reference	112
4.59.1	Detailed Description	113
4.60	CxVideoSource Class Reference	113
4.60.1	Detailed Description	114
4.60.2	Member Function Documentation	114
4.60.2.1	currentOutputImageInfo(SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL)	114
4.60.2.2	run()	114
4.60.2.3	saveSettings(QDomDocument &aDoc, QDomElement &aElm) const	115
4.61	CxVideoToFile Class Reference	115
4.61.1	Detailed Description	116

4.61.2	Member Function Documentation	116
4.61.2.1	init(QString *pErrMsg=NULL)	116
4.62	CxViewImageData Class Reference	116
4.62.1	Detailed Description	117
4.62.2	Constructor & Destructor Documentation	117
4.62.2.1	CxViewImageData(CxImageData *pOrigData, SxPicBuf *pViewPicBufRgb32)	117
4.62.3	Member Data Documentation	117
4.62.3.1	m_pOrigData	117
4.63	CxViewImageReceiver Class Reference	118
4.63.1	Detailed Description	119
4.63.2	Member Function Documentation	119
4.63.2.1	configurationWidget()	119
4.63.2.2	processData(IxChainData *pReceivedData)	120
4.63.2.3	queryImagesCountInMemoryPool(const SxPicBufInfo &picInfo, TxMemPoolHandle hPool)	120
4.64	CxViewReceiverConfWidget Class Reference	120
4.65	CxHsiCameraParams::CxVirtualBand Class Reference	121
4.65.1	Detailed Description	121
4.66	IxAppDelegate Class Reference	121
4.66.1	Detailed Description	123
4.66.2	Member Enumeration Documentation	123
4.66.2.1	ExSound	123
4.66.3	Member Function Documentation	124
4.66.3.1	capturePictureWithGUI(QWidget *pParent, XICORE_HANDLE hCamera, const QString &sWndTitle, con	
4.67	IxChainData Class Reference	124
4.67.1	Detailed Description	124
4.68	IxChainDataLoader Class Reference	125
4.68.1	Detailed Description	125
4.68.2	Member Function Documentation	125
4.68.2.1	loadedData(qint32 iSeqIdx=0)=0	125
4.69	IxChainDataSaver Class Reference	125
4.69.1	Detailed Description	126
4.70	IxImageDataLoader Class Reference	126
4.70.1	Detailed Description	126
4.71	IxImageDataSaver Class Reference	127
4.71.1	Detailed Description	127
4.72	IxVideoAccessor Class Reference	127
4.72.1	Detailed Description	128
4.73	IxVideoDecoder Class Reference	128
4.73.1	Detailed Description	128
4.73.2	Member Function Documentation	128

4.73.2.1	decode(qint32 iFrameIndex, IxVideoAccessor *pVideo, SxPicBuf &rDecodedPic)=0	128
4.74	IxVideoEncoder Class Reference	128
4.74.1	Detailed Description	129
4.74.2	Member Function Documentation	129
4.74.2.1	encode(const SxPicBuf &rPic, void *pEncoded, size_t &ruiSize, bool *pbKeyFrame=NULL)=0	129
4.75	CxHsiCameraParams::SxBandIdx Struct Reference	129
4.76	SxCameraInfo Struct Reference	129
4.76.1	Detailed Description	130
4.77	CxVideoCodecManager::SxCodecItem Struct Reference	130
4.78	CxCameraParam::SxEnumItem Struct Reference	130
4.78.1	Detailed Description	131
4.79	CxMemoryManager::SxLimitDef Struct Reference	131
4.80	SxPicBuf Struct Reference	131
4.80.1	Detailed Description	132
4.80.2	Member Enumeration Documentation	132
4.80.2.1	ExDataOwner	132
4.80.3	Member Data Documentation	132
4.80.3.1	m_hMemoryPoolOwner	132
4.80.3.2	m_iFrameNo	132
4.80.3.3	m_iXiApiFrameNo	133
4.80.3.4	m_pData	133
4.80.3.5	m_pMemoryPoolClient	133
4.80.3.6	m_uiAllocSize	133
4.80.3.7	m_uiTimeStamp	133
4.81	SxPicBufInfo Class Reference	133
4.81.1	Detailed Description	134
4.81.2	Member Data Documentation	134
4.81.2.1	m_eColorFilterArray	134
4.81.2.2	m_eDataType	134
4.81.2.3	m_eFormat	134
4.81.2.4	m_uiBpc	134
4.81.2.5	m_uiChannels	134
4.81.2.6	m_uiHeight	134
4.81.2.7	m_uiStride	134
4.81.2.8	m_uiWidth	134
5	File Documentation	135
5.1	Chainable.h File Reference	135
5.1.1	Macro Definition Documentation	135
5.1.1.1	CHAINABLE_RUN_END_WHILE	135

5.1.1.2	CHAINABLE_RUN_WHILE_TRUE	136
---------	------------------------------------	-----

Chapter 1

Hierarchical Index

1.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

CxHsiCameraParams::CxBand	9
CxBinaryLayerContainer	11
CxCameraParam	15
CxLut::CxCompContrast	52
CxHsiCameraParams::CxCorrectionMatrix	53
CxFFT	57
CxHsiCameraParams::CxFilterZone	60
CxImageFormatManager::CxFormatItem	60
CxGuiFunctions	61
CxHistoData	61
CxHsiCameraParams	62
CxHsiCorrectionImage	64
CxImageFormat	66
CxImageMetadata	69
CxShadingPreset::CxItem	73
CxLut	73
CxHsiCameraParams::CxOpticalComponent	82
CxPicBufAPI	83
CxPicBufMeas	88
CxHistoData::CxPixelStats	89
CxPlugin	89
CxCorePlugin	53
CxProfileData	90
CxLut::CxPseudoColorTable	91
CxShadingPreset	96
CxUnits	101
CxValue	102
CxValueRegInvalidator	103
CxValueRegister	104
CxVecObjContainer	104
CxHsiCameraParams::CxVirtualBand	121
IxAppDelegate	121
IxChainDataLoader	125
IxImageDataLoader	126
CxCircularBuffer	48
IxChainDataSaver	125
IxImageDataSaver	127

IxVideoAccessor	127
IxVideoDecoder	128
IxVideoEncoder	128
QGraphicsObject	
CxVectorObject	105
CxVectorObjectLine	108
CxVectorObjectRect	110
CxVectorObjectEllipse	107
QLabel	
CxClickableLabel	51
CxColorPickerLabel	51
QObject	
CxBinaryLayer	9
CxCameraManager	12
CxCameraParameters	17
CxChain	24
CxChainable	30
CxCircularBuffer	48
CxImageProvider	70
CxDemosaic	54
CxImageSource	72
CxCameraSource	20
CxStaticImageSource	97
CxVideoSource	113
CxRgb32Convertor	91
CxViewImageReceiver	118
CxVideoToFile	115
CxChainManager	43
CxCorePlugin	53
CxImageFormatManager	67
CxMemoryManager	77
CxRTTI	94
CxVideoCodecManager	112
IxChainData	124
CxImageData	65
CxFFTImageData	59
CxViewImageData	116
QToolButton	
CxExportToolButton	56
QWidget	
CxMdiView	76
CxSegmentedControl	95
CxTuningBar	99
CxViewReceiverConfWidget	120
CxHsiCameraParams::SxBandIdx	129
SxCameraInfo	129
CxVideoCodecManager::SxCodecItem	130
CxCameraParam::SxEnumItem	130
CxMemoryManager::SxLimitDef	131
SxPicBufInfo	133
SxPicBuf	131

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

CxHsiCameraParams::CxBand	9
CxBinaryLayer	9
CxBinaryLayerContainer	11
CxCameraManager	
The class that provides all necessary methods to access cameras and to set up camera parameters	12
CxCameraParam	15
CxCameraParameters	17
CxCameraSource	
The implementation of image source for Ximea cameras	20
CxChain	
This class implements a processing chain which is together with the CxChainable class (chainable objects) the keystones of the XIMEA CamTool	24
CxChainable	
A base class of all chainable objects	30
CxChainManager	
The class that creates, manages and manipulates with chains	43
CxCircularBuffer	48
CxClickableLabel	51
CxColorPickerLabel	51
CxLut::CxCompContrast	52
CxCorePlugin	
The core plugin that registers all core classes (image loaders/savers etc.)	53
CxHsiCameraParams::CxCorrectionMatrix	53
CxDemosaic	
The class performs a debayering on images in RAW format	54
CxExportToolButton	56
CxFFT	
The class covers several 1D and 2D FFT functions	57
CxFFTImageData	
Implementation of CxImageData that can keep the FFT on an image. It consists of real and imaginary parts both of type <code>extypeInt32</code>	59
CxHsiCameraParams::CxFilterZone	60
CxImageFormatManager::CxFormatItem	60
CxGuiFunctions	61
CxHistoData	
Class holding measured histogram and pixel statistics data	61
CxHsiCameraParams	62

CxHsiCorrectionImage	64
CxImageData	
The class implements an image data object exchanged by CxChainable objects	65
CxImageFormat	66
CxImageFormatManager	
Keeps list of supported file formats using CxRTTI , creates loader and saver objects	67
CxImageMetadata	
Metadata stored next to image pixel data	69
CxImageProvider	
The base class for any chainable object which provides images	70
CxImageSource	
The base class of any image provider which is an image source, i.e. it has no predecessors and it provides images	72
CxShadingPreset::CxItem	73
CxLut	73
CxMdiView	
The base of all MDI subwindows in the application	76
CxMemoryManager	
The singleton class for management of memory pools	77
CxHsiCameraParams::CxOpticalComponent	
< Calibration data for used rejection filter (on camera, may not be present)	82
CxPicBufAPI	83
CxPicBufMeas	88
CxHistoData::CxPixelStats	
Class for measurement of data in the ROI/image per channel	89
CxPlugin	
Base class of all plugins. For more informations about plugin development see the plugins	89
CxProfileData	
Class holding measured line profile data	90
CxLut::CxPseudoColorTable	91
CxRgb32Convertor	
The class converts all incoming image data to RGB32 format."	91
CxRTTI	
Our RTTI builds upon Qt meta object system, and just adds class discovery	94
CxSegmentedControl	
Set of exclusively checkable buttons next to each other, inspired by OS X	95
CxShadingPreset	96
CxStaticImageSource	
The image source of static image	97
CxTuningBar	
Our own implementation of slides, working with CxValue types etInt and etFloat	99
CxUnits	
Class the defines all known units and measures, and includes conversions between them	101
CxValue	
Class for storing values similar to variant. Used in camera parameters	102
CxValueRegInvalidator	
Class describing invalidator connection	103
CxValueRegister	
This class is responsible for the communitaion with the camera	104
CxVecObjContainer	
Our special class above QGraphicsScene maintains list of our vector objects and assigns IDs to them	104
CxVectorObject	
Base class for all our vector objects in QGraphicsScene. Needed to send signals away	105
CxVectorObjectEllipse	
Object representing the ellipse, using the functionality implemented in rectangle	107
CxVectorObjectLine	
Object representing the single line, from start point to end point	108

CxVectorObjectRect	Object representing the rectangular object. Implements mouse resizing and caption	110
CxVideoCodecManager	The class that manages video codecs used by video loaders and savers	112
CxVideoSource	Provides image sequences	113
CxVideoToFile	The chainable class that allows to save images coming from camera to a video file	115
CxViewImageData	Implementation of CxImageData intended for displaying in views	116
CxViewImageReceiver	This class is primarily used by images views to receive a image data from a chain	118
CxViewReceiverConfWidget		120
CxHsiCameraParams::CxVirtualBand	< Band record inside <virtual_bands>	121
IxAppDelegate	Class for plugins to communicate with the application. Use IxAppDelegate::instance() to get the pointer	121
IxChainData	The base class for all data exchanged by CxChainable objects	124
IxChainDataLoader	The base class for any kind of chain data loading	125
IxChainDataSaver	The base class for any kind of class for saving IxChainData	125
IxImageDataLoader	The base class for image data loading	126
IxImageDataSaver	The base class for any kind of image data saver	127
IxVideoAccessor	The class for accessing the raw video data buffers	127
IxVideoDecoder	The base class of all video decoders	128
IxVideoEncoder	The base class of all video encoders	128
CxHsiCameraParams::SxBandIdx		129
SxCameraInfo	Information about camera device	129
CxVideoCodecManager::SxCodecItem		130
CxCameraParam::SxEnumItem	Struct holding all info for each item of enumerated parameter (ecpEnum)	130
CxMemoryManager::SxLimitDef		131
SxPicBuf	Picture buffer	131
SxPicBufInfo	Information about picture buffer	133

Chapter 3

File Index

3.1 File List

Here is a list of all documented files with brief descriptions:

AppDelegate.h	??
BinaryLayer.h	??
BinaryLayerContainer.h	??
CameraManager.h	??
CameraParameters.h	??
CameraSource.h	??
Chain.h	??
Chainable.h	135
ChainManager.h	??
CircularBuffer.h	??
ClickableLabel.h	??
CorePlugin.h	??
Data.h	??
DataLoader.h	??
DataSaver.h	??
Demosaic.h	??
ExportToolButton.h	??
FFT.h	??
GuiFunctions.h	??
HsiCameraParams.h	??
ImageData.h	??
ImageFormat.h	??
ImageFormatManager.h	??
ImageMetadata.h	??
Lut.h	??
MdiView.h	??
MemoryManager.h	??
PicBuf.h	??
PicBufAPI.h	??
PicBufMacros.h	??
PicBufMeas.h	??
Plugin.h	??
Rgb32Convertor.h	??
SegmentedControl.h	??
StaticImageSource.h	??
TuningBar.h	??
Units.h	??
ValueTypes.h	??

VectorContainer.h	??
VectorObject.h	??
VideoAccessor.h	??
VideoCodecManager.h	??
VideoDecoder.h	??
VideoEncoder.h	??
VideoSource.h	??
VideoToFile.h	??
ViewImageReceiver.h	??
xiCore.h	??
xiCoreGlobal.h	??
xiCoreGuiGlobal.h	??
xiCoreLog.h	??
xiCoreVersion.h	??
xiCrash.h	??
xiRtti.h	??

Chapter 4

Class Documentation

4.1 CxHsiCameraParams::CxBand Class Reference

Public Member Functions

- bool **isValid** () const
- double **peakInRange** (int iMinNm, int iMaxNm) const
Returns peak value in a specified range.
- QString **name** () const
Returns band name as "123 nm".

Public Attributes


- bool **m_bSelected**
Flags if the signal of this band is within specifications (true) or not (false). Usage of non-selected bands will result in wrong spectra.
- QVector< double > **m_vecPeakWaveLengths**
Central wavelength of peaks [nm]. At least one peak should be defined.
- QVector< double > **m_vecCalibData**
Band response data, tag <response_composition> (not required for images)

The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.2 CxBinaryLayer Class Reference

Inheritance diagram for CxBinaryLayer:



Signals

- void **changed** (TxBinaryLayerID idLayer)
- void **visibilityChanged** (TxBinaryLayerID idLayer)

Public Member Functions

- **CxBinaryLayer** (uint uiWidth, uint uiHeight)
- TxBinaryLayerID **id** () const
ID of the layer.
- void **setId** (TxBinaryLayerID ild)
Used by [CxBinaryLayerContainer](#).
- QString **name** () const
A custom name of the layer (default is "Layer #id").
- void **setName** (const QString &sName)
Sets the custom layer name.
- const [SxPicBuf](#) & **picbuf** () const
- [SxPicBuf](#) & **picbuf** ()
Use with caution!! Change only the buffer's data (m_pData) but nothing more!!!!
- void **reset** ()
Resets the layer's data to zeros.
- void **release** ()
Releases the layer's data.
- bool **alloc** (uint uiWidth, uint uiHeight)
Allocs the layer's data. If the data exists then the allocation is not done and the method returns false.
- const QColor **color** () const
- void **setColor** (const QColor &color)
- int **alpha** () const
- void **setAlpha** (int iAlpha)
- bool **visible** () const
- void **setVisible** (bool bVisible)
- void **notifyChanged** ()
When the layer's data is modified, call this method.

Static Public Member Functions

- static bool **MakeColoredBinary** (const [CxBinaryLayer](#) *pLayer, [SxPicBuf](#) &picRgb32)
Converts the binary layer to colored RGB32 image with premultiplied rgb values (see QImage::Format_ARGB32_↔ Premultiplied)
- static bool **BlendRGBWithBinary** (const [SxPicBuf](#) &picRgb, const [CxBinaryLayer](#) *pLayer, [SxPicBuf](#) &pic↔ BlendedRgb)
Blends the given RGB image (RGB32, RGB24, RGB565) with the binary layer. The resulting picBlendedRgb is the image of the same format as picRgb.

Protected Member Functions

- [CxBinaryLayer](#) & **operator=** (const [CxBinaryLayer](#) &rOther)
Prevent to use the operator=.
- void **setAsShallowCopyOf** (const [CxBinaryLayer](#) &rOther)

Private Attributes

- TxBinaryLayerID **m_ild**
- QString **m_sName**
- SxPicBuf **m_binary**
- QColor **m_color**
- bool **m_bVisible**
- bool **m_bChanged**

The documentation for this class was generated from the following file:

- BinaryLayer.h

4.3 CxBinaryLayerContainer Class Reference

Public Member Functions

- TxBinaryLayerID **addLayer** (CxBinaryLayer *pLayer, bool bSetAutoColor=true)
Add new layer and assign ID.
- void **deleteLayer** (TxBinaryLayerID idLayer)
Deletes (and deallocates) the layer with given ID.
- void **deleteAllLayers** ()
Deletes all layers from container.
- CxBinaryLayer * **removeLayer** (TxBinaryLayerID idLayer)
Remove the layer from container and returns its pointer.
- QList< TxBinaryLayerID > **listLayers** () const
List of all layers.
- CxBinaryLayer * **layerWithId** (TxBinaryLayerID idLayer) const
Returns a pointer to the layer with given id. If it not exists then it returns NULL.
- bool **hasVisibleLayer** () const
Returns true if at least one layer is visible.
- QList< TxBinaryLayerID > **listVisibleLayers** () const
List of all visible layers.
- void **saveLayersVisibilityState** ()
Saves the per-layer visibility state for later.
- void **restoreLayersVisibilityState** ()
Restores the layers visibility state as saved in `saveLayersVisibilityState()`;
- void **hideAllLayers** ()
Sets all layers as invisible.

Private Member Functions

- QColor **autoColor** ()

Private Attributes

- TxBinaryLayerMap **m_mapLayers**
- QColor **m_colCurrent**
For auto color.
- int **m_idMax**
ID generator.
- QMap< TxBinaryLayerID, bool > **m_mapVisibilityStates**

4.3.1 Member Function Documentation

4.3.1.1 void CxBinaryLayerContainer::deleteAllLayers ()

Deletes all layers from container.

See also

[deleteLayer\(\)](#)


The documentation for this class was generated from the following file:

- BinaryLayerContainer.h

4.4 CxCameraManager Class Reference

The class that provides all necessary methods to access cameras and to set up camera parameters.

Inheritance diagram for CxCameraManager:



Public Slots

- void [finishAll](#) ()
Stops all cameras and closes them. This method should be called on application exit.

Signals

- void [acquisitionStarted](#) (XICORE_HANDLE hCamera)
The signal is emitted when the acquisition on the given camera was started.
- void [acquisitionFinished](#) (XICORE_HANDLE hCamera)
The signal is emitted when the acquisition on the given camera was finished.
- void [cameraOpened](#) (XICORE_HANDLE hCamera)
The signal is emitted when camera is opened.
- void [cameraClosed](#) (XICORE_HANDLE hCamera)
The signal is emitted when the camera was closed.
- void [cameraLost](#) (XICORE_HANDLE hCamera)
The signal is emitted when camera disconnected during app session.
- void [cameraAcqError](#) (XICORE_HANDLE hCamera, int iErrorCode, const QString &slnFunction)
The signal is emitted when camera acquisition fails.
- void [imageFormatChanging](#) (XICORE_HANDLE hCamera)
The signal is emitted when format of camera image (dimension, components..) is to be changed.
- void [imageFormatChanged](#) (XICORE_HANDLE hCamera)
The signal is emitted when format of the camera image (dimension, komponents..) changed.

Public Member Functions

- [~CxCameraManager \(\)](#)
Destructor.
- uint [cameraCount \(\)](#) const
Gets the number of available cameras connected to PC.
- bool [cameraInfo](#) (uint uiCameraIdx, [SxCameraInfo](#) &aInfo) const
Gets the info about camera with given index uiCameraIdx.
- bool [cameraSerialNumber](#) (uint uiCameraIdx, QString &sSerialNumber) const
- XICORE_HANDLE [openCamera](#) (uint uiCameraIdx, [SxCameraInfo](#) *pPreloadedInfo=NULL)
Opens a camera with given index and returns its handle.
- uint [openedCameraCount \(\)](#)
Gets the number of opened cameras.
- QList< XICORE_HANDLE > [openedCameras \(\)](#)
Gets a list of handles to the opened cameras.
- bool [isCameraOpened](#) (XICORE_HANDLE hCamera)
- bool [closeCamera](#) (XICORE_HANDLE hCamera)
- bool [startAcquisition](#) (XICORE_HANDLE hCamera)
- bool [stopAcquisition](#) (XICORE_HANDLE hCamera)
- uint [runningCameraCount \(\)](#)
- bool [isCameraRunning](#) (XICORE_HANDLE hCamera)
- QList< XICORE_HANDLE > [runningCameras \(\)](#)
- [CxCameraParameters](#) * [cameraParameters](#) (XICORE_HANDLE hCamera)
Gets the camera parameters.
- bool [readCameraHsiParams](#) (XICORE_HANDLE hCamera, const [SxCameraInfo](#) &info, [CxHsiCameraParams](#) *pHsiParams) const
Reads HSI calibration XML from camera.
- void [notifyCameraFormatChanging](#) (XICORE_HANDLE hCamera)
Call this method to notify the rest of the application that the camera format (image dimensions, components, ...) are to be changed.
- void [notifyCameraFormatChanged](#) (XICORE_HANDLE hCamera)
Call this method to notify rest of the application that the camera format (image dimensions, components, ...) has been changed.
- void [notifyCameraLost](#) (XICORE_HANDLE hCamera)
*Call this method to notify rest of the application that the camera just disconnected */.*
- void [notifyCameraAcqError](#) (XICORE_HANDLE hCamera, int iErrorCode, const QString &sInFunction)
*Call this method to notify rest of the application that the camera acquisition is failing */.*
- bool [cameraPicbufInfo](#) (XICORE_HANDLE hCamera, [SxPicBufInfo](#) &picInfo)
Gets the picbuf info currently produced by the given camera.
- XICORE_HANDLE [cameraHandleBySerialId](#) (const QString &sSerialId)
Gets the camera handle by the serial id of the camera.
- bool [newProcessChainEnabled \(\)](#) const
State of XI_PRM_NEW_PROCESS_CHAIN_ENABLE parameter.
- void [setNewProcessChainEnabled](#) (bool bEnable)
- bool [newProcessChainSupported \(\)](#) const
- void [saveSettings](#) (QSettings *pSettings)
store some options to app settings
- void [loadSettings](#) (QSettings *pSettings)
load some options from app settings
- void [createCameraSimulator \(\)](#)
- void [disconnectCameraSimulator \(\)](#)
- void [createHsiCameraSimulator \(\)](#)

Static Public Member Functions

- static [CxCameraManager](#) * [instance](#) ()
Gets the instance if this object.

Private Member Functions

- [CxCameraManager](#) ()
Private constructor (the class is singleton).

Private Attributes

- QMap< XICORE_HANDLE, [CxCameraParameters](#) * > [m_mapCameras](#)
- QSet< XICORE_HANDLE > [m_setRunningCameras](#)
- QMutex [m_lock](#)
- XICORE_HANDLE [m_hLastCameraLost](#)
- XICORE_HANDLE [m_hLastCameraAcqError](#)
- QString [m_sLastCameraAcqErrorFunc](#)
- bool [m_bEnableNewProcessChain](#)
- bool [m_bNewProcessingChainSupported](#)

4.4.1 Detailed Description

The class that provides all necessary methods to access cameras and to set up camera parameters.

4.4.2 Member Function Documentation

4.4.2.1 void [CxCameraManager::cameraClosed](#) (XICORE_HANDLE *hCamera*) [signal]

The signal is emitted when the camera was closed.

Remarks

The camera handle *hCamera* is invalid in the moment of emitting of this signal, so any attempt to communication with the camera fails.

4.4.2.2 bool [CxCameraManager::cameraInfo](#) (uint *uiCameraIdx*, [SxCameraInfo](#) & *ainfo*) const

Gets the info about camera with given index *uiCameraIdx*.

Remarks

The index *uiCameraIdx* is related to cameras count returned by the [cameraCount\(\)](#) method (NOT by the [openedCameraCount\(\)](#) method!).

4.4.3 Member Data Documentation

4.4.3.1 XICORE_HANDLE [CxCameraManager::m_hLastCameraAcqError](#) [private]

Do not emit the cameraAcqError signal more then once

4.4.3.2 XICORE_HANDLE CxCameraManager::m_hLastCameraLost [private]

Do not emit the cameraLost signal more then once

4.4.3.3 QMutex CxCameraManager::m_lock [private]

The lock of the manager.

4.4.3.4 QMap<XICORE_HANDLE, CxCameraParameters*> CxCameraManager::m_mapCameras [private]

Camera parameters by the camera handle

4.4.3.5 QSet<XICORE_HANDLE> CxCameraManager::m_setRunningCameras [private]

Set of camera where the acquisition is running.

The documentation for this class was generated from the following file:

- CameraManager.h

4.5 CxCameraParam Class Reference

Classes

- struct [SxEnumItem](#)
Struct holding all info for each item of enumerated parameter (ecpEnum)

Public Types

- enum **ExType** {
 ecpInvalid, **ecpInt**, **ecpFloat**, **ecpBool**,
 ecpString, **ecpEnum**, **ecpCommand** }
- enum **ExRepresentation** { **erepDefault**, **erepLogarithmic** }
- enum [ExVisibility](#) {
 evisInvisible = 0, **evisBeginner** = 1, **evisExpert** = 2, **evisGuru** = 3,
 evisMaxLevel = evisGuru }

Camera parameter visibility level (selected in app GUI)

Public Member Functions

- **CxCameraParam** (ExType eType, const QString &sXiApiName, const QString &sCaption, const QString &s↵
 UIGroup, [ExVisibility](#) eVisibility)
- [CxCameraParam](#) ()
This default constructor should not be used.
- void [readInfo](#) (XICORE_HANDLE hCamera)
Reads m_vMin, m_vMax and m_vIncrement from xiApi, should be called only once.
- bool [readFromCamera](#) (XICORE_HANDLE hCamera)
Reads current value from the camera.
- bool [writeToCamera](#) (XICORE_HANDLE hCamera) const
Sends the current value to camera.
- void [filterEnumValues](#) (XICORE_HANDLE hCamera)

- Filter out enum values.*

 - void [limitEnumsWorkaround](#) (XICORE_HANDLE hCamera)
- Filter out enum values (old hack)*

 - bool [updateIfInvalidatedWithParam](#) (XICORE_HANDLE hCamera, const QString &sInvalidator)

Returns true if it has sInvalidator among its invalidators.
- bool [isVisibleInGUI](#) () const

Parameter should be visible in Camera settings.
- bool [isVisibleForLevel](#) ([ExVisibility](#) eLevel) const

Same as isVisibleInGUI, but checks visibility level too.
- bool [belongsToGroup](#) (const QString &sUIGroup, [ExVisibility](#) eLevel=evisMaxLevel) const

Same as isVisibleForLevel, but verifies the group membership too.
- bool [isInvalidatedByParam](#) (const QString &sParam) const

Returns true if this parameter is invalidated by another parameter.
- [CxValue](#) [paramValue](#) () const

returns current value. When enum, sends the value, not selection index.
- int [indexOfEnumValue](#) (const [CxValue](#) &vValue) const

Returns index if value is among this param's enum values.
- bool [setParamValue](#) (const [CxValue](#) &vValue)

sets parameter value. When enum, sets the correct selection index.

Public Attributes

- QString [m_sXiApiName](#)

Parameter name as in xiApi (e.g. XI_PRM_EXPOSURE)
- QString [m_sCaption](#)

User-friendly name of this parameter.
- QString [m_sTooltip](#)

Some longer text about this parameter.
- QString [m_sUIGroup](#)

Name of the section in camera settings this parameter should belong to.
- bool [m_bReadOnly](#)

If we can write this parameter to camera.
- [ExVisibility](#) [m_eVisible](#)

Tells that this camera parameter should be present in GUI.
- bool [m_bNeedsRegularRefresh](#)

We should poll this parameter to get new updates automatically (even if there is no invalidator)
- bool [m_bRequiresAcqStop](#)

Acquisition loop has to be stopped before settings this parameter.
- [ExRepresentation](#) [m_eRepresentation](#)

How the parameter should look like in GUI.
- [ExType](#) [m_eType](#)

Type of the parameter: int, float, bool, string, enum, command...
- [CxValue](#) [m_vValue](#)

Current value of the parameter, type is selected according to m_eType.
- [CxValue](#) [m_vMin](#)

Minimum value as read from xiApi or defined in XML.
- [CxValue](#) [m_vMax](#)

Maximum value as read from xiApi or defined in XML.
- [CxValue](#) [m_vIncrement](#)

Value increment as read from xiApi or defined in XML.

- QVector< [SxEnumItem](#) > [m_vecEnumValues](#)
Special case for enumerated parameter, selection stored in m_vValue.intVal. Filtered using XI_PRM_INFO_SETTABLE.
- CxValue [m_vAppDefault](#)
Default value as defined in XML (may be undefined - type is etInvalid)
- bool [m_bCameraError](#)
There was an error when getting the value of this parameter. Probably not supported by camera.

Private Attributes

- CxValueRegister [m_aValueRegister](#)
- CxValueRegister [m_aMinRegister](#)
- CxValueRegister [m_aMaxRegister](#)
- CxValueRegister [m_alncRegister](#)
- QVector< [SxEnumItem](#) > [m_vecEnumValuesUnfiltered](#)
All enum values as defined in XML. Filtered values according to camera caps are in m_vecEnumValues.

Friends


- class **CxCameraCapsLoader**
- class **CxCameraParameters**

The documentation for this class was generated from the following file:

- CameraParameters.h

4.6 CxCameraParameters Class Reference

Inheritance diagram for CxCameraParameters:



Signals

- void [parametersNeedRefresh](#) (const QStringList &lstParamNames)
Signal for GUI that it should refresh controls for these parameters.

Public Member Functions

- **CxCameraParameters** (XICORE_HANDLE hCamera, const [SxCameraInfo](#) &aCamInfo)
- XICORE_HANDLE [cameraHandle](#) () const
Returns camera handle to which this object is connected.
- bool [loadFromXMLDescription](#) ()
Clears everything and uses the definition stored in the external xml file for this camera.
- void [saveDescriptionToXML](#) (const QString &sFile)
For testing purposes, inverse to loadFromXMLDescription.

- const [SxCameraInfo](#) & **cameraInfo** () const
- const [CxHsiCameraParams](#) & **hsiParams** () const
- const QSize & [sensorResolution](#) () const
Physical resolution of the sensor.
- const QSize & [downsampledResolution](#) () const
Maximum resolution at the current downsampling (when no ROI active)
- const QRectF & [lastROI](#) () const
ROI relatively to resolution as float vals (i.e. 1.0=full width)
- void [setLastROI](#) (const QRectF &rc)
Tells us to remember provided ROI geometry for next time.
- bool [setCameraUserId](#) (const QString &sUserId)
Write custom User ID to camera flash memory.
- void [saveParamValues](#) (const QString &sStorageFolder)
Save values of all parameters to file, so that we can load it next time (called on app exit)
- void [saveParamValuesToFile](#) (const QString &sFilename)
Save values of all parameters to custom file.
- QString [saveParamValuesToString](#) ()
Save values of all parameters to string.
- bool [loadParamValues](#) (const QString &sStorageFolder, QStringList *plstLoadedParams=NULL)
Load saved values, called at the start of the application.
- bool [loadParamValuesFromFile](#) (const QString &sFilename, QStringList *plstLoadedParams=NULL, bool b↔
NotifyGUI=false)
Load saved value from custom file.
- bool [loadParamValuesFromString](#) (const QString &sString, QStringList *plstLoadedParams=NULL, bool b↔
NotifyGUI=false)
Load saved value from string.
- void [saveAdditionalSettingsToXml](#) (QDomDocument *pXmlDoc, QDomElement *pRootElement)
- void [loadAdditionalSettingsFromXml](#) (QDomElement *pRootElement)
- bool [addParam](#) (const [CxCameraParam](#) &aParam)
Adds new parameter to our list.
- [CxCameraParam](#) * [findParam](#) (const QString &sXiApiParamName)
Returns NULL when not found.
- [CxValue](#) [currentParamValue](#) (const QString &sXiApiParamName)
Refresh current value from camera and return it.
- const QStringList & [paramsOrderList](#) () const
Returns list with parameter names in the order they appeared in xml.
- bool [commitSingleParamValue](#) (const QString &sXiApiParamName, bool bRefreshCommittedParamsIn↔
Gui=false)
Adds parameter to batch and commits it right away.
- void [addParamCommitToBatch](#) (const QString &sXiApiParamName)
Marks the parameter to be sent to camera in this batch.
- bool [endCommitBatch](#) (bool bRefreshCommittedParamsInGui=false)
Sends all parameters in batch to camera, and checks invalidators.
- bool [setParamDirect](#) (const QString &sXiApiParamName, const [CxValue](#) &vValue)
Use commitSingleParamValue instead of this.
- bool [getParamDirect](#) (const QString &sXiApiParamName, [CxValue](#) &vValue) const
Note: vValue.m_eType must be correctly filler.
- void [setROIRectParams](#) (int x, int y, int w, int h, bool bCheckMinMax=false)
Sends new ROI to camera, while checking min/max/inc optionally.
- void [setROIRectParams](#) (const QRect &rc, bool bCheckMinMax=false)
Sends new ROI to camera, while checking min/max/inc optionally.

- QRect [currentROIRectParams](#) ()
Returns current ROI as rectangle.
- void [setImageFormatKeepRoi](#) (int eXImgFormat)
Switches XI_PRM_IMAGE_DATA_FORMAT and forces to keep current ROI.
- bool [rereadHsiParamsForRange](#) (int iStartNm, int iEndNm)
switch to different filter range (when lens changed, e.g. for 5x5)
- bool **setWhitelImage** (const [SxPicBuf](#) &aPic, const [CxImageMetadata](#) &aMetadata)
- bool **setDarkImage** (const [SxPicBuf](#) &aPic, const [CxImageMetadata](#) &aMetadata)
- bool **setDarkBiasImage** (const [SxPicBuf](#) &aPic, const [CxImageMetadata](#) &aMetadata)
- [CxHsiCorrectionImage](#) & **whitelImage** ()
- [CxHsiCorrectionImage](#) & **darkImage** ()
- [CxHsiCorrectionImage](#) & **darkBiasImage** ()
- ExHsiCorrection **hsiCorrection** ()
- void **setHsiCorrection** (ExHsiCorrection eCorr)
- bool **isCorrectionReady** ()
- void **notifyShadingChanged** ()

Static Public Member Functions

- static bool [queryCameraStoredInParamFile](#) (const QString &sFilename, QString *psCamera, QString *psSerial)
Returns camera model and serial stored in xml file with values.
- static bool [queryCameraStoredInParamString](#) (const QString &sString, QString *psCamera, QString *psSerial)
Returns camera model and serial stored in xml string with values.

Public Attributes

- TxShadingPresetList **m_lstShadingPresets**

Private Slots

- void [on_timerRefreshParams_update](#) ()
Timer for auto-refreshed parameters (temperature,...)

Private Member Functions

- void **addParamIfNotPresent** (const QString &sXiApiParamName, CxCameraParam::ExType eType)

Private Attributes

- QMap< QString, [CxCameraParam](#) > **m_mapParams**
Map of all supported parameters by this camera, keys are their XiApiNames.
- QStringList **m_lstParamsOrder**
Parameter names in the order they appeared in xml.
- XICORE_HANDLE **m_hCamera**
Handle to camera.
- [SxCameraInfo](#) **m_aCamInfo**
Basic camera informations.
- [CxHsiCameraParams](#) **m_aHsiParams**

- Hyperspectral information.*

 - QSize [m_sizeSensor](#)
Physical resolution of the sensor.
 - QSize [m_sizeDownsampledRes](#)
Remember the maximum resolution we can get at the current downsampling mode. It will enable us reset ROI in one go.
 - QRectF [m_rcLastROI](#)
Save ROI relatively to DownsampledRes as float vals (i.e. 1.0=full width)
 - QStringList [m_lstParamsToCommit](#)
Parameter xiApiNames in current commit batch.
 - CxHsiCorrectionImage [m_picWhite](#)
Exposure time T_{ref} .
 - CxHsiCorrectionImage [m_picDark](#)
Exposure time T_{ref} .
 - CxHsiCorrectionImage [m_picDarkBias](#)
Exposure time same as the current live image being corrected.
 - ExHsiCorrection [m_eHsiCorrection](#)
 - QMutex [m_lock](#)


The documentation for this class was generated from the following file:

- CameraParameters.h

4.7 CxCameraSource Class Reference

The implementation of image source for Ximea cameras.

Inheritance diagram for CxCameraSource:



Public Member Functions

- virtual QString [title](#) () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual int [buffersCountInMemoryPool](#) () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*
- virtual [CxChainable](#) * [clone](#) ()
Creates a deep copy of this object.

- virtual bool [saveSettings](#) (QDomDocument &aDoc, QDomElement &aElm) const
The method saves the settings of this object to the DOM element.
- virtual bool [loadSettings](#) (const QDomElement &aDom)
Loads settings from DOM.
- virtual void [run](#) ()
This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).
- virtual bool [currentOutputImageInfo](#) (SxPicBufInfo &picInfo, [CxImageMetadata](#) *pMetadata=NULL)
Gets the current output image info (format) provided by this object.
- [CxCameraParameters](#) * [parameters](#) ()
- XICORE_HANDLE [cameraHandle](#) () const
- void [invalidateCameraHandle](#) ()
- double [cameraSourceFrameRate](#) ()
The current receiving frame rate of the camera source (it is not the same as camera frame rate!).

Protected Member Functions

- virtual int [recommendedAdditionalBuffersInMemoryPool](#) () const
Recommended additional number of buffers that should be present in memory pool, but not reserved for this chainable only. Used in CameraSource to add more then minimum number of frames to pool.
- virtual void [setupStart](#) ()
Sets up the object before the object is started.
- virtual void [setupStop](#) ()
Sets up the object after the object is stopped.

Private Slots

- void [onImageFormatChanged](#) (XICORE_HANDLE hCamera)
Connected to camera manager.

Private Member Functions

- [CxCameraSource](#) (XICORE_HANDLE hCamera)
The constructor is private. Call the [CxChainManager::createCameraSource\(\)](#) for new instance of the class.
- int [bufferPolicy](#) ()
Gets the current buffer policy.
- bool [createFakeTransportPicture](#) (SxPicBuf *pPic, const [SxPicBufInfo](#) *pPicFormat, void *pSrcXilmg)
- void [checkCameraSkippingFrames](#) (QTime &timer, int &iApiCounter, int &iTransportCounter)
- bool [isUnpackingInMyChain](#) () const

Private Attributes

- XICORE_HANDLE [m_hCamera](#)
- qint32 [m_iFrameNo](#)
- QQueue< int > [m_queLastFrameTimes](#)
- QTime [m_timerLastFrames](#)
- [SxPicBufInfo](#) [m_aCurrentOutputFormat](#)
- bool [m_bCurrentOutputFormatValid](#)
- bool [m_bFakeTransportFormat](#)
- int [m_iPayloadSize](#)
The payload size of the current image.

Friends

- class **CxChainManager**

Additional Inherited Members

4.7.1 Detailed Description

The implementation of image source for Ximea cameras.

4.7.2 Member Function Documentation

4.7.2.1 virtual bool CxCameraSource::currentOutputImageInfo (SxPicBufInfo & *picInfo*, CxImageMetadata * *pMetadata* = NULL) [virtual]

Gets the current output image info (format) provided by this object.

Remarks

The method should call `getCurrentOutputImageInfo()` on its predecessors and calculate the output image info provided by this object.

Reimplemented from [CxImageProvider](#).

4.7.2.2 virtual void CxCameraSource::run () [virtual]

This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).

See also

[CxChainable](#)

Reimplemented from [CxChainable](#).

4.7.2.3 virtual bool CxCameraSource::saveSettings (QDomDocument & *aDoc*, QDomElement & *aElm*) const [virtual]

The method saves the settings of this object to the DOM element.

Remarks

The element `aElm` is already named by a caller, do not change the name (ie. do not call the `aElm.setTag←Name()`). You can add attributes and child elements.

Reimplemented from [CxChainable](#).

4.7.2.4 virtual void CxCameraSource::setupStart () [protected],[virtual]

Sets up the object before the object is started.

Remarks

The method is internally used by the [start\(\)](#) method immediately before the [run\(\)](#) method is invoked. For example, the [CxCameraSource](#) sets up internal frame counter to zero and starts camera acquisition. The default implementation does nothing.

Reimplemented from [CxChainable](#).

4.7.2.5 virtual void CxCameraSource::setupStop () [protected],[virtual]

Sets up the object after the object is stopped.

Remarks

The method is internally used by the [stop\(\)](#) method immediately after the [run\(\)](#) method is finished. For example, the [CxCameraSource](#) stops camera acquisition here. The default implementation does nothing.

Reimplemented from [CxChainable](#).

4.7.3 Member Data Documentation**4.7.3.1 SxPicBufInfo CxCameraSource::m_aCurrentOutputFormat [private]**

Cached current output format.

4.7.3.2 bool CxCameraSource::m_bCurrentOutputFormatValid [private]

Whether the cached output format is valid.

4.7.3.3 bool CxCameraSource::m_bFakeTransportFormat [private]

Camera is in Transport data format, generate and send fake images

4.7.3.4 XICORE_HANDLE CxCameraSource::m_hCamera [private]

Handle to camera.

4.7.3.5 qint32 CxCameraSource::m_iFrameNo [private]

Internal frame counter.

4.7.3.6 QQueue<int> CxCameraSource::m_queLastFrameTimes [private]

Acquisition times of the last frames.

4.7.3.7 QTime CxCameraSource::m_timerLastFrames [private]

Measuring of the frame's acquisition times.


The documentation for this class was generated from the following file:

- CameraSource.h

4.8 CxChain Class Reference

This class implements a processing chain which is together with the [CxChainable](#) class (chainable objects) the keystones of the XIMEA CamTool.

Inheritance diagram for CxChain:



Public Types

- enum **ExChainState** { **eStateStopped**, **eStateStopping**, **eStateStarted**, **eStateStarting** }
- enum **ExChainErrorState** { **eErrStateNoErr**, **eErrStateWarning**, **eErrStateError** }

Signals

- void [chainStarted](#) ([CxChain](#) *pChain)
The signal is emitted when the chain was started.
- void [chainStopped](#) ([CxChain](#) *pChain)
The signal is emitted when the chain was stopped.
- void [aboutToDelete](#) ([CxChain](#) *pSender)
The signal is emitted in destructor of the object.
- void [aboutToRemove](#) ([CxChain](#) *pSender, [CxChainable](#) *pRemoving)
The signal is emitted from [CxChain::remove\(\)](#) just before the chainable object is removed from the chain.
- void [chainErrorState](#) ([CxChain::ExChainErrorState](#) eState, [QString](#) sErrMsg)
The signal is emitted from [CxChain::start\(\)](#) if the chain start failed or some problem appeared while starting.

Public Member Functions

- [~CxChain](#) ()
Destructor.
- bool [start](#) ()
Starts processing on the chain.
- void [stop](#) ()
Stops the chain. Each chainable stops its task immediately, all unprocessed data will be lost. The method leaves its execution when all chainables are stopped (it's a blocking operation).
- bool [close](#) ()
Closes (but not deletes) the chain. It deletes all chainables and views in the chain.
- bool [started](#) ()
Returns true if the chain is running.
- bool [stopped](#) ()
Returns true if the chain is stopped.
- [ExChainState](#) [state](#) ()
Returns a state of the chain.
- [QString](#) [fullChainTitle](#) ()
Get the chain's full title.
- [QString](#) [chainTitle](#) ([CxChainable](#) *pChainable)

- Gets the chain's title for the given chainable.*
- QSet< CxChainable * > [allChainables](#) ()
Gets a set of all chainables in the chain. It does not matter if the chain contains branches.
- QSet< CxImageSource * > [imageSources](#) ()
Gets a set of all image sources in the chain.
- QSet< CxChainable * > [consumers](#) ()
Gets a set of all data consumers in the chain.
- QSet< QObject * > [nonChainableListeners](#) ()
Gets a set of all objects listening this chain.
- QSet< CxViewImageReceiver * > [viewReceivers](#) ()
Gets all view receivers in the chain.
- QSet< XICORE_HANDLE > [cameras](#) ()
Gets all cameras associated with this chain. Internally it searches all image sources which are the [CxCameraSource](#).

Private Member Functions

- [CxChain](#) ()
Constructor.
- bool [insertBefore](#) (CxChainable *pChainable, CxChainable *pBefore)
Inserts new element before the specified element.
- bool [branchAfter](#) (CxChainable *pChainable, CxChainable *pAfter)
Inserts new element after the specified element as a new branch.
- bool [insertAfter](#) (CxChainable *pChainable, CxChainable *pAfter)
Insert new element after the specified element.
- bool [remove](#) (CxChainable *pChainable)
The method removes given pChainable from this chain.
- bool [connectNonChainable](#) (CxChainable *pChainable, QObject *pListener, const char *onDataAvailable↔ Slot, const char *onFreePoolBuffersSlot)
Connects an arbitrary object pListener to the pChainable.
- bool [disconnectNonChainable](#) (CxChainable *pChainable, QObject *pListener, const char *onData↔ AvailableSlot, const char *onFreePoolBuffersSlot)
Disconnects an arbitrary object pListener from the pChainable which was connected using the connect↔ Chainable() method.
- CxChain * [duplicateChainWithViews](#) (CxChainable *pRoot, CxChainable *pNewRoot, bool bIncludeRoot)
Duplicates the chain starting from pRoot. All views (CxMdiView) will be duplicated too.
- CxChain * [duplicateBranch](#) (CxChainable *pRoot, CxChainable *pLast, CxChainable *pNewRoot, bool b↔ IncludeRoot)
Duplicates the branch of the starting from pRoot to pLast. If the pLast is an instance of CxViewImageReceiver then the eventual image view attached to the receiver will be duplicated as well.
- bool [addClonesWithViews](#) (CxChain *pChain, CxChainable *pChainable, CxChainable *pPrevClone, bool bAddChainableClone, bool bInsert)
Internal method used by duplicateChainWithViews().
- CxMdiView * [receiverView](#) (CxViewImageReceiver *pChainable)
Returns a view connected the the chainable.
- void [deleteAll](#) (CxChainable *pChainable)
Deletes the pChainable and all its successors.
- QSet< CxChainable * > [leaves](#) (CxChainable *pRoot=NULL)
Returns a set leaves (chainables a the end of each branche).
- void [removeBranch](#) (CxChainable *pLeaf)
Removes a branch given by its leaf.
- bool [saveSettings](#) (QDomDocument &aXml, QDomElement &aRoot)

- Exports the chain to DOM document.*

 - bool **loadSettings** (const QDomElement &aDom, QString *pErrorMsg=NULL)

Loads the chain from DOM.
 - bool **saveChainableSettings** (const CxChainable *pChainable, QDomDocument &aXml, QDomElement &aDom)

Returns UniqueID from metadata. If it is not present then it returns class name.
 - CxChainable * **loadChainable** (const QDomElement &aElm, CxChainable *pAfter)

Recursive method that creates and loads a chainable and all successors from DOM.
 - QSet< CxChainable * > **disconnectSuccessors** (CxChainable *pChainable)

All chainables in the disconnected branches will be removed from chain (their chain will be set to NULL). For reconnecting use [connectBranch\(\)](#)
 - void **connectBranch** (CxChainable *pBranch, CxChainable *pRoot)

Connect the pBranch to the pRoot. It is opposite method to [disconnectSuccessors\(\)](#)
 - void **deleteSuccessors** (CxChainable *pChainable)
 - void **deleteViewReceiver** (CxViewImageReceiver *pReceiver)
 - void **viewReceivers** (CxChainable *pChainable, QSet< CxViewImageReceiver * > &setReceivers)
- Recursive method that fills the setReceivers by all view receivers after the pChainable.*
- CxCameraSource * **getExistingNotOpenedCamera** (const QString &sCamSerial, const QString &sCamName, const QString &sCamType, QString *pErrorMsg)

Static Private Member Functions

- static bool **isImageSource** (const QString &sClass)

Private Attributes

- QMutex **m_lock**
- CxChainable * **m_pRoot**
- ExChainState **m_eState**

Friends

- class **CxChainManager**

4.8.1 Detailed Description

This class implements a processing chain which is together with the [CxChainable](#) class (chainable objects) the keystones of the XIMEA CamTool.

4.8.2 Chains and chainable objects

Processing chain (a pipeline) defines a data flow and processing of the data.

Note

Currently the chains are implemented as a tree with single root (data source).

In the XIMEA CamTool the chain consists of these chainables:

- Data source - currently implemented by [CxImageSource](#).
- Processing elements - mostly implemented by [CxImageProvider](#).

- Data consumers - a chainables where the chain ends i.e. consumers are leafs in the chain tree. The xiCore implements these consumers: [CxViewImageReceiver](#), [CxCircularBuffer](#), [CxVideoToFile](#).

The data flow and interaction between chainables is shown on the image below. The chain on the image consists of:

- CameraSource ([CxCameraSource](#)) which is the implementation of the [CxImageSource](#).
- RAW coloring ([CxBayerPseudoColorsCnbl](#)) - the implementation of the [CxImageProvider](#).
- VideoToFile ([CxVideoToFile](#)) - the implementation of [CxChainable](#) as a data consumer. The [VideoToFile](#) has pointer to ImageSaver ([IxImageDataSaver](#)).

The data that flows through the chain must inherit the [IxChainData](#) interface. Currently the [CxImageData](#) and [CxViewImageData](#) classes based on that interface are implemented.

Any chain can consist of an arbitrary number of branches. Each branch gets its own copy of the data coming from data source. The copy of data is created automatically in the [CxChainable::data\(\)](#) method.

The [CxChain](#) class implements several methods for adding/removing the [CxChainable](#) to/from the chain. These methods are private. Chains are managed by [CxChainManager](#) which is a friend class of the [CxChain](#).

4.8.2.1 Connecting non-chainable object to the chain

The current implementation of chains allows to connect an arbitrary object to any chainable element. This object then receives its own copy of chain data. The arbitrary object can be connected to the chainable when it satisfies these conditions:

- It must be a descendant of the [QObject](#).
- It must implement a slot to be able to listen the [CxChainable::dataAvailable\(\)](#) signal.
- It must implement a slot to be able to listen the [CxMemoryManager::freePoolBuffers\(\)](#) signal.

The object can be connected to the chainable in the chain by [CxChainManager::connectNonChainable\(\)](#) method. When the object is connected to the chain then it behaves similar to [CxChainable](#). It means that it receives the [CxChainable::dataAvailable\(\)](#) signal and can call [CxChainable::data\(\)](#) or [CxChainable::discardData\(\)](#) etc.

The mechanism of connected non-chainable objects is used by views ([CxMdiView](#)) in the GUI frontend. So that the chain architecture allows to connect any number of views to chainables in the chain and it is guaranteed that each of the connected objects (views) receives its own copy of the data.

4.8.3 Member Function Documentation

4.8.3.1 `bool CxChain::branchAfter ( CxChainable * pChainable, CxChainable * pAfter ) [private]`

Inserts new element after the specified element as a new branch.

Remarks

If the `pAfter` has successors then the `pChainable` will be added as a new branch. To be able to insert after the element, these conditions must be satisfied:

- both of the parameters must be at the state of [CxChainable::eStopped](#),
- the current chain of the `pChainable` must be `NULL` (`chain()`),
- the current chain of the `pAfter` must be this chain,
- the `pAfter` must not be a data consumer (the [CxChainable::isDataConsumer\(\)](#) must return true,
- if this chain is empty then the `pAfter` must be `NULL`.
- if the `pAfter` is not `NULL` then the `pChainable` must not be [CxImageSource](#).

4.8.3.2 `bool CxChain::connectNonChainable ( CxChainable * pChainable, QObject * pListener, const char * onDataAvailableSlot, const char * onFreePoolBuffersSlot ) [private]`

Connects an arbitrary object `pListener` to the `pChainable`.

Remarks

the `onDataAvailableSlot` and `onFreePoolBuffersSlot` must be generated by the SLOT macro.

See also

[disconnectNonChainable\(\)](#)

4.8.3.3 `bool CxChain::disconnectNonChainable ( CxChainable * pChainable, QObject * pListener, const char * onDataAvailableSlot, const char * onFreePoolBuffersSlot ) [private]`

Disconnects an arbitrary object `pListener` from the `pChainable` which was connected using the `connectNonChainable()` method.

Remarks

the `onDataAvailableSlot` and `onFreePoolBuffersSlot` must be generated by the SLOT macro.

4.8.3.4 `CxChain* CxChain::duplicateBranch ( CxChainable * pRoot, CxChainable * pLast, CxChainable * pNewRoot, bool bIncludeRoot ) [private]`

Duplicates the branch of the starting from `pRoot` to `pLast`. If the `pLast` is an instance of [CxViewImageReceiver](#) then the eventual image view attached to the receiver will be duplicated as well.

Remarks

All the duplicated members of the chain will be appended to the `pNewRoot`. If the `bIncludeRoot` is true, then the duplicated `pRoot` will be included in the duplicated chain as well.

4.8.3.5 `CxChain* CxChain::duplicateChainWithViews ( CxChainable * pRoot, CxChainable * pNewRoot, bool bIncludeRoot ) [private]`

Duplicates the chain starting from `pRoot`. All views ([CxMdiView](#)) will be duplicated too.

Remarks

All the duplicated members of the chain will be appended to the `pNewRoot`. If the `bIncludeRoot` is true, then the duplicated `pRoot` will be included in the duplicated chain as well.

4.8.3.6 `bool CxChain::insertAfter ( CxChainable * pChainable, CxChainable * pAfter ) [private]`

Insert new element after the specified element.

Remarks

The current successors of the `pAfter` will be reconnected to the `pChainable` and the `pChainable` becomes the new successor of the `pAfter`.

4.8.3.7 `bool CxChain::insertBefore ( CxChainable * pChainable, CxChainable * pBefore ) [private]`

Inserts new element before the specified element.

Remarks

If the `pBefore` has a predecessor, then it will be reconnected so that it became the predecessor of the `pChainable`. To be able to insert before the element, these conditions must be satisfied:

- both of the parameters must be at the state of `CxChainable::eStopped`,
- the current chain of the `pChainable` must be `NULL` (`chain()`),
- the current chain of the `pBefore` must be this chain,
- the `pBefore` must not be a [CxImageSource](#),
- the `pChainable` must not be a data consumer (the [CxChainable::isDataConsumer\(\)](#) must return false.
- if this chain is empty then the `pBefore` must be `NULL`.

4.8.3.8 `bool CxChain::loadSettings ( const QDomElement & aDom, QString * pErrorMsg = NULL ) [private]`

Loads the chain from DOM.

Remarks

If the chain is empty all the chain from DOM will be imported. If the current chain is not empty (it contains a image source) then the image source will be kept and the rest of the chain will be replaced by the chain from DOM and the image source from DOM will be excluded from import.

4.8.3.9 `bool CxChain::saveChainableSettings ( const CxChainable * pChainable, QDomDocument & aXml, QDomElement & aDom ) [private]`

Returns UniqueID from metadata. If it is not present then it returns class name.

Recursive method that saves a single chainable settings and all its successors.

4.8.3.10 `bool CxChain::start ( )`

Starts processing on the chain.

Returns

The method return false if the chain is locked or already started.

4.8.4 Member Data Documentation

4.8.4.1 `ExChainState CxChain::m_eState [private]`

State of this chain.

4.8.4.2 `QMutex CxChain::m_lock [private]`

Lock of this object

4.8.4.3 CxChainable* CxChain::m_pRoot [private]

The root of the chain


The documentation for this class was generated from the following file:

- Chain.h

4.9 CxChainable Class Reference

A base class of all chainable objects.

Inheritance diagram for CxChainable:



Public Types

- enum [ExChainableState](#) { **eUndefined** = 0, **eStopped**, **eRunning** }
Value representing a state of the [CxChainable](#) object.
- enum [ExDataSourceType](#) { **eSourceTypeNone** = 0, **eSourceTypeCamera**, **eSourceTypeStaticImage**, **eSourceTypeVideo** }
Values representing an data source type.
- enum [ExChainableFlag](#) {
eHasOwnMemoryPool = 0x1, **eRunsInfiniteLoop** = 0x2, **eSupportsFpsRegulation** = 0x4, **eCanChangeEnableState** = 0x8,
eTakeLatestImage = 0x10 }
Flags of the chainable object.

Signals

- void [dataAvailable](#) ([CxChainable](#) *pSender, int iFrameNo)
The signal is emitted when this objects has data available for next chainable objects.
- void [errorStateChanged](#) ([CxChainable](#) *pSender, bool bError)
The signal is emitted when some error appears or when an error is cleared. It is emitted by method [setErrorMessag](#)();.
- void [aboutToDelete](#) ([CxChainable](#) *pSender)
The signal is emitted from the destructor of this object.
- void [startRun](#) ()
The signal is internally connected to the [run\(\)](#) slot.

Public Member Functions

- **CxChainable** (ExChainableFlags eFlags)
- virtual QString **title** () const =0
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual bool **acceptsDataFrom** (CxChainable *pPredecessor)=0
Query if this object can accept data from the given chainable object.
- virtual int **bufferCountInMemoryPool** () const =0
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*
- virtual CxChainable * **clone** ()=0
Creates a deep copy of this object.
- virtual bool **isDataConsumer** () const
Returns true if the object is data consumer, i.e. the object has no successors and all the received data are "consumed" here e.g. they are displayed in image view, saved to disk etc. In other words this object is the a leaf in the chain.
- virtual ExChainableState **state** ()
- virtual IxChainData * **data** (QObject *pReceiver, qint32 iFrameNo)
Gets the pointer to data processed by this object.
- virtual void **discardData** (QObject *pReceiver, qint32 iFrameNo)
This method discards data for the given receiver.
- virtual IxChainData * **latestData** (QObject *pReceiver)
Gets the latest data from the output queue. It is called automatically if ExChainableFlag::eTakeLatestImage is set.
- virtual void **discardLatestData** (QObject *pReceiver)
Discards the latest data for the given receiver. It is called automatically if ExChainableFlag::eTakeLatestImage is set.
- virtual QWidget * **configurationWidget** ()
Displays the object's configuration dialog.
- virtual bool **hasConfigurationWidget** () const
Returns true when the the object has configuration widget (ie. the [configurationWidget\(\)](#) returns not NULL).
- virtual bool **saveSettings** (QDomDocument &aDoc, QDomElement &aElm) const
The method saves the settings of this object to the DOM element.
- virtual bool **loadSettings** (const QDomElement &aDom)
Loads settings from DOM.
- virtual bool **hasHelp** () const
Returns true if the object has help. The default is false.
- virtual QString **help** (QString &slmagesPath)
Returns a string with help.
- virtual void **setMaxFps** (double dFps)
Sets max. FPS (when zero or negative, then the FPS will not be regulated)
- virtual double **maxFps** () const
Returns max. FPS (zero means that the max FPS is not set).
- const QSet< CxChainable * > & **precedessors** ()
Gets the set of precedessors of this objects in the chain.
- const QSet< CxChainable * > & **successors** () const
Gets the set of successors of this objects in the chain.
- CxChain * **chain** () const
Gets a chain in which this object is contained.
- bool **supportsFpsRegulation** ()
Whether this object supports FPS regulation.
- QSet< QObject * > **nonChainableListeners** ()
Returns a set of all non-CxChainable objects connected to this object.
- TxMemPoolHandle **memoryPool** ()

- Returns a handle to the object's memory pool.*
- bool **enabled** () const
Whether the object is enabled.
- void **setEnabled** (bool bEnabled)
Can be applied when [CxChainable::eCanChangeEnableState](#) flag is set. The chain of this object must be stopped.
- ExChainableFlags **flags** () const
- QList< [ExDataSourceType](#) > **chainDataSourceTypes** () const
Gets the type of data sources of this chainable object.
- QList< [CxChainable](#) * > **dataSources** ()
Gets list of data sources of this chainable object.
- QString **errorMessage** ()
Error message of this object.

Protected Slots

- virtual void **onFreePoolBuffers** (TxMemPoolHandle hMemoryPool)
Releases all buffers allocated from the given memory pool.
- virtual void **onDataAvailable** ([CxChainable](#) *pSender, int iFrameNo)
Reaction on the [dataAvailable\(\)](#) emitted by predecessors of this object. The slot is connected automatically by the [addSuccessor\(\)](#).
- virtual void **run** ()
This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).
- virtual void **onChainStarted** ([CxChain](#) *pChain)
Handle the [CxChain::chainStarted\(\)](#). The default implementation does nothing.
- virtual void **onChainStopped** ([CxChain](#) *pChain)
Handle the [CxChain::chainStarted\(\)](#). The default implementation does nothing.

Protected Member Functions

- virtual bool **init** (QString *pErrMsg=NULL)=0
Initializes the object. When calling this method, the object MUST be stopped.
- virtual bool **initialized** ()=0
- virtual void **setupStart** ()
Sets up the object before the object is started.
- virtual void **setupStop** ()
Sets up the object after the object is stopped.
- bool **isMemoryPoolNeeded** ()
Returns true if this object requires memory pool. The default implementation just tests the [m_eFlags](#), but for example the camera source should take into account the actual buffer policy of the camera.
- virtual int **recommendedAdditionalBuffersInMemoryPool** () const
Recommended additional number of buffers that should be present in memory pool, but not reserved for this chainable only. Used in CameraSource to add more then minimum number of frames to pool.
- virtual [IxChainData](#) * **processData** ([IxChainData](#) *pReceivedData)
The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.
- virtual bool **addSuccessor** ([CxChainable](#) *pSuccessor)
Adds new successor to this object and connects the [dataAvailable\(\)](#) signal to the successor's [onDataAvailable\(\)](#) slot. It is internally used by [CxChain](#) class.
- virtual bool **addPredecessor** ([CxChainable](#) *pPredecessor)

- Adds new predecessor to this object. It is internally used by [CxChain](#) class.*

 - virtual bool [removeSuccessor](#) (CxChainable *pSuccessor)

Removes the successor from this object and discards the connection of the [dataAvailable\(\)](#) signal from the successor's [onDataAvailable\(\)](#) slot. It is internally used by [CxChain](#) class.
 - virtual bool [removePredecessor](#) (CxChainable *pPredecessor)

Removes the predecessor from this object. It is internally used by [CxChain](#) class.
 - virtual bool [clearSuccessors](#) ()

Removes and disconnect all the successors of this object. It internally calls the [removeSuccessor\(\)](#).
 - virtual bool [clearPredecessors](#) ()

Removes all the predecessors of this object. It internally calls the [removePredecessor\(\)](#).
 - virtual bool [addToOutputData](#) (IxChainData *pData)

Adds the given data to output data and emits the [dataAvailable\(\)](#) signal if necessary.
 - void [deleteOldestOutputDataNoLock](#) ()

Deletes the output data. The method is not thread safe, it is necessary to lock the [m_lockOutputData](#)!
 - void [setErrorMessag](#) (const QString &errMsg)

Set the error message and emits [errorStateChanged\(\)](#) signal.
 - bool [start](#) (CxChain::ExChainErrorState *pErrState=NULL, QString *pErrMsg=NULL)

Starts this object and sets the correct state of the object.
 - bool [stop](#) ()

Stops running and resets this object (release appropriate buffers etc.). The method stops receiving input data and processing incoming data.
 - void [waitForStarted](#) ()

The method waits until the [run\(\)](#) method is started.
 - void [waitForStopped](#) ()

The method waits until the [run\(\)](#) method is finished.
 - void [freeOutputData](#) ()

The method frees all output data [CxChainable::m_mapOutputData](#).
 - void [connect](#) (QObject *pReceiver, const char *method)

This method registers an arbitrary non-chainable receiver to receive a data provided by this object.
 - void [disconnect](#) (QObject *pReceiver, const char *method)

*Discards the connection between the [pReceiver](#) and this object created by [connect\(const QObject *pReceiver, const char *method\)](#) method.*
 - int [objectBuffersInMemoryPool](#) ()

*Returns the number of buffers required by **this** object to have in memory pool.*

Protected Attributes

- QSet< CxChainable * > [m_setPredecessors](#)
- QSet< CxChainable * > [m_setSuccessors](#)
- int [m_iMaxOutputDataCount](#)
- ExChainableFlags [m_eFlags](#)
- double [m_dMaxFps](#)
- QQueue< int > [m_queFpsCounter](#)
- QTime [m_timerFps](#)
- ExChainableState [m_eState](#)
- ExChainableState [m_eReqState](#)
- QMutex [m_lock](#)
- QMutex [m_lockOutputData](#)
- QMap< int, IxChainData * > [m_mapOutputData](#)
- QMap< QObject *, QSet< int > > [m_mapObjectOutputData](#)
- XICORE_HANDLE [m_hMemoryPool](#)
- QSemaphore [m_semStartStop](#)
- bool [m_bDataConsumer](#)

Private Member Functions

- bool [setChain](#) ([CxChain](#) *pChain, bool bSetSuccessors=false)
Sets the 'parent' chain to this object.
- double [calcFps](#) ()
Calculates current FPS (if FPS regulation is not supported the it returns 0)
- bool [addToOutputDataNoLock](#) ([IxChainData](#) *pData)
- [IxChainData](#) * [dataNoLock](#) (QObject *pReceiver, qint32 iFrameNo)
- void [discardDataNoLock](#) (QObject *pReceiver, qint32 iFrameNo)
- QThread * [myThread](#) () const
Gets the internal thread of this object.

Private Attributes

- [CxChain](#) * [m_pMyChain](#)
- QThread * [m_pThread](#)
- QString [m_sErrorMessage](#)
- bool [m_bLoosingFrames](#)
- int [m_iFirstLostFrame](#)
- bool [m_bEnabled](#)

Friends

- class [CxChain](#)
- class [CxChainManager](#)

4.9.1 Detailed Description

A base class of all chainable objects.

Remarks

The [CxChainable](#) class together with the [CxChain](#) is the key stone of the xiCore architecture. It is recommended to read the [CxChain](#) documentation before you start to study this class.

4.9.2 Implementation of chainable object

The [CxChainable](#) class is an abstract class that implements:

- starting and stopping the object ([start\(\)](#), [stop\(\)](#));
- receiving a data from its predecessors in the chain ([onDataAvailable\(\)](#), [data\(\)](#), [discardData\(\)](#));
- notifying its successors when it has a processed data ready ([dataAvailable\(\)](#)).

This class generally can process any kind of data ([IxChainData](#)). When implementing a custom [CxChainable](#) class then in most cases the [CxImageProvider](#) is the class you should inherit from. Currently the [CxImageProvider](#) is the only class that can be registered by the application and thus that can be added to a processing chain by the chain editor (see plugins).

Memory management

Any chainable object can use its own memory pool ([CxMemoryManager](#)). This is specified by the [CxChainable::eHasOwnMemoryPool](#) flag passed to the constructor. The memory pool should be initialized in the pure virtual method [init\(\)](#). This method is implemented by the [CxImageProvider](#) so that when you

inherit from it you do not need to implement this method as well as the `initialized()` method.

Each chainable object in the chain has its `onFreePoolBuffers()` slot connected to the `CxMemoryManager::freePoolBuffers()`.

When the signal is emitted then the chainable object **always must** immediately release all buffers allocated from the given memory pool otherwise the application crashes on assertion!

Data ownership

There are simple rules describing the data ownership:

1. If the object created (allocated) a data (`lxChainData`) then it is the owner of the data until some of its successors take the data by calling `data()` method.
2. If the object takes the data by calling `data()` method on its predecessor then it become the owner of the data until some of its successors take the data by calling `data()` method.

When the object is the owner of the data then there are two scenarios how the data can be managed:

1. The chainable object processes the data directly. It means that the format and the size of the data is unchanged and the data are processed in-place. After the processing the data are passed to the object's successors. In this case the chainable object cannot delete the data until the successors pick up the data by calling the `data()` method.
2. The chainable object changes the data size, format etc. so that it needs to create (alloc) new data which are passed to the object's successors. In this case the object must delete the original (received) data after the processing finished.

Sometimes it can happen that the chainable object cannot accept a data from its predecessor. It means that its predecessor notified the object by `dataAvailable()` signal but the object cannot process them i.e. it cannot call the `data()` method. In this case the object should call the `discardData()` on the predecessor instead of the `data()`.

4.9.3 Member Enumeration Documentation

4.9.3.1 enum CxChainable::ExChainableFlag

Flags of the chainable object.

Enumerator

eHasOwnMemoryPool Object uses own memory pool.

eRunsInfiniteLoop The `run()` method is implemented as an infinite loop.

eSupportsFpsRegulation The object allows to regulate a framerate by discarding input data.

eCanChangeEnableState The object can be disabled so that it will not process a data and the data flow is bypassed

eTakeLatestImage The object always gets the latest image from its predecessor's output queue.

4.9.3.2 enum CxChainable::ExDataSourceType

Values representing an data source type.

Enumerator

eSourceTypeNone The object is not a data source.

eSourceTypeCamera Data come from camera (`CxCameraSource`).

eSourceTypeStaticImage Data come from static image (`CxStaticImageSource`).

eSourceTypeVideo Data come from video (`CxVideoSource`).

4.9.4 Member Function Documentation

4.9.4.1 `QList<ExDataSourceType> CxChainable::chainDataSourceTypes ( ) const`

Gets the type of data sources of this chainable object.

Remarks

The current implementation of the [CxChain](#) supports only one data source so that the method return a list containing at most one element.

See also

[dataSources\(\)](#)

4.9.4.2 `virtual QWidget* CxChainable::configurationWidget ( ) [virtual]`

Displays the object's configuration dialog.

Remarks

The configuration widget should not contain the Ok and Cancel button as they will be added automatically by the application.

Warning

When this chainable object changes the image format, don't forget to stop the current chain before applying new value from the dialog!
This method MUST be called from GUI (main) thread only!!

Returns

The default value is NULL. When overridden then it should return the pointer to the configuration widget.

Reimplemented in [CxViewImageReceiver](#).

4.9.4.3 `void CxChainable::connect ( QObject * pReceiver, const char * method ) [protected]`

This method registers an arbitrary non-chainable receiver to receive a data provided by this object.

Remarks

The method creates connection between its [dataAvailable\(\)](#) signal and the receiver's slot given by `method` parameter.

Warning

Always use the SLOT macro to specify the method name given by `method` parameter.
Each object which wants to receive the data must be connected using this method otherwise the behaviour is undefined.

See also

`disconnect(const QObject *pReceiver, const char *method)`

4.9.4.4 `virtual IxChainData* CxChainable::data ( QObject * pReceiver, qint32 iFrameNo ) [virtual]`

Gets the pointer to data processed by this object.

See also

[discardData\(\)](#)
[CxChainable](#)

Reimplemented in [CxStaticImageSource](#).

4.9.4.5 `void CxChainable::dataAvailable ( CxChainable * pSender, int iFrameNo ) [signal]`

The signal is emitted when this objects has data available for next chainable objects.

See also

[CxChainable](#)

4.9.4.6 `QList<CxChainable*> CxChainable::dataSources ( )`

Gets list of data sources of this chainable object.

Remarks

The current implementation of the [CxChain](#) supports only one data source so that the method return a list containing at most one element.

See also

[chainDataSourceTypes\(\)](#)

4.9.4.7 `virtual void CxChainable::discardData ( QObject * pReceiver, qint32 iFrameNo ) [virtual]`

This method discards data for the given receiver.

See also

[data\(\)](#)
[CxChainable](#)

4.9.4.8 `void CxChainable::disconnect ( QObject * pReceiver, const char * method ) [protected]`

Discards the connection between the `pReceiver` and this object created by `connect(const QObject *pReceiver, const char *method)` method.

Warning

Allways use the SLOT macro to specify the method name given by the `method` parameter.

See also

`connect(const QObject *pReceiver, const char *method)`

4.9.4.9 bool CxChainable::enabled () const

Whether the object is enabled.

See also

[CxChainable::eCanChangeEnableState](#)
[setEnabled\(\)](#)

4.9.4.10 QString CxChainable::errorMessage ()

Error message of this object.

See also

[setErrorMessage\(\)](#)

4.9.4.11 virtual QString CxChainable::help (QString & sImagesPath) [virtual]

Returns a string with help.

Remarks

The text can be plain text or HTML. If the HTML help contains images then sImagesPath should contain an path to the images. The images path must be either absolute or relative to the application path (qApp->applicationDirPath()) and it also can be a prefix to embedded resource images (for example ":/images").

Reimplemented in [CxDemosaic](#).

4.9.4.12 virtual bool CxChainable::init (QString * pErrMsg = NULL) [protected],[pure virtual]

Initializes the object. When calling this method, the object MUST be stopped.

Remarks

Methods [init\(\)](#) and [initialized\(\)](#) are internally used by the [start\(\)](#) method. It should initialize the chainable object so that it is ready to start ([start\(\)](#)) . If this object is the [CxImageProvider](#) then the method initializes internal memory pool.

Implemented in [CxImageProvider](#), [CxVideoToFile](#), and [CxCircularBuffer](#).

4.9.4.13 TxMemPoolHandle CxChainable::memoryPool ()

Returns a handle to the object's memory pool.

Remarks

If the object has no internal memory pool then the method returns NULL.

4.9.4.14 int CxChainable::objectBuffersInMemoryPool () [protected]

Returns the number of buffers required by **this** object to have in memory pool.

Remarks

The method calls internally the [buffersCountInMemoryPool\(\)](#) and add a number of listeners of this object.

4.9.4.15 `virtual void CxChainable::onFreePoolBuffers ( TxMemPoolHandle hMemoryPool ) [protected], [virtual], [slot]`

Releases all buffers allocated from the given memory pool.

Remarks

This slot is always connected to [CxMemoryManager::freePoolBuffers\(\)](#) signal. The chainable object **MU↔ST** release all buffers allocated from the given memory pool. If the buffers will not be released then the application's behaviour is undefined (the app. most probably crashes or freezes). The default implementation of the method checks the internal output data ([CxChainable::m_mapOutputData](#)) and releases all the data from the given memory pool.

4.9.4.16 `virtual IxChainData* CxChainable::processData ( IxChainData * pReceivedData ) [protected], [virtual]`

The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Remarks

If the implementation creates (allocates) a new data here then the original (received) data should be deleted using the delete operator. The default implementation of this method returns a pointer to the received data without any processing.

Warning

If the method returns NULL, then the *pReceivedData* will be deleted automatically!

Reimplemented in [CxViewImageReceiver](#), [CxDemosaic](#), and [CxRgb32Convertor](#).

4.9.4.17 `virtual void CxChainable::run ( ) [inline], [protected], [virtual], [slot]`

This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).

See also

[CxChainable](#)

Reimplemented in [CxCameraSource](#), [CxStaticImageSource](#), and [CxVideoSource](#).

4.9.4.18 `virtual bool CxChainable::saveSettings ( QDomDocument & aDoc, QDomElement & aElm ) const [virtual]`

The method saves the settings of this object to the DOM element.

Remarks

The element *aElm* is already named by a caller, do not change the name (ie. do not call the *aElm.setTag↔Name()*). You can add attributes and child elements.

Reimplemented in [CxCameraSource](#), [CxStaticImageSource](#), and [CxVideoSource](#).

4.9.4.19 `bool CxChainable::setChain ( CxChain * pChain, bool bSetSuccessors = false ) [private]`

Sets the 'parent' chain to this object.

Remarks

The method succeeds when this object is not assigned to other chain. If the `pChain` is NULL, then the method succeeds when the current object's chain does not have this chainable object connected i.e. the `CxChain::isConnectedInChain()` method returns false.

4.9.4.20 `virtual void CxChainable::setupStart ( ) [protected],[virtual]`

Sets up the object before the object is started.

Remarks

The method is internally used by the `start()` method immediately before the `run()` method is invoked. For example, the `CxCameraSource` sets up internal frame counter to zero and starts camera acquisition. The default implementation does nothing.

Reimplemented in `CxCameraSource`, and `CxDemosaic`.

4.9.4.21 `virtual void CxChainable::setupStop ( ) [protected],[virtual]`

Sets up the object after the object is stopped.

Remarks

The method is internally used by the `stop()` method immediately after the `run()` method is finished. For example, the `CxCameraSource` stops camera acquisition here. The default implementation does nothing.

Reimplemented in `CxCameraSource`, `CxDemosaic`, and `CxRgb32Convertor`.

4.9.4.22 `bool CxChainable::start ( CxChain::ExChainErrorState * pErrState = NULL, QString * pErrMsg = NULL ) [protected]`

Starts this object and sets the correct state of the object.

Remarks

When the object is started then it can process and accept data.

4.9.4.23 `void CxChainable::startRun ( ) [signal]`

The signal is internally connected to the `run()` slot.

Remarks

It is internally emitted in the `start()` method.

Warning

Do not emit or receive the signal!

4.9.4.24 `bool CxChainable::stop ( ) [protected]`

Stops running and resets this object (release appropriate buffers etc.). The method stops receiving input data and processing incoming data.

Remarks

When the chainable object is stopped then it must be guaranteed that:

- the object cannot emit `dataAvailable()`,
- when the `data()` is invoked, then the method must return NULL,
- when the `onDataAvailable()` slot is invoked then it cannot call the `data()` on the sender. I.e. when the object is stopped it cannot provide or receive any data. The default implementation of the `CxChainable` satisfies the conditions above.

4.9.4.25 `void CxChainable::waitForStarted ( ) [protected]`

The method waits until the `run()` method is started.

Remarks

The method is internally used by the `start()` method.

Warning

Do not use it in your code!

4.9.4.26 `void CxChainable::waitForStopped ( ) [protected]`

The method waits until the `run()` method is finished.

Remarks

The method is internally used by the `stop()` method.

Warning

Do not use it in your code!

4.9.5 Member Data Documentation

4.9.5.1 `bool CxChainable::m_bDataConsumer [protected]`

Whether the object is data consumer.

4.9.5.2 `bool CxChainable::m_bEnabled [private]`

Enable state is valid when `CxChainable::eCanChangeEnableState` flag is set

4.9.5.3 `bool CxChainable::m_bLoosingFrames [private]`

It is internally used to signaling the error state.

Warning

Do not touch!

4.9.5.4 double CxChainable::m_dMaxFps [protected]

Max. FPS processed by this object. (-1 for unlimited)

4.9.5.5 ExChainableFlags CxChainable::m_eFlags [protected]

Flags of this object.

4.9.5.6 ExChainableState CxChainable::m_eReqState [protected]

Required state of this object. It is used when `run()` is an infinite loop. Do not touch!

4.9.5.7 ExChainableState CxChainable::m_eState [protected]

State of this object.

4.9.5.8 XICORE_HANDLE CxChainable::m_hMemoryPool [protected]

Handle to memory pool.

4.9.5.9 int CxChainable::m_iFirstLostFrame [private]

It is internally used to signalizing the error state.

Warning

Do not touch!

4.9.5.10 int CxChainable::m_iMaxOutputDataCount [protected]

Max number of output data per successor.

4.9.5.11 QMutex CxChainable::m_lock [protected]

Lock of this object.

4.9.5.12 QMutex CxChainable::m_lockOutputData [protected]

Lock of the output data.

4.9.5.13 QMap<QObject*, QSet<int> > CxChainable::m_mapObjectOutputData [protected]

The map of the output frame numbers. Key ~ the receiver, value ~ set of frame numbers

4.9.5.14 QMap<int, IxChainData*> CxChainable::m_mapOutputData [protected]

The output data, key ~ the frame number, value ~ the processed data.

4.9.5.15 CxChain* CxChainable::m_pMyChain [private]

Chain which this object belongs to.

4.9.5.16 QThread* CxChainable::m_pThread [private]

The thread in which this object lives (if NULL, then the object lives in the GUI thread)

4.9.5.17 QQueue<int> CxChainable::m_queFpsCounter [protected]

Queue of times when data was received - for FPS calculation.

4.9.5.18 QSemaphore CxChainable::m_semStartStop [protected]

Semaphore signaling that the [run\(\)](#) started or stopped.

Warning

Do not touch!

4.9.5.19 QString CxChainable::m_sErrorMessage [private]

The error message of this object.

4.9.5.20 QSet<CxChainable*> CxChainable::m_setPredecessors [protected]

Set of predecessors of this object in the chain.

4.9.5.21 QSet<CxChainable*> CxChainable::m_setSuccessors [protected]

Set of successors of this object in the chain.

4.9.5.22 QTimer CxChainable::m_timerFps [protected]

Timer for FPS measurement


The documentation for this class was generated from the following file:

- [Chainable.h](#)

4.10 CxChainManager Class Reference

The class that creates, manages and manipulates with chains.

Inheritance diagram for CxChainManager:



Signals

- void **chainErrorState** (CxChain *pChain, CxChain::ExChainErrorState eState, QString sErrMsg)

Public Member Functions

- **~CxChainManager** ()
Destructor.
- CxChain * **createChain** ()
Creates an instance of new chain.
- CxCameraSource * **createCameraSource** (XICORE_HANDLE hCamera)
Creates an instance of camera source.
- CxCameraSource * **cameraSource** (XICORE_HANDLE hCamera)
Gets the pointer to the existing camera source.
- void **cameraDisconnected** (XICORE_HANDLE hCamera)
For internal use only. Keeps the existing cameraSource, but sets its handle to NULL.
- CxImageSource * **createImageSource** (IxImageDataLoader *pLoader)
Creates an instance of an image source.
- CxStaticImageSource * **createStaticImageSource** (CxImageData *pImageData, const QString &sTitle)
Creates an instance of static image source.
- bool **closeChain** (CxChain *pChain)
Closes the chain.
- bool **deleteChainable** (CxChainable *pChainable)
Removes and deletes the chainable from the manager. The chainable must not be connected to any chain.
- bool **deleteBranches** (CxChain *pChain, CxChainable *pRoot)
Deletes all branches (including views) of the given chainable pRoot.
- CxChain * **cameraChain** (XICORE_HANDLE hCamera)
Gets a camera's chain.
- bool **insertBefore** (CxChain *pChain, CxChainable *pChainable, CxChainable *pBefore)
Inserts the pChainable before the pBefore object in the chain given by pChain.
- bool **branchAfter** (CxChain *pChain, CxChainable *pChainable, CxChainable *pAfter)
Inserts the pChainable after the pAfter in the chain given by pChain as a new branch.
- bool **insertAfter** (CxChain *pChain, CxChainable *pChainable, CxChainable *pAfter)
Inserts the pChainable after the pAfter in the chain given by pChain.
- bool **remove** (CxChain *pChain, CxChainable *pChainable, bool bDeleteBranchWhenNoListener=false)
Removes the pChainable from the pChain.
- bool **connectNonChainable** (CxChain *pChain, CxChainable *pChainable, QObject *pListener, const char *onDataAvailableSlot, const char *onFreePoolBuffersSlot)
The method registers a custom object as a data receiver of the pChainable.
- bool **disconnectNonChainable** (CxChain *pChain, CxChainable *pChainable, QObject *pListener, const char *onDataAvailableSlot, const char *onFreePoolBuffersSlot, bool bDeleteWhenNoListener=false)
Discards the connection created by connectNonChainable()
- CxChain * **duplicateChainWithViews** (CxChain *pChain, CxChainable *pRoot, CxChainable *pNewRoot, bool bIncludeRoot)
Duplicates the pChain starting from pRoot. All views (CxMdiView) will be duplicated too.
- CxChain * **duplicateChainBranch** (CxChain *pChain, CxChainable *pRoot, CxChainable *pLast, CxChainable *pNewRoot, bool bIncludeRoot)
Duplicates the branch of the pChain starting from pRoot to pLast. If the pLast is an instance of CxView/ImageReceiver then the eventual image view attached to the receiver will be duplicated as well.
- bool **loadChainSettings** (const QString &sFileName, QDomDocument &aDocChain, QString *pErrMsg=NULL)
Loads the chain settings from the given file to a DOM document.

- bool [importToExistingChain](#) ([CxChain](#) *pChain, const QDomDocument &aDocChain)
Replaces the existing chain by the chain imported from the given DOM document.
- [CxChain](#) * [loadChain](#) (const QDomDocument &aDocChain, QString *pErrorMsg=NULL)
Loads and creates a chain from the DOM document.
- const QMap< XICORE_HANDLE, [CxCameraSource](#) * > [cameraSources](#) () const
Gets all the camera sources.
- [CxChainable](#) * [findChainableByUniqueID](#) ([CxChain](#) *pChain, const QString &sUniqueID)
Returns a chainable with given UniqueId. If it is not present in the pChain then it returns NULL.

Static Public Member Functions

- static [CxChainManager](#) * [instance](#) ()
Gets the class instance.
- static bool [exportChain](#) ([CxChain](#) *pChain, QDomDocument &aXml, QDomElement &aRoot)
Exports the chain to DOM document.
- static bool [exportChain](#) ([CxChain](#) *pChain, QString &sXml)
Export the chain to XML string.
- static bool [saveChain](#) ([CxChain](#) *pChain, const QString &sFileName)
Saves the chain to XML file.
- static void [saveChainableSettings](#) ([CxChainable](#) *pObject, QSettings *pSettings)
Saved the object settings to QSettings.
- static void [loadChainableSettings](#) ([CxChainable](#) *pObject, QSettings *pSettings)
Loads the object settings from QSettings.

Private Slots

- void [onChainErrorState](#) (CxChain::ExChainErrorState eState, QString sErrMsg)

Private Member Functions

- [CxChainManager](#) ()
Private constructor (the class is singleton).
- void [removeBlindBranches](#) ([CxChain](#) *&pChain)
Removes branches with no listeners. If the chain is empty, then it will be destroyed.
- void [removeBlindBranches](#) ([CxChain](#) *&pChain, [CxChainable](#) *pChainable)
Removes the blind branches the pChainable is connected to.

Private Attributes

- QSet< [CxChain](#) * > [m_setChains](#)
- QMap< XICORE_HANDLE, [CxCameraSource](#) * > [m_mapCameraSources](#)
- QSet< [CxChainable](#) * > [m_setChainables](#)
- QMutex [m_lock](#)

4.10.1 Detailed Description

The class that creates, manages and manipulates with chains.

4.10.2 Member Function Documentation

4.10.2.1 `bool CxChainManager::branchAfter ( CxChain * pChain, CxChainable * pChainable, CxChainable * pAfter )`

Inserts the `pChainable` after the `pAfter` in the chain given by `pChain` as a new branch.

Remarks

The `pAfter` must already be connected in the `pChain`.

4.10.2.2 `CxCameraSource* CxChainManager::cameraSource ( XICORE_HANDLE hCamera )`

Gets the pointer to the existing camera source.

Returns

If the camera source exists then it returns pointer to it, otherwise it returns NULL.

4.10.2.3 `bool CxChainManager::closeChain ( CxChain * pChain )`

Closes the chain.

Remarks

If the chain is started it will be stopped. The method removes all the chainable object from the chain and then deletes the chain.

4.10.2.4 `bool CxChainManager::connectNonChainable ( CxChain * pChain, CxChainable * pChainable, QObject * pListener, const char * onDataAvailableSlot, const char * onFreePoolBuffersSlot )`

The method registers a custom object as a data receiver of the `pChainable`.

Remarks

Always use the SLOT macro to define the `onDataAvailableSlot` and `onFreePoolBuffersSlot`. These slots will be connected to [CxChainable::dataAvailable\(\)](#) signal and [CxMemoryManager::freePoolBuffers\(\)](#).

See also

[disconnectNonChainable\(\)](#)

4.10.2.5 `CxCameraSource* CxChainManager::createCameraSource ( XICORE_HANDLE hCamera )`

Creates an instance of camera source.

Remarks

Each camera can have only one instance of the camera source. If the method is called more than once for the valid camera handle then it fails (it returns NULL).

4.10.2.6 CxImageSource* CxChainManager::createImageSource (IxImageDataLoader * pLoader)

Creates an instance of an image source.

Remarks

The function take the ownership of the pLoader, so that it takes care about its deletion.

4.10.2.7 CxStaticImageSource* CxChainManager::createStaticImageSource (CxImageData * pImageData, const QString & sTitle)

Creates an instance of static image source.

Remarks

If the function succedes then it takes the ownership of the pImageData, so that it takes care about its deletion.

4.10.2.8 CxChain* CxChainManager::duplicateChainBranch (CxChain * pChain, CxChainable * pRoot, CxChainable * pLast, CxChainable * pNewRoot, bool bIncludeRoot)

Duplicates the branch of the pChain starting from pRoot to pLast. If the pLast is an instance of [CxViewImageReceiver](#) then the eventual image view attached to the receiver will be duplicated as well.

Remarks

All the duplicated members of the chain will be appended to the pNewRoot. If the bIncludeRoot is true, then the duplicated pRoot will be included in the duplicated chain as well.

4.10.2.9 CxChain* CxChainManager::duplicateChainWithViews (CxChain * pChain, CxChainable * pRoot, CxChainable * pNewRoot, bool bIncludeRoot)

Duplicates the pChain starting from pRoot. All views ([CxMdiView](#)) will be duplicated too.

Remarks

All the duplicated members of the chain will be appended to the pNewRoot. If the bIncludeRoot is true, then the duplicated pRoot will be included in the duplicated chain as well.

4.10.2.10 bool CxChainManager::importToExistingChain (CxChain * pChain, const QDomDocument & aDocChain)

Replaces the existing chain by the chain imported from the given DOM document.

Remarks

. The chain being imported will be appended to the existing image source of the pChain. The image source from the from the DOM will be excluded.

4.10.2.11 bool CxChainManager::insertAfter (CxChain * pChain, CxChainable * pChainable, CxChainable * pAfter)

Inserts the pChainable after the pAfter in the chain given by pChain.

Remarks

The pAfter must already be connected in the pChain. All successors of the pAfter will be disconnected fro it and reconnected to the pChainable.

4.10.2.12 `bool CxChainManager::insertBefore ( CxChain * pChain, CxChainable * pChainable, CxChainable * pBefore )`

Inserts the `pChainable` before the `pBefore` object in the chain given by `pChain`.

Remarks

The `pBefore` must already be connected in the `pChain`. All predecessors of the `pBefore` will be disconnected from it and reconnected to `pChainable`.

4.10.2.13 `bool CxChainManager::remove ( CxChain *& pChain, CxChainable * pChainable, bool bDeleteBranchWhenNoListener = false )`

Removes the `pChainable` from the `pChain`.

Remarks

If the `bDeleteBranchWhenNoListener` and the branch from which the `pChainable` was removed has no listener then the branch will be removed from the chain. If the chain has only one branch then the chain will be removed and the `pChain` will be set to `NULL`.

4.10.3 Member Data Documentation

4.10.3.1 `QMutex CxChainManager::m_lock [private]`

The lock of this object.

4.10.3.2 `QMap<XICORE_HANDLE, CxCameraSource*> CxChainManager::m_mapCameraSources [private]`

Map of camera sources

4.10.3.3 `QSet<CxChainable*> CxChainManager::m_setChainables [private]`

Set of chainables.

4.10.3.4 `QSet<CxChain*> CxChainManager::m_setChains [private]`


The set of chains.

The documentation for this class was generated from the following file:

- ChainManager.h

4.11 CxCircularBuffer Class Reference

Inheritance diagram for `CxCircularBuffer`:



Signals

- void **frameReceived** (qint32 iCirclePosition)

Public Member Functions

- virtual QString **title** () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual bool **acceptsDataFrom** (CxChainable *pPredecessor)
Query if this object can accept data from the given chainable object.
- virtual int **buffersCountInMemoryPool** () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the CxChainable.*
- virtual CxChainable * **clone** ()
Creates a deep copy of this object.
- virtual bool **currentOutputImageInfo** (SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL)
- virtual bool **queryOutputImageInfo** (const SxPicBufInfo &picInfoInput, SxPicBufInfo &picInfoOutput, const CxImageMetadata *pMetadataInput, CxImageMetadata *pMetadataOutput=NULL)
- virtual qint32 **sequenceSize** ()
In case the source has multiple data items, their count is returns here.
- virtual lxChainData * **loadedData** (qint32 iSeqIdx=0)
*Gets the pointer to the current data (after the **load()** method was called).*
- virtual bool **load** ()
Loads data from some source.
- virtual QString **inputFile** ()
If the loader load data from file then this method returns the the file's name.
- virtual void **setInputFile** (const QString &sFileName)
- virtual double **loadedFps** ()
In case the source is a sequence, this should get its fps. Return 0 if uncertain.
- virtual TxImageFormatList **supportedImageFormats** ()
Image formats capable of loading.
- void **setCircleSize** (qint32 iSize)
Sets the number of buffers in the circular buffer.
- qint32 **circleSize** () const

Protected Slots

- virtual void **onDataAvailable** (CxChainable *pSender, int iFrameNo)

Protected Member Functions

- virtual bool **init** (QString *pErrMsg=NULL)
*Initializes the object. When calling this method, the object **MUST** be stopped.*
- virtual bool **initialized** ()

Protected Attributes

- qint32 [m_iCircleSize](#)
number of allocated pictures (cannot be all filled - when stopped in the first pass)
- qint32 [m_iFramesCount](#)
number of frames captured to circular buffer
- qint32 [m_iLastFrame](#)
frame that we stored as last at this moment
- double [m_dFps](#)
fps of the incoming images. Need to tell when providing the data further
- QVector< [CxImageData](#) * > [m_vecCircBuffer](#)
this is the allocated circular buffer

Private Member Functions

- void **freeBuffers** ()

Additional Inherited Members

4.11.1 Detailed Description

Circular buffer has two modes:

1. as a [CxChainable](#), serves as a storage for frames (primary task)
2. as a [IxImageDataLoader](#), serves as a in-memory storage for [CxVideoSource](#)

4.11.2 Member Function Documentation

4.11.2.1 virtual bool CxCircularBuffer::init (QString * *pErrMsg* = NULL) [protected],[virtual]

Initializes the object. When calling this method, the object MUST be stopped.

Remarks

Methods [init\(\)](#) and [initialized\(\)](#) are internally used by the [start\(\)](#) method. It should initialize the chainable object so that it is ready to start ([start\(\)](#)) . If this object is the [CxImageProvider](#) then the method initializes internal memory pool.

Implements [CxChainable](#).

4.11.2.2 virtual IxChainData* CxCircularBuffer::loadedData (qint32 *iSeqIdx* = 0) [virtual]

Gets the pointer to the current data (after the [load\(\)](#) method was called).

Remarks

The caller of this method is the owner of the returned data so that it has to take care about its deletion using the `delete` operator.


Implements [IxChainDataLoader](#).

The documentation for this class was generated from the following file:

- CircularBuffer.h

4.12 CxClickableLabel Class Reference

Inheritance diagram for CxClickableLabel:



Signals

- void **clicked** ()

Public Member Functions

- **CxClickableLabel** (QWidget *parent=0)

Protected Member Functions


- virtual void **mouseReleaseEvent** (QMouseEvent *event)

The documentation for this class was generated from the following file:

- ClickableLabel.h

4.13 CxColorPickerLabel Class Reference

Inheritance diagram for CxColorPickerLabel:



Signals

- void **colorChanged** ()

Public Member Functions

- **CxColorPickerLabel** (QWidget *parent=0)
- unsigned **pickerColor** ()
- void **setPickerColor** (unsigned clColor)

Static Public Member Functions

- static QCursor **droppperCursor** ()

Protected Member Functions

- virtual void **paintEvent** (QPaintEvent *)

Protected Attributes

- unsigned **m_color**

Private Slots

- void **on_pickerClicked** ()

The documentation for this class was generated from the following file:

- ClickableLabel.h

4.14 CxLut::CxCompContrast Class Reference

Public Member Functions

- bool **isReset** (quint32 uiSrcBpc)
does it applies contast?
- QString **saveToString** () const
- bool **loadFromString** (const QString &sString)

Public Attributes

- int **m_iSrcMinVal**
- int **m_iSrcMaxVal**
increase contrast
- int **m_iDstMinVal**
- int **m_iDstMaxVal**
decrease contrast
- double **m_dGamma**
gamma during contrast


The documentation for this class was generated from the following file:

- Lut.h

4.15 CxCorePlugin Class Reference

The core plugin that registers all core classes (image loaders/savers etc.)

Inheritance diagram for CxCorePlugin:



Public Member Functions

- virtual bool [init](#) ()
Called on application start.
- virtual QString [name](#) () const
Plugin name displayed in the Plugin Manager view.
- virtual QString [author](#) () const
Name of the author (company) of this plugin.
- virtual QString [description](#) () const
Short description of the plugin, notes, warnings etc.
- virtual void [version](#) (int &major, int &minor, int &build) const
Version of the plugin.
- virtual QString [website](#) () const
The author's website adress.
- virtual QString [copyrightNotice](#) () const
In case plugin needs to write some legal info to Help - About - Legal...

Static Public Member Functions

- static [CxPlugin](#) * [instance](#) ()

Private Attributes

- bool [m_bInitialized](#)

4.15.1 Detailed Description

The core plugin that registers all core classes (image loaders/savers etc.)

The documentation for this class was generated from the following file:

- CorePlugin.h

4.16 CxHsiCameraParams::CxCorrectionMatrix Class Reference

Public Attributes

- QString [m_sName](#)
Name of the correction matrix.

- int **m_iStartNm**
- int **m_iEndNm**
Filter range.
- **ExCorrectionMatrixType** **m_eType**
Type of correction matrix.
- **ExCorrectionMatrixAlgorithm** **m_eAlgorithm**
Algorithm used to create the matrix.
- QVector< **CxVirtualBand** > **m_vecVirtualBands**
size is equal or lower than vecBands.count, may be present even in images (correction of spectrum). tag <spectral↔_correction_info>
- QVector< **CxOpticalComponent** > **m_vecOpticalComponents**
List of optical components used in this matrix.


The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.17 CxDemosaic Class Reference

The class performs a debayering on images in RAW format.

Inheritance diagram for CxDemosaic:



Public Member Functions

- **CxDemosaic** (const **SxPicBufInfo** &rOutputFormat)
- virtual QString **title** () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual bool **acceptsDataFrom** (**CxChainable** *pPredecessor)
Query if this object can accept data from the given chainable object.
- virtual int **buffersCountInMemoryPool** () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the **CxChainable**.*
- virtual **CxChainable** * **clone** ()
Creates a deep copy of this object.
- virtual bool **hasHelp** () const
Returns true if the object has help. The default is false.
- virtual QString **help** (QString &sImagesPath)
Returns a string with help.
- virtual bool **queryOutputImageInfo** (const **SxPicBufInfo** &picInfoInput, **SxPicBufInfo** &picInfoOutput, const **CxImageMetadata** *pMetadataInput, **CxImageMetadata** *pMetadataOutput=NULL)
Query the output image info (format) for given input image format.

Protected Member Functions

- virtual void [setupStart](#) ()
Sets up the object before the object is started.
- virtual void [setupStop](#) ()
Sets up the object after the object is stopped.
- virtual [IxChainData](#) * [processData](#) ([IxChainData](#) *pReceivedData)
The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Private Types

- enum **ExOfflineProcessing** { **eopNone**, **eopUnpacking**, **eopDebayering** }

Private Attributes

- bool **m_bUsePredefOutput**
- [SxPicBufInfo](#) **m_aPredefOutput**
- void * **m_pProc**
- ExOfflineProcessing **m_eProcType**

Additional Inherited Members

4.17.1 Detailed Description

The class performs a debayering on images in RAW format.

4.17.2 Member Function Documentation

4.17.2.1 virtual QString CxDemosaic::help (QString & *sImagesPath*) [virtual]

Returns a string with help.

Remarks

The text can be plain text or HTML. If the HTML help contains images then *sImagesPath* should contain an path to the images. The images path must be either absolute or relative to the application path (`qApp->applicationDirPath()`) and it also can be a prefix to embedded resource images (for example `"/images"`).

Reimplemented from [CxChainable](#).

4.17.2.2 virtual IxChainData* CxDemosaic::processData (IxChainData * *pReceivedData*) [protected], [virtual]

The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Remarks

If the implementation creates (allocates) a new data here then the original (received) data should be deleted using the delete operator. The default implementation of this method returns a pointer to the received data without any processing.

Warning

If the method returns NULL, then the `pReceivedData` will be deleted automatically!

Reimplemented from [CxChainable](#).

4.17.2.3 `virtual void CxDemosaic::setupStart ( ) [protected],[virtual]`

Sets up the object before the object is started.

Remarks

The method is internally used by the [start\(\)](#) method immediately before the [run\(\)](#) method is invoked. For example, the [CxCameraSource](#) sets up internal frame counter to zero and starts camera acquisition. The default implementation does nothing.

Reimplemented from [CxChainable](#).

4.17.2.4 `virtual void CxDemosaic::setupStop ( ) [protected],[virtual]`

Sets up the object after the object is stopped.

Remarks

The method is internally used by the [stop\(\)](#) method immediately after the [run\(\)](#) method is finished. For example, the [CxCameraSource](#) stops camera acquisition here. The default implementation does nothing.


Reimplemented from [CxChainable](#).

The documentation for this class was generated from the following file:

- Demosaic.h

4.18 CxExportToolButton Class Reference

Inheritance diagram for CxExportToolButton:

**Public Types**

- enum **ExExportType** { `extImageToClipboard` = 0x01, `extImageToFile` = 0x02, `extDataToClipboard` = 0x04, `extDataToTextFile` = 0x08 }

Signals

- void **exportClicked** (int eType)

Public Member Functions

- **CxExportToolButton** (ExExportTypes aSupportedExportTypes, QWidget *parent=0)
- void **addSupportedExportType** (ExExportType eType)
- void **setExportType** (ExExportType eType)
- ExExportType **exportType** ()

Static Public Member Functions

- static QString **textForType** (ExExportType eType)
- static QIcon **iconForType** (ExExportType eType)

Private Slots

- void **on_exportMenuClicked** ()
- void **on_selfClicked** ()

Private Attributes

- QMenu * **m_pPopupMenu**
- ExExportType **m_eSelExportType**

The documentation for this class was generated from the following file:

- ExportToolButton.h

4.19 CxFFT Class Reference

The class covers several 1D and 2D FFT functions.

Static Public Member Functions

- static bool **Alloc2D** (const [SxPicBufInfo](#) &picSrc, void *&pCfg, void *&pFFTin, void *&pFFTout)
Preallocates a buffers needed for 2D FFT. The buffers must be freed using [Free2D\(\)](#)
- static void **Free2D** (void *pCfg, void *pFFTin, void *pFFTout)
Frees the temporary buffers allocated by [Alloc2D\(\)](#)
- static bool **Calc2D** (const [SxPicBuf](#) &picSrcMono, [SxPicBuf](#) &picRe, [SxPicBuf](#) &picIm, void *pCfg, void *p←FFTin, void *pFFTout)
Performs the 2D FFT on the `picSrcMono` image.
- static bool **Calc2D** (const [SxPicBuf](#) &picSrcMono, [SxPicBuf](#) &picRe, [SxPicBuf](#) &picIm)
Performs the 2D FFT on the `picSrcMono` image.
- static bool **Alloc1D** (int iLength, void *&pCfg, void *&pFFTin, void *&pFFTout)
Preallocates a buffers needed for 1D FFT. The buffers must be freed using [Free1D\(\)](#)
- static void **Free1D** (void *pCfg, void *pFFTin, void *pFFTout)
Frees the temporary buffers allocated by [Alloc1D\(\)](#)
- static bool **Calc1D** (float *pInput, float *pReal, float *pImaginary, void *pCfg, void *pFFTin, void *pFFTout)
Performs the 1D FFT on the `pInput` array.
- static bool **Calc1D** (float *pInput, float *pReal, float *pImaginary, int iLength)
Performs the 1D FFT on the `pInput` array.

4.19.1 Detailed Description

The class covers several 1D and 2D FFT functions.

4.19.2 Member Function Documentation

4.19.2.1 `static bool CxFFT::Calc1D ( float * pInput, float * pReal, float * plmaginary, void * pCfg, void * pFFTin, void * pFFTout ) [static]`

Performs the 1D FFT on the *pInput* array.

Remarks

The *pReal* and *pImaginary* are array of the same length as *pInput* and the length must be the same as passed to the [Alloc1D\(\)](#). The *pCfg*, *pFFTin* and *pFFTout* are buffers preallocated by [Alloc1D\(\)](#).

4.19.2.2 `static bool CxFFT::Calc1D ( float * pInput, float * pReal, float * plmaginary, int iLength ) [static]`

Performs the 1D FFT on the *pInput* array.

Remarks

The *pReal* and *pImaginary* are array of the same length as *pInput* given by *iLength*.

Warning

This function allocates the temporary buffers internally for each call so that it is less effective.

4.19.2.3 `static bool CxFFT::Calc2D ( const SxPicBuf & picSrcMono, SxPicBuf & picRe, SxPicBuf & picIm, void * pCfg, void * pFFTin, void * pFFTout ) [static]`

Performs the 2D FFT on the *picSrcMono* image.

Remarks

The *picRe* and *picIm* are images of the same size as *picSrc*. The image type of these images must be `extypeInt32`. The *pCfg*, *pFFTin* and *pFFTout* are buffers preallocated by [Alloc2D\(\)](#).

4.19.2.4 `static bool CxFFT::Calc2D ( const SxPicBuf & picSrcMono, SxPicBuf & picRe, SxPicBuf & picIm ) [static]`

Performs the 2D FFT on the *picSrcMono* image.

Remarks

The *picRe* and *picIm* are images of the same size as *picSrc*. The image type of these images must be `extypeInt32`.

Warning

This function allocates the temporary buffers internally for each call so that it is less effective. It is suitable for static images.


The documentation for this class was generated from the following file:

- `FFT.h`

4.20 CxFFTImageData Class Reference

Implementation of [CxImageData](#) that can keep the FFT on an image. It consists of real and imaginary parts both of type `extypeInt32`.

Inheritance diagram for CxFFTImageData:



Public Member Functions

- virtual [IxChainData](#) * [clone](#) (TxMemPoolHandle hMemPool, QObject *pMemPoolClient) const
Creates the deep copy of this object. If the `hMemPool` is not NULL and it is a valid memory pool, then the data will be allocated from the memory pool. If the `hMemPool` is NULL then the data will be allocated from heap (operator new).
- virtual void [setPicBuf](#) (const [SxPicBuf](#) &picBuf)
Set the real part of the FFT image. It must be one-channel image of type `extypeInt32`.
- bool [setPicBufIm](#) (const [SxPicBuf](#) &picIm)
Sets the imaginary part of the FFT image. It must be one-channel image of type `extypeInt32`.
- [SxPicBuf](#) & [picBufIm](#) ()
Reference to the imaginary part of the FFT image. Can be modified.
- const [SxPicBuf](#) & [picBufIm](#) () const
Const variant that should not change the returned buffer.

Protected Attributes

- [SxPicBuf](#) `m_aPicIm`

4.20.1 Detailed Description

Implementation of [CxImageData](#) that can keep the FFT on an image. It consists of real and imaginary parts both of type `extypeInt32`.

Remarks

The [setPicBuf\(\)](#) and [picBuf\(\)](#) methods sets/gets the real part of the FFT.

Warning

This class can keep only one-channel images!

The documentation for this class was generated from the following file:

- ImageData.h

4.21 CxHsiCameraParams::CxFilterZone Class Reference

Public Member Functions

- bool **isValid** () const
- [CxBand](#) * **bandAt** (int xPos, int yPos)
For xPos < m_iPatternWidth and yPos < m_iPatternHeight.
- const [CxBand](#) * **bandAt** (int xPos, int yPos, int *pidx=NULL) const
For xPos < m_iPatternWidth and yPos < m_iPatternHeight.
- const [CxBand](#) * **bandAtPicPos** (int xPos, int yPos, int *pidx=NULL) const
Use m_rcActiveArea to get the pattern pos.
- bool **isPicPosInZone** (int xPos, int yPos) const
Use m_rcActiveArea to determine if picture position belongs into this zone.
- [QVector](#)< int > **sortedBandIdxs** (bool bOnlySelected=false) const
Indexes to band vector sorted by band peak wavelength.
- int **selectedBandsCount** () const
Number of selected bands in the spectrum.
- [QVector](#)< int > **selectedBandsMapping** () const
Mapping of band index to selected band index. If a band is not selected then the mapping index is -1.

Public Attributes

- [ExPattern](#) **m_eFilterPattern**
Mosaic?
- int **m_iPatternWidth**
- int **m_iPatternHeight**
*Pattern size (4x4, 5x5 for mosaic, or 1*h for line scan)*
- int **m_iFilterWidth**
- int **m_iFilterHeight**
*Filter size (1x1 for mosaic, w*X for line scan) - repetition of same band inside pattern.*
- [QRect](#) **m_rcActiveArea**
Filter active area (optional, zero when N/A)
- int **m_iSpectralRangeStartNm**
- int **m_iSpectralRangeEndNm**
Active sensor range in nm (optional, zero when N/A)
- [QVector](#)< [CxBand](#) > **m_vecBands**

The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.22 CxImageFormatManager::CxFormatItem Class Reference

Public Attributes

- [QString](#) **m_sRttiClass**
CxRTTI class name.
- [TxImageFormatList](#) **m_lstFormats**
list of supported formats

The documentation for this class was generated from the following file:

- ImageFormatManager.h

4.23 CxGuiFunctions Class Reference

Static Public Member Functions

- static void [makeToolBarButtonsSameHeight](#) (QToolBar *pToolBar)
ToolBar buttons without icon may be smaller then those that have one - looks horrible.
- static QMessageBox::StandardButton [warningMsgBox](#) (bool &bShowNextTime, QWidget *parent, const QString &title, const QString &text, QMessageBox::StandardButtons buttons=QMessageBox::StandardButtons(QMessageBox::Ok), QMessageBox::StandardButton defaultButton=QMessageBox::NoButton)
Displays a warning message box with "Do not show next time this message" checkbox.
- static QMessageBox::StandardButton [informationMsgBox](#) (bool &bShowNextTime, QWidget *parent, const QString &title, const QString &text, QMessageBox::StandardButtons buttons=QMessageBox::StandardButtons(QMessageBox::Ok), QMessageBox::StandardButton defaultButton=QMessageBox::NoButton)
Displays an information message box with "Do not show next time this message" checkbox.
- static QMessageBox::StandardButton [criticalMsgBox](#) (bool &bShowNextTime, QWidget *parent, const QString &title, const QString &text, QMessageBox::StandardButtons buttons=QMessageBox::StandardButtons(QMessageBox::Ok), QMessageBox::StandardButton defaultButton=QMessageBox::NoButton)
Displays a critical message box with "Do not show next time this message" checkbox.
- static QPixmap [dpiScaledPixmap](#) (const QString &sPath, int iOrigWidth=0)
Image scaled to current DPI.

The documentation for this class was generated from the following file:

- GuiFunctions.h

4.24 CxHistoData Class Reference

Class holding measured histogram and pixel statistics data.

Classes

- class [CxPixelStats](#)
Class for measurement of data in the ROI/image per channel.

Public Types

- enum [ExHistoDataFlag](#) { [eMeasureHisto](#) = 0x01, [eMeasureStats](#) = 0x02, [eSkipRows](#) = 0x04, [eTreatRawAsMono](#) = 0x08 }

Public Member Functions

- [CxHistoData](#) & **operator=** (const [CxHistoData](#) &rhs)
- void **allocChannels** (int iTotalChannels, quint32 uiMaxVal, bool bIncludeStats)
- void **freeChannelData** ()
- bool **isChannelDataCompatible** (const [CxHistoData](#) &other) const
- void **calcPixelStats** (bool bEnableAlloc=false)
- void **sumWithData** (const [CxHistoData](#) &other)
- void **maxWithData** (const [CxHistoData](#) &other)
- void **divideBy** (quint32 uiVal)

Public Attributes

- quint32 [m_uiMaxVal](#)
max value according to bpc (i.e. 1024 for 10 bit images)
- ExColorFilterArray [m_eCfa](#)
RAW format of the image.
- QVector< quint32 * > [m_vecChannelData](#)
measured histogram per channel
- QVector< QColor > [m_vecChannelColor](#)
suggested channel colors
- QVector< QString > [m_vecChannelName](#)
suggested channel names
- QVector< [CxPixelStats](#) > [m_vecChannelStats](#)
measured pixel statistics
- QRect [m_rc](#)
rectangle that was measured

4.24.1 Detailed Description

Class holding measured histogram and pixel statistics data.

4.24.2 Member Enumeration Documentation

4.24.2.1 enum CxHistoData::ExHistoDataFlag

Enumerator

- eMeasureHisto** measure histogram (fill [m_vecChannelData](#))
- eMeasureStats** measure pixel statistics (fill [m_vecChannelStats](#))
- eSkipRows** skip rows in large images
- eTreatRawAsMono** RAW image are usually split into channels, keep as mono.

The documentation for this class was generated from the following file:

- PicBufMeas.h

4.25 CxHsiCameraParams Class Reference

Classes

- class [CxBand](#)
- class [CxCorrectionMatrix](#)
- class [CxFilterZone](#)
- class [CxOpticalComponent](#)
< Calibration data for used rejection filter (on camera, may not be present)
- class [CxVirtualBand](#)
< Band record inside <virtual_bands>
- struct [SxBandIdx](#)

Public Types

- enum **ExPattern** { epNone, epMosaic, epTiled, epWedge }
- enum **ExCorrectionMatrixType** { ecmtHyperspectral = 0, ecmtRgb = 1, ecmtPanchromatic = 2 }
- enum **ExCorrectionMatrixAlgorithm** { ecmaPitchfork = 0, ecmaTacktile = 1, ecmaBcls = 2 }
- enum **ExOpticalComponentType** {
 eoctShortpassFilter = 0, eoctLongpassFilter = 1, eoctBandpassFilter = 2, eoctNotchFilter = 3,
 eoctLinearPolarizer = 4, eoctCircularPolarizer = 5, eoctSource = 6 }
- typedef QVector< **SxBandIdx** > **TxBandIdxVector**

Public Member Functions

- bool **isValid** () const
- bool **isSimpleMosaicFilter** () const
Returns true if there is only one zone, and is mosaic.
- int **totalBandCountInAllZones** () const
Returns the sum of number of bands in all zones.
- const **CxFilterZone** * **filterZoneWithTotalBandIdx** (int iTotBandIdx, **SxBandIdx** *pIdx=NULL) const
Find the zone the band belongs to when iterating from 0 to totalBandCountInAllZones()
- const **CxFilterZone** * **filterZoneAtPicPos** (int xPos, int yPos) const
use m_rcActiveArea to get the pattern pos
- const **CxBand** * **bandAtPicPos** (int xPos, int yPos, **SxBandIdx** *pIdx=NULL) const
use m_rcActiveArea to get the pattern pos
- **TxBandIdxVector** **sortedBandIdxs** (bool bOnlySelected=false) const
Indexes to band vector sorted by band peak wavelength across all filter zones.
- void **removeDataOutsideFilterRange** ()
Removes band peaks and correction matrixes outside current filter range.
- void **removeBandCalibrationData** ()
Clears m_vecCalibData for all bands when no longer needed.
- bool **isCorrectionMatrixPresent** () const
Returns true if the correction to virtual bands is available.
- const **CxCorrectionMatrix** * **currentCorrectionMatrix** () const
Returns matrix if the correction to virtual bands for current filter range is available.
- bool **loadFromCameraXmlFile** (const QString &sFilename, QString *psErrorMsg=NULL)
Load from camera XML file.
- bool **loadFromXml** (QDomElement *pRootElement)
Load from XML in metadata or camera XML file.
- bool **saveToXml** (QDomDocument *pXmlDoc, QDomElement *pRootElement) const
Save to XML compatible with camera XML.
- void **loadFromHsiApiStruct** (const XI_HSI_CAMERA_PARAMS *pParams)

Static Public Member Functions

- static bool **isHsiCamera** (XICORE_HANDLE hCamera)

Public Attributes

- int **m_iCalibrationRangeStartNm**
- int **m_iCalibrationRangeEndNm**
Range of band calibration data (optional, zero when N/A)
- int **m_iCalibrationResNm**
Calibration data resolution (is usually 1 nm) (optional, zero when N/A)
- int **m_iFilterRangeStartNm**
- int **m_iFilterRangeEndNm**
Range of the active filter on camera (matches the filename, not read from XML)
- QVector< CxFilterZone > **m_vecFilterZones**
List of zones on the sensor with active areas, bands definitions (>= 1)
- QVector< CxCorrectionMatrix > **m_vecCorrectionMatrices**
Possible filter options for the camera sensor (as read from camera XMLs)

4.25.1 Member Enumeration Documentation

4.25.1.1 enum CxHsiCameraParams::ExCorrectionMatrixAlgorithm

Enumerator

- ecmaTacktile** correction matrix must be applied prior to reflectance calculation
- ecmaBcls** correction matrix must be applied after reflectance calculation

4.25.1.2 enum CxHsiCameraParams::ExCorrectionMatrixType

Enumerator

- ecmtHyperspectral** matrix generates corrected hyperspectral samples
- ecmtRgb** matrix generates true color RGB samples
- ecmtPanchromatic** matrix generates panchromatic samples

The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.26 CxHsiCorrectionImage Class Reference

Public Member Functions

- void **reset** ()
- bool **isEmpty** ()
- bool **setPicbuf** (const SxPicBuf &aPic, const CxImageMetadata &aMetadata, bool isWhiteImage)
The ownership of the aPic is passed to this object.
- bool **initDefault** (const SxPicBufInfo &aPicInfo, const CxImageMetadata &aMetadata)
- CxImageMetadata **metadata** ()
- bool **load** (const QString &sPath, bool blsWhite)
- bool **save** (const QString &sPath)
- const double * **meansPerBand** ()
- void **lock** ()
- void **unlock** ()

Static Public Member Functions

- static bool **apply** (CxHsiCorrectionImage &aWhite, CxHsiCorrectionImage &aDark, CxHsiCorrectionImage &aDarkBias, SxPicBuf &aPic, const CxImageMetadata &aPicMetadata, ExHsiCorrection eCorrType, QString *pErrMsg=NULL)
- static bool **validate** (CxHsiCorrectionImage &aWhite, CxHsiCorrectionImage &aDark, CxHsiCorrectionImage &aDarkBias, bool &bDarkValid, bool &bBiasValid, QString *pErrMsg=NULL)
- static bool **validate** (CxHsiCorrectionImage &aWhite, CxHsiCorrectionImage &aDark, CxHsiCorrectionImage &aDarkBias, SxPicBufInfo &aPicInfo, const CxImageMetadata &aPicMetadata, bool &bDarkValid, bool &bBiasValid, QString *pErrMsg=NULL)

Private Member Functions

- void **calculateMeansPerBand** ()

Static Private Member Functions

- static bool **validateNoLock** (CxHsiCorrectionImage &aWhite, CxHsiCorrectionImage &aDark, CxHsiCorrectionImage &aDarkBias, bool &bDarkValid, bool &bBiasValid, QString *pErrMsg=NULL)
- static bool **validateNoLock** (CxHsiCorrectionImage &aWhite, CxHsiCorrectionImage &aDark, CxHsiCorrectionImage &aDarkBias, SxPicBufInfo &aPicInfo, const CxImageMetadata &aPicMetadata, bool &bDarkValid, bool &bBiasValid, QString *pErrMsg=NULL)

Private Attributes

- SxPicBuf **m_pic**
- CxImageMetadata **m_mdata**
- double * **m_pMeansPerBand**
Mean values per each band in the center of the WHITE image. For dark images it is NULL.
- QMutex **m_lock**


The documentation for this class was generated from the following file:

- CameraParameters.h

4.27 CxImageData Class Reference

The class implements an image data object exchanged by CxChainable objects.

Inheritance diagram for CxImageData:



Public Member Functions

- [CxImageData](#) ()
Constructor.
- virtual [~CxImageData](#) ()
Destructor.
- virtual qint32 [frameNo](#) () const
Gets the frame number of the data. For "static" data (e.g. static image) it is always zero.
- virtual [IxChainData](#) * [clone](#) (TxMemPoolHandle hMemPool, QObject *pMemPoolClient) const
Creates the deep copy of this object. If the hMemPool is not NULL and it is a valid memory pool, then the data will be allocated from the memory pool. If the hMemPool is NULL then the data will be allocated from heap (operator new).
- virtual [IxChainData](#) * [clone](#) (QObject *pMemPoolClient) const
Creates the deep copy of this object from the same memory source as the original data.
- virtual bool [isFromMemoryPool](#) (TxMemPoolHandle hMemPool) const
Checks if the data is from given memory pool.
- virtual [SxPicBuf](#) & [picBuf](#) ()
Reference to the picture. Can be modified.
- virtual const [SxPicBuf](#) & [picBuf](#) () const
Const variant should not change the returned buffer.
- virtual void [setPicBuf](#) (const [SxPicBuf](#) &[picBuf](#))
- virtual const [CxImageMetadata](#) & [imageMetadata](#) () const
- virtual void [setImageMetadata](#) (const [CxImageMetadata](#) &aMetadata)

Protected Attributes

- [SxPicBuf](#) [m_aPic](#)
- [CxImageMetadata](#) [m_aMetadata](#)

4.27.1 Detailed Description

The class implements an image data object exchanged by [CxChainable](#) objects.

4.27.2 Member Data Documentation

4.27.2.1 [SxPicBuf](#) [CxImageData::m_aPic](#) [protected]

The image.

The documentation for this class was generated from the following file:

- [ImageData.h](#)

4.28 CxImageFormat Class Reference

Public Types

- enum [EType](#) { [eSingle](#) = 0x01, [eSequence](#) = 0x02 }

Public Member Functions

- [CxImageFormat](#) (const QString &sFormat, int iType)

Public Attributes

- QString [m_sFormat](#)
Should be in format "Name (.ext)".*
- int [m_iType](#)
Combination of EType flags.

4.28.1 Member Enumeration Documentation

4.28.1.1 enum CxImageFormat::EType

Enumerator

- eSingle** Format can be used for saving single image.
- eSequence** Format can be used for saving sequence of images.


The documentation for this class was generated from the following file:

- ImageFormat.h

4.29 CxImageFormatManager Class Reference

Keeps list of supported file formats using [CxRTTI](#), creates loader and saver objects.

Inheritance diagram for CxImageFormatManager:



Classes

- class [CxFormatItem](#)

Public Member Functions

- QStringList [formatsForLoad](#) (int iFormatTypes=[CxImageFormat::eSingle](#))
All formats we are capable of loading, items as "Name (.ext)".*
- QStringList [formatsForSave](#) (int iFormatTypes=[CxImageFormat::eSingle](#))
All formats we are capable of saving, items as "Name (.ext)".*
- QString [stringWithAllFormatForLoad](#) (int iFormatTypes=[CxImageFormat::eSingle](#), bool bIncludeAllImages↔
Filter=true)
All formats we are capable of loading in format suitable for QFileDialog.
- QString [stringWithAllFormatForSave](#) (int iFormatTypes=[CxImageFormat::eSingle](#))
All formats we are capable of saving in format suitable for QFileDialog.
- [IxImageDataLoader](#) * [loaderForFile](#) (const QString &sFileName)
Creates new image loader instance according to the file extension.
- [IxImageDataSaver](#) * [saverForFile](#) (const QString &sFileName)
Creates new image saver instance according to the file extension.

- TxImageFormatOptions [formatOptions](#) (const QString &sExt)
get defaults for file format
- void [setFormatOptions](#) (const QString &sExt, const TxImageFormatOptions &aOpt)
set defaults for file format
- void [saveSettings](#) (QSettings *pSettings)
store format options to app settings
- void [loadSettings](#) (QSettings *pSettings)
load format options from app settings

Static Public Member Functions

- static [CxImageFormatManager](#) * [instance](#) ()
Singleton class, this is the way to obtain the pointer.
- static QString [fileExtFromFormatSignature](#) (const QString &sFormatItem)
returns "ext" from "name (.ext)"*
- static bool [loadPicBufFromFile](#) (const QString &sFileName, [SxPicBuf](#) &rPic, [CxImageMetadata](#) *pMetadata=0)
Helper function to load one picture. The result is stored to rPic, caller is responsible freeing it.
- static bool [savePicBufToFile](#) (const QString &sFileName, const [SxPicBuf](#) *pPic, const [CxImageMetadata](#) *pMetadata=0)
Helper function to store one picture.

Private Slots

- void [onRttiltemEnableStateChanged](#) (TxRttiltem hItem, bool bEnabled)
Handle the [CxRTTI](#) changes.

Private Member Functions

- void [resetFormatsForLoad](#) ()
- void [resetFormatsForSave](#) ()
- void [loadAvailableLoaders](#) ()
Asks our [CxRTTI](#) to return all available loader classes.
- void [loadAvailableSavers](#) ()
Asks our [CxRTTI](#) to return all available saver classes.

Static Private Member Functions

- static QString [findClassForFormat](#) (const QList< [CxFormatItem](#) > &lst, const QString &sExt)

Private Attributes

- bool [m_bLoaderFormatsLoaded](#)
Did we call [loadAvailableLoaders\(\)](#)?
- bool [m_bSaverFormatsLoaded](#)
Did we call [loadAvailableSavers\(\)](#)?
- QList< [CxFormatItem](#) > [m_lstFormatsForLoad](#)
- QList< [CxFormatItem](#) > [m_lstFormatsForSave](#)
- QMap< QString, TxImageFormatOptions > [m_mapFormatOptions](#)
maps "ext" to its options (not class!)

4.29.1 Detailed Description

Keeps list of supported file formats using [CxRTTI](#), creates loader and saver objects.

The documentation for this class was generated from the following file:

- [ImageFormatManager.h](#)

4.30 CxImageMetadata Class Reference

Metadata stored next to image pixel data.

Public Member Functions

- bool **isColorCorrMatrixEmpty** () const
- void **fillFromCameraSettings** ([CxCameraParameters](#) *pPars)
Reads metadata values from current (cached) camera parameters.
- void **saveToXml** (QDomDocument *pXmlDoc, QDomElement *pRootElement) const
Saves all values to XML document.
- bool **loadFromXml** (QDomElement *pRootElement)
Loads all values from XML document.
- QString **saveToXmlString** () const
Calls saveToXml, and saves the XML doc to string.
- bool **loadFromXmlString** (const QString &sString)
Loads all values from XML doc stored in a string.

Static Public Member Functions

- static void **saveTag** (QDomDocument *pXmlDoc, QDomElement *pRootElement, const QString &sTag, const QString &sVal)

Public Attributes

- QString **m_sCameraModel**
String got from XI_PRM_DEVICE_NAME.
- int **m_iCameraModelId**
Integer got from XI_PRM_DEVICE_MODEL_ID.
- QString **m_sCameraSerialNumber**
String got from XI_PRM_DEVICE_SN.
- QString **m_sSensorSerialNumber**
String got from XI_PRM_DEVICE_SENS_SN.
- QString **m_sCameraUserName**
String got from XI_PRM_DEVICE_USER_ID.
- QString **m_sApiContextList**
String got from XI_PRM_API_CONTEXT_LIST.
- bool **m_bAutoExposure**
Is auto exposure switched on?
- int **m_iExposure**
In microseconds - as in xiApi.
- float **m_fGain**

Gain in dB, as in xiApi.

- bool [m_bAutoWB](#)

Is AWB switched on?

- float **m_fWhiteR**
- float **m_fWhiteG**
- float [m_fWhiteB](#)

White balance coefficients.

- float [m_fColorCorrMatrix](#) [16]

Current color correction matrix.

- float **m_fGammaY**
- float **m_fGammaC**
- float [m_fSharpness](#)

Other color correction parameters.

- int [m_iSensorTaps](#)

Number of sensor taps used when capturing the image.

- int [m_iDownsamplingType](#)

Value of XI_PRM_DOWNSAMPLING_TYPE.

- float **m_fLensAperture**
- float [m_fLensFocalLength](#)

Lens parameters (when present, 0 for none)

- float **m_fTempHousing**
- float [m_fTempSensor](#)

Housing and sensor temperatures in deg Celsius.

- ExColorFilterArray [m_eColorFilterArray](#)

Type of Bayer matrix, excfaNone when not RAW image.

- uint [m_uTransportDataFormat](#)

The GenTL image format (see XI_GenTL_Image_Format_e in xiApi.h).

- CxHsiCameraParams [m_aHsiParams](#)

Hyper Spectral Image information.

- QDateTime [m_aAcqDateTime](#)

Date and time of exposure start.

4.30.1 Detailed Description

Metadata stored next to image pixel data.


The documentation for this class was generated from the following file:

- ImageMetadata.h

4.31 CxImageProvider Class Reference

The base class for any chainable object which provides images.

Inheritance diagram for CxImageProvider:



Public Member Functions

- **CxImageProvider** (ExChainableFlags eFlags)
- virtual bool [queryOutputImageInfo](#) (const [SxPicBufInfo](#) &picInfoInput, [SxPicBufInfo](#) &picInfoOutput, const [CxImageMetadata](#) *pMetadataInput, [CxImageMetadata](#) *pMetadataOutput=NULL)=0
Query the output image info (format) for given input image format.
- virtual bool [currentOutputImageInfo](#) ([SxPicBufInfo](#) &picInfo, [CxImageMetadata](#) *pMetadata=NULL)
Gets the current output image info (format) provided by this object.
- virtual bool [currentInputImageInfo](#) ([SxPicBufInfo](#) &picInfo, [CxImageMetadata](#) *pMetadata=NULL)
Gets the current input image info (format) which this object receives on its input.
- virtual int [queryImagesCountInMemoryPool](#) (const [SxPicBufInfo](#) &picInfo, TxMemPoolHandle hPool)
Gets the required number of images reserved in the memory pool by its successors.

Protected Member Functions

- virtual bool [init](#) (QString *pErrMsg=NULL)
Initializes the object. When calling this method, the object MUST be stopped.
- virtual bool [initialized](#) ()

Private Member Functions

- int [totalBuffersCountInMemoryPool](#) (TxMemPoolHandle hPool)
Returns the total buffers count required by this object and its successors. If this object does not have its own memory pool then it returns zero.

Additional Inherited Members

4.31.1 Detailed Description

The base class for any chainable object which provides images.

4.31.2 Member Function Documentation

- 4.31.2.1 virtual bool CxImageProvider::currentOutputImageInfo ([SxPicBufInfo](#) & *picInfo*, [CxImageMetadata](#) * *pMetadata* = NULL) [virtual]

Gets the current output image info (format) provided by this object.

Remarks

The method should call `getCurrentOutputImageInfo()` on its predecessors and calculate the output image info provided by this object.

Reimplemented in [CxCameraSource](#), [CxStaticImageSource](#), and [CxVideoSource](#).

4.31.2.2 `virtual bool CxImageProvider::init ( QString * pErrMsg = NULL ) [protected],[virtual]`

Initializes the object. When calling this method, the object MUST be stopped.

Remarks

Methods `init()` and `initialized()` are internally used by the `start()` method. It should initialize the chainable object so that it is ready to start (`start()`). If this object is the [CxImageProvider](#) then the method initializes internal memory pool.

Implements [CxChainable](#).

4.31.2.3 `virtual int CxImageProvider::queryImagesCountInMemoryPool ( const SxPicBufInfo & picInfo, TxMemPoolHandle hPool ) [virtual]`

Gets the required number of images reserved in the memory pool by its successors.

Remarks

If this object uses its own memory pool then the method allways returns zero!

Reimplemented in [CxViewImageReceiver](#).


The documentation for this class was generated from the following file:

- [Chainable.h](#)

4.32 CxImageSource Class Reference

The base class of any image provider which is an image source, i.e. it has no predecessors and it provides images.

Inheritance diagram for `CxImageSource`:



Public Member Functions

- **CxImageSource** (ExChainableFlags eFlags)
- virtual bool [acceptsDataFrom](#) ([CxChainable](#) *pPredecessor)
Query if this object can accept data from the given chainable object.
- virtual bool [addPredecessor](#) ([CxChainable](#) *pPredecessor)
Adds new predecessor to this object. It is internally used by [CxChain](#) class.
- virtual bool [queryOutputImageInfo](#) (const [SxPicBuflInfo](#) &picInfoInput, [SxPicBuflInfo](#) &picInfoOutput, const [CxImageMetadata](#) *pMetadataInput, [CxImageMetadata](#) *pMetadataOutput=NULL)
Query the output image info (format) for given input image format.

Additional Inherited Members

4.32.1 Detailed Description

The base class of any image provider which is an image source, i.e. it has no predecessors and it provides images.
The documentation for this class was generated from the following file:

- [Chainable.h](#)

4.33 CxShadingPreset::CxItem Class Reference

Public Attributes

- QString [m_sXiApiName](#)
Parameter name.
- [CxValue](#) [m_vVal](#)
Parameter value.

The documentation for this class was generated from the following file:

- CameraParameters.h

4.34 CxLut Class Reference

Classes

- class [CxCompContrast](#)
- class [CxPseudoColorTable](#)

Public Types

- enum **ExDemosaicMode** { **exdemNone** = 0, **exdemMin** = 1, **exdemMax** = 2, **exdemAvg** = 3 }

Public Member Functions

- **CxLut** (quint32 uiInBpc, quint32 uiInComps)
- bool [alloc](#) ()
allocate the mapping table

- bool **allocRgb** (quint32 uiRgbBits)
allocate the mapping table to RGB24, RGB32 or RGB565 format
- CxLut * **transformToFormat** (quint32 uiInBpc, quint32 uiInComps) const
create new Lut based on current setting
- bool **isCompatibleWith** (const SxPicBuf *pDstPic, const SxPicBuf *pSrcPic)
- bool **apply** (SxPicBuf *pDstPic, const SxPicBuf *pSrcPic, const CxImageMetadata *pSrcMetadata)

Static Public Member Functions

- static quint32 **fullMask** (quint32 uiComps)
- static quint32 **numCompsInMask** (quint32 uiMask)
- static int **firstComplnMask** (quint32 uiMask)
- static bool **addPseudoColorLut** (const QString &sFile)
- static int **pseudoColorLutCount** ()
- static QString **pseudoColorLutName** (int idx)
- static int **indexOfPseudoColorLutName** (const QString &sName)
- static void **setCustomPseudoColorLutFolder** (const QString &sFolderName)
- static QString **customPseudoColorLutFolder** ()
- static void **reloadCustomPseudoColorLuts** ()

Public Attributes

- quint32 **m_uiSrcBpc**
bits per component on input
- quint32 **m_uiSrcComps**
number of input components
- quint32 **m_uiCompMask**
binary mask of input comps idxs to map
- quint32 **m_uiDstBpc**
bits per component on output
- quint32 **m_uiDstComps**
number of output components
- bool **m_bSingleComplnColor**
show masked channel in its color or in gray
- int **m_iPseudoColor**
in case one channel is mapped, use this pseudocolor map (-1 for none)
- quint32 **m_clMonoCompColor**
when mono image with m_bSingleComplnColor, use this color
- QVector< CxCompContrast > **m_vecContrast**
curve mapping per component. If size==1, same applied to all comps. Pre-created with size 1
- bool **m_bMarkOverexposed**
- bool **m_bMarkUnderexposed**
mark over- or underexposed pixels?
- quint32 **m_clOverexposed**
- quint32 **m_clUnderexposed**
color for over- and underexposed pixels (default: red, blue)
- enum CxLut::ExDemosaicMode **m_eDemosaicMode**

Protected Types

- enum [ExMapMode](#) { **exmmUnset**, **exmmToRgb**, **exmmSeparateComps**, **exmmSeparate3Comps**, **exmmSeparateCompsSameTable**, **exmmSingleChannelToRgb** }
- enum [ExLutDataType](#) { **exldtNull**, **exldtUint8Ptr**, **exldtUint16Ptr**, **exldtUint32Ptr** }

Protected Member Functions

- void [init](#) ()
called from all constructors
- void **freeLutData** ()
- bool **applyOnOneCompToRgb** ([SxPicBuf](#) *pDstPic, const [SxPicBuf](#) *pSrcPic, int iComp, const [CxImageMetadata](#) *pSrcMetadata)
- void **markOverExposedRgb32** ([SxPicBuf](#) *pPic)
- void **markOverExposedRgb565** ([SxPicBuf](#) *pPic)

Protected Attributes

- bool **m_bWasAllocatedForRgb**
- enum [CxLut::ExMapMode](#) **m_eMapMode**
- void * **m_pLutVoidData**
- enum [CxLut::ExLutDataType](#) **m_eLutDataType**
- [SxPicBuf](#) * **m_pPicDemosaicCache**

4.34.1 Member Enumeration Documentation

4.34.1.1 enum [CxLut::ExLutDataType](#) [protected]

Enumerator

- exldtUint8Ptr** mapping components to 8-bit component
- exldtUint16Ptr** typically mapping single component to RGB16
- exldtUint32Ptr** typically mapping single component to RGB32

4.34.1.2 enum [CxLut::ExMapMode](#) [protected]

Enumerator

- exmmToRgb** try to convert the picture to RGB just using generating functions (no tables)
- exmmSeparateComps** map each component separately, each component has its own table
- exmmSeparate3Comps** like **exmmSeparateComps**, but use only first 3 components (for RGB32)
- exmmSeparateCompsSameTable** map each component separately, while using same table for all
- exmmSingleChannelToRgb** color is table for one comp, producing RGB image


The documentation for this class was generated from the following file:

- [Lut.h](#)

4.35 CxMdiView Class Reference

The base of all MDI subwindows in the application.

Inheritance diagram for CxMdiView:



Public Types

- enum [ExViewCap](#) {
[evcCanZoomInOut](#) = 0x01, [evcCanZoomFit](#) = 0x02, [evcCanSaveToFile](#) = 0x04, [evcCanCapture](#) = 0x10,
[evcCanRecordVideo](#) = 0x20 }

Public Slots

- virtual void [zoomIn](#) ()
Change zoom as user clicked on the toolbar button.
- virtual void [zoomOut](#) ()
Change zoom as user clicked on the toolbar button.
- virtual void [zoomFitToScreen](#) ()
Change zoom as user clicked on the toolbar button.
- virtual void [updateViewTitle](#) ()

Signals

- void [receivedFrame](#) (qint32 iFrame)
New data arrived through chain to the view.
- void [requestAdjustWindowSize](#) ([CxMdiView](#) *pSender)
Emit when parent MdiWindow should adjust its size to this view's "sizeHint() const".
- void [requestCloseWindow](#) ([CxMdiView](#) *pSender)
For internal use only, called from closeMdiView()

Public Member Functions

- CxMdiView** ([CxChainable](#) *pDataProvider, [CxChainable](#) *pDataReceiver)
- void [closeMdiView](#) ()
- int [viewId](#) ()
Returns ID of this view.
- void [setViewId](#) (int iViewId)
Do not use, main window sets this automatically.
- virtual [ExViewCaps](#) [viewCaps](#) ()
- virtual bool [containsData](#) ()=0
Check whether view shows any data (maybe image view did not receive any imageData)
- virtual bool [canZoomIn](#) ()
- virtual bool [canZoomOut](#) ()
- virtual bool [save](#) (QString &sLastUsedFolder)

- Called from toolbar button, you should show the SaveAs dialog, or save the data directly.*
- virtual bool **saveToFile** (const QString &sFilename)
Save contents to file.
 - virtual bool **recordVideoToFile** (const QString &sFilename)
Start recording the video.
 - virtual CxMdiView * **clone** (CxChainable *pNewDataProvider, CxChainable *pNewDataReceiver)=0
 - CxChainable * **dataProvider** () const
Returns chainable object view is connected to.
 - CxChainable * **dataReceiver** () const
Returns view's helper object (that is connected to dataProvider)
 - CxImageSource * **imageSource** () const
Image source in the chain view is connected to.
 - virtual QString **viewTitle** () const
Returns the view title excluding additional decorations (i.e. w/ zoom factor, fps)

Protected Attributes

- CxChainable * **m_pDataReceiver**
The chainable object as a data receiver.
- int **m_iViewId**

4.35.1 Detailed Description

The base of all MDI subwindows in the application.

4.35.2 Member Enumeration Documentation

4.35.2.1 enum CxMdiView::ExViewCap

Enumerator

- evcCanZoomInOut** Enable zooming using + / - toolbar buttons.
- evcCanZoomFit** Enable zooming using bestFit toolbar button.
- evcCanSaveToFile** This view can save some data to file, enable toolbar button, and calling.
- evcCanCapture** Is reasonable to make snaphots of the data inside the view.
- evcCanRecordVideo** This view can record video.


The documentation for this class was generated from the following file:

- MdiView.h

4.36 CxMemoryManager Class Reference

The singleton class for management of memory pools.

Inheritance diagram for CxMemoryManager:



Classes

- struct [SxLimitDef](#)

Public Types

- enum **ExAppOperation** { **eaopStartProgram**, **eaopStartAcquisition**, **eaopStartVideoRecord** }

Public Slots

- void [finishAll](#) ()
Destroys all memory pools.

Signals

- void [freePoolBuffers](#) (TxMemPoolHandle hMemPool)
The signal is emitted when the memory pool is freed and it still contains some allocated (non-released) buffers. The receiver of this signal MUST immediately release these buffers by calling [freePicBuf\(\)](#). It is very important to connect this signal with direct connection (see `QObject::connect()` and `Qt::DirectConnection`).

Public Member Functions

- [~CxMemoryManager](#) ()
Destructor.
- TxMemPoolHandle [createPool](#) ()
Creates new empty memory pool.
- bool [initPool](#) (TxMemPoolHandle hPool, const [SxPicBufInfo](#) &picInfo, int iCount, int &rAllocatedCount)
Initialization of the existing empty pool.
- void [freePool](#) (TxMemPoolHandle hPool)
Releases the pool's internal buffers. It means that the memory pool itself will not be deleted. After calling this method the memory pool will be empty and it won't be possible to allocate any buffer from it. The method succeeds if the pool has all buffers free.
- bool [deletePool](#) (TxMemPoolHandle hPool)
Tries to delete the memory pool.
- int [buffersCount](#) (TxMemPoolHandle hPool)
Gets total number of buffers kept by this pool.
- int [freeBuffersCount](#) (TxMemPoolHandle hPool)
Gets free buffers kept by this pool.
- bool [poolSize](#) (TxMemPoolHandle hPool, size_t &uiSize)
Gets to total allocated size of the pool..
- bool [bufferInfo](#) (TxMemPoolHandle hPool, [SxPicBufInfo](#) &picInfo)
Gets buffer info for which the pool was created.
- void [reserveBuffersInPool](#) (TxMemPoolHandle hPool, QObject *pClient, int iRequestCount)
Objects (chainables) can reserve some buffer places for themselves, so that cloning imageData would not fail.
- void [invalidateBufferReservations](#) (TxMemPoolHandle hPool)
Clear all reservations made in pool.
- bool [allocPicBuf](#) (TxMemPoolHandle hPool, QObject *pClient, [SxPicBuf](#) &picBuf, bool bReset)
Allocates a buffer from given pool.
- bool [freePicBuf](#) ([SxPicBuf](#) &picBuf)
Tries to releases the image buffer.
- void [freeAllBuffers](#) (TxMemPoolHandle hPool)

The method releases all buffers in the memory pool but will not dealocate the memory pool itself. If the pool is non empty then it emits the [freePoolBuffers\(\)](#) signal.

- bool [isImageCompatibleWithPool](#) (TxMemPoolHandle hMemPool, const [SxPicBufInfo](#) &picInfo)
Checks whether the given image info is compatible with the memory pool. It means whether it is possible to allocate the a buffer with given image info from the memory pool.
- bool [loadLimits](#) ()
Loads the limits from XML file.
- bool [checkLimits](#) (ExAppOperation eOperation, XICORE_HANDLE hCam, QString &sMessage, const Q←String *pParamName=NULL, const [CxValue](#) *pParValue=NULL)
Check a resources for given operation and camera.

Static Public Member Functions

- static [CxMemoryManager](#) * [instance](#) ()
Gets the singleton instance of the class.
- static quint64 [availableMemory](#) ()
Returns the amount of available RAM in bytes.
- static quint32 [L2CacheSize](#) ()
Returns the size of L2 cache in bytes.

Private Slots

- void [onFreePoolBuffers](#) (TxMemPoolHandle hMemPool)
Reaction on signal [CxMemoryPool::freeBuffersRequest\(\)](#)

Private Member Functions

- [CxMemoryManager](#) ()
Private constructor.

Private Attributes

- QMap< ExAppOperation, QList< [SxLimitDef](#) > > [m_mapLimitDefs](#)
- QSet< TxMemPoolHandle > [m_setMemPools](#)
- QMutex [m_lock](#)

4.36.1 Detailed Description

The singleton class for management of memory pools.

4.36.2 Member Function Documentation

- 4.36.2.1 bool [CxMemoryManager::allocPicBuf](#) (TxMemPoolHandle *hPool*, QObject * *pClient*, [SxPicBuf](#) & *picBuf*, bool *bReset*)

Allocates a buffer from given pool.

Remarks

The `picBuf` parameter must fit the `picInfo` parameter passed to `createPool()`. The method actualizes these members of the `picBuf` :

- `SxPicBuf::m_pData`
- `SxPicBuf::m_uiAllocSize`
- `SxPicBuf::m_bDataOwner`

4.36.2.2 `bool CxMemoryManager::bufferInfo ( TxMemPoolHandle hPool, SxPicBufInfo & picInfo )`

Gets buffer info for which the pool was created.

See also

`createPool()`

4.36.2.3 `int CxMemoryManager::buffersCount ( TxMemPoolHandle hPool )`

Gets total number of buffers kept by this pool.

Returns

If the method fails then it return -1.

4.36.2.4 `bool CxMemoryManager::checkLimits ( ExAppOperation eOperation, XICORE_HANDLE hCam, QString & sMessage, const QString * pParamName = NULL, const CxValue * pParValue = NULL )`

Check a resources for given operation and camera.

Remarks

If it is necessary to check an parameter and its value that is not set on the given camera yer, pass it via `pParamName` and `pParValue`. It is suitable for checking resources before parameter change.

4.36.2.5 `TxMemPoolHandle CxMemoryManager::createPool ( )`

Creates new empty memory pool.

See also

`deletePool()`
`initPool()`
`freePool()`

4.36.2.6 `bool CxMemoryManager::deletePool ( TxMemPoolHandle hPool )`

Tries to delete the memory pool.

Remarks

The method internally calls `freePool()` and if it succeeds then the pool is deleted and the handle `hPool` became invalid.

4.36.2.7 void CxMemoryManager::finishAll () [slot]

Destroys all memory pools.

Remarks

The method forces the memory pools destruction without checking if the pools contain only free (nonoccupied) buffers so that the method should be used very carefully. It should be called when the application is closing.

4.36.2.8 int CxMemoryManager::freeBuffersCount (TxMemPoolHandle *hPool*)

Gets free buffers kept by this pool.

Returns

If the method fails then it returns -1.

4.36.2.9 bool CxMemoryManager::freePicBuf (SxPicBuf & *picBuf*)

Tries to release the image buffer.

Remarks

The function releases the buffer only when the SxPicBuf::m_hDataOwner contains a valid handle to memory pool.

4.36.2.10 void CxMemoryManager::freePool (TxMemPoolHandle *hPool*)

Releases the pool's internal buffers. It means that the memory pool itself will not be deleted. After calling this method the memory pool will be empty and it won't be possible to allocate any buffer from it. The method succeeds if the pool has all buffers free.

Remarks

The method does not delete the pool itself. The handle stays valid. It should be called before initialization of the already initialized pool.

See also

[deletePool\(\)](#)
[initPool\(\)](#)

4.36.2.11 bool CxMemoryManager::initPool (TxMemPoolHandle *hPool*, const SxPicBufInfo & *picInfo*, int *iCount*, int & *rAllocatedCount*)

Initialization of the existing empty pool.

Remarks

The parameter `picInfo` must have set these members correctly:

- [SxPicBufInfo::m_uiWidth](#)
- [SxPicBufInfo::m_uiHeight](#)
- [SxPicBufInfo::m_uiStride](#)
- [SxPicBufInfo::m_uiChannels](#).

4.36.3 Member Data Documentation

4.36.3.1 QMutex CxMemoryManager::m_lock [private]

Lock of this object.

4.36.3.2 QSet<TxMemPoolHandle> CxMemoryManager::m_setMemPools [private]

Map of memory pools and its threads.

The documentation for this class was generated from the following file:

- MemoryManager.h

4.37 CxHsiCameraParams::CxOpticalComponent Class Reference

< Calibration data for used rejection filter (on camera, may not be present)

Public Member Functions

- bool **isValid** () const
- void **clear** ()
- QString **typeString** () const

Public Attributes

- ExOpticalComponentType **m_eType**
- QString **m_sManufacturer**
Manufacturer of the optical component.
- QString **m_sPartId**
Textual description of the component.
- QString **m_sDescription**
Unique part identifier as given by the manufacturer.
- QString **m_sTag**
Custom tag uniquely identifying the component.
- int **m_iStartNm**
- int **m_iEndNm**
Range of band calibration data.
- int **m_iResNm**
Calibration data resolution (is usually 1 nm)
- QVector< double > **m_vecCalibData**
Band response data, tag <response_composition> (not required for images)

4.37.1 Detailed Description

< Calibration data for used rejection filter (on camera, may not be present)

The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.38 CxPicBufAPI Class Reference

Static Public Member Functions

- static `size_t ImageTypeSize` (`ExImageDataType eDataType`)
Returns number of bytes occupied by given image data type.
- static `bool AllocPicBufFromHeap` (`SxPicBuf &pic`, `quint32 uiWidth`, `quint32 uiHeight`, `quint32 uiChannels`, `ExImageDataType eDataType`, `quint32 uiStride=0`, `quint32 uiColorBitDepth=0`, `ExDataStorageFormat eFormat=exdataStoragePixel`, `ExColorFilterArray eCfa=excfaNone`, `bool bResetData=false`, `uint uiGenTL=0`)
Allocates and initializes the `SxPicBuf` structure.
- static `bool AllocPicBufFromHeap` (`SxPicBuf &aPic`, `const SxPicBufInfo &aInfo`, `bool bResetData=false`)
Allocates and initializes the `SxPicBuf` using the initialized `SxPicBufInfo` structure. The allocation is preformed using the `malloc()` function from heap.
- static `bool AllocPicBufFromPool` (`TxMemPoolHandle hPool`, `QObject *pMemPoolClient`, `SxPicBuf &aPic`, `const SxPicBufInfo &aInfo`, `bool bResetData=false`)
Allocates the `aPic` from the memory pool given by `hPool` handle using the initialized `aInfo` structure.
- static `void FreePicBuf` (`SxPicBuf &pic`, `bool bReset=true`)
Deallocates the given image depending on `SxPicBuf::m_bDataOwner`.
- static `bool CopyPicBuf` (`SxPicBuf &dstPic`, `const SxPicBuf &srcPic`, `bool bIncludeMetadata=true`)
Copies the `srcPic` to `dstPic`. The destination image must be correctly preallocated otherwise the function fails.
- static `bool CopyRect` (`SxPicBuf &dstPic`, `quint32 dx`, `quint32 dy`, `const SxPicBuf &srcPic`, `quint32 sx`, `quint32 sy`, `quint32 w`, `quint32 h`)
Copies the rectangle from `srcPic` to `dstPic`.
- static `bool CheckPicValidData` (`const SxPicBuf &pic`)
Checks if the data in N-bit image are really N-bit. E.g. returns false when value 1234 is stored in 10-bit image, where max value is 1023.
- static `bool ConvertPicToRGB` (`const SxPicBuf &picSrc`, `SxPicBuf &picDst`)
Attemp to convert source image to RGB, depending on the src and dst format.
- static `bool ConvertGray8ToRGB32` (`const SxPicBuf &picMono8`, `SxPicBuf &picRgb32`)
Converts the source MONO8 or RAW8 picbuf to picbuf in RGB32 format (0x00RRGGBB).
- static `bool ConvertGray8ToRGB24` (`const SxPicBuf &picMono8`, `SxPicBuf &picRgb24`)
Converts the source MONO8 or RAW8 picbuf to picbuf in RGB24 format (0xRRGGBB).
- static `bool ConvertGray8ToRGB565` (`const SxPicBuf &picMono8`, `SxPicBuf &picRgb565`)
Converts the source MONO8 or RAW8 picbuf to picbuf in RGB565 format.
- static `bool ConvertGray16ToRGB32` (`const SxPicBuf &picMono16`, `SxPicBuf &picRgb32`)
Converts the source MONO16 or RAW16 picbuf to picbuf in RGB32 format (0x00RRGGBB).
- static `bool ConvertGray16ToRGB24` (`const SxPicBuf &picMono16`, `SxPicBuf &picRgb24`)
Converts the source MONO16 or RAW16 picbuf to picbuf in RGB24 format (0xRRGGBB).
- static `bool ConvertGray16ToGray8` (`const SxPicBuf &picMono16`, `SxPicBuf &picMono8`)
Converts the source MONO16 or RAW16 picbuf to MONO8 or RAW8 picbuf.
- static `bool ConvertGray16ToRGB565` (`const SxPicBuf &picMono16`, `SxPicBuf &picRgb565`)
Converts the source MONO16 or RAW16 picbuf to picbuf in RGB565 format.
- static `bool ConvertRGB24ToRGB32` (`const SxPicBuf &picRgb24`, `SxPicBuf &picRgb32`)
Converts the source RGB24 (0xRRGGBB) picbuf to picbuf in RGB32 format (0x00RRGGBB).
- static `bool ConvertRGB24ToRGB565` (`const SxPicBuf &picRgb24`, `SxPicBuf &picRgb565`)
Converts the source RGB24 (0xRRGGBB) picbuf to picbuf in RGB565 format.
- static `bool ConvertRGB32ToRGB24` (`const SxPicBuf &picRgb32`, `SxPicBuf &picRgb24`)
Converts the source RGB32 (0x00RRGGBB) picbuf to picbuf in RGB24 format (0xRRGGBB).
- static `bool ConvertRGB32ToRGB565` (`const SxPicBuf &picRgb32`, `SxPicBuf &picRgb565`)
Converts the source RGB32 (0x00RRGGBB) picbuf to picbuf in RGB565 format.
- static `bool ConvertRGBPlanarToRGB32` (`const SxPicBuf &picRgbPlanar`, `SxPicBuf &picRgb32`)
Converts the source RGB Planar picbuf to picbuf in RGB32 format (0x00RRGGBB).

- static bool [ConvertRGBPlanarToRGB24](#) (const [SxPicBuf](#) &picRgbPlanar, [SxPicBuf](#) &picRgb24)
Converts the source RGB Planar picbuf to picbuf in RGB24 format (0xRRGGBB).
- static bool [ConvertRGBPlanarToRGB565](#) (const [SxPicBuf](#) &picRgbPlanar, [SxPicBuf](#) &picRgb565)
Converts the source RGB Planar picbuf to picbuf in RGB565 format.
- static bool [ConvertFFTTToRGB](#) (const [SxPicBuf](#) &picRe, const [SxPicBuf](#) &picIm, [SxPicBuf](#) &picRGB, [SxPicBuf](#) *pTmp=NULL)
Converts the source FFT images to a picbuf in the format of RGB.
- static bool [ConvertFFTTToRGB32](#) (const [SxPicBuf](#) &picRe, const [SxPicBuf](#) &picIm, [SxPicBuf](#) &picRGB32, [SxPicBuf](#) *pTmp=NULL)
Converts the source FFT images to a picbuf in the format of RGB32.
- static bool [ConvertFFTTToRGB24](#) (const [SxPicBuf](#) &picRe, const [SxPicBuf](#) &picIm, [SxPicBuf](#) &picRGB24, [SxPicBuf](#) *pTmp=NULL)
Converts the source FFT images to a picbuf in the format of RGB24.
- static bool [ConvertFFTTToRGB565](#) (const [SxPicBuf](#) &picRe, const [SxPicBuf](#) &picIm, [SxPicBuf](#) &picRGB565, [SxPicBuf](#) *pTmp=NULL)
Converts the source FFT images to a picbuf in the format of RGB565.
- static bool [ConvertPicToDifBpc](#) (const [SxPicBuf](#) &picSrc, [SxPicBuf](#) &picDst, bool bStretchValues=false)
Converts the source picbuf to picbuf with uint32 or uin16 type and same number of components.
- static void [SwapRedAndBlue](#) ([SxPicBuf](#) &picRgb)
Swaps the red and green channels.
- static void [SwapAllChannels](#) ([SxPicBuf](#) &pic)
Swaps all channels (i.e. ARGB -> BGRA).
- static void [SwapBytes](#) ([SxPicBuf](#) &pic)
From little to big endian or back efficiently.
- static void [FillChannelWithValue](#) ([SxPicBuf](#) &pic, quint32 uiChannel, quint16 uiValue)
Sets color to one component.
- static void [ResetRGB32AlphaByte](#) ([SxPicBuf](#) &pic)
Sets alpha byte to 0xFF, as this is required by Qt in format QImage::Format_RGB32.
- static bool [ReverseRowOrder](#) ([SxPicBuf](#) &pic)
The function reverses the row order of the image.
- static bool [ReverseMemRowOrder](#) (quint8 *pData, uint uiHeight, uint uiStride)
The function reverses the row order of the memory block.
- static bool [ReverseColOrder](#) ([SxPicBuf](#) &pic)
The function reverses the order of columns of the image.
- static bool [DetectColorBitDepth](#) (const [SxPicBuf](#) &picBuf, uint &uiBitDepth)
The function detect the color bit depth of the image (per channel). Possible values that can be returned: 1, 8, 10, 12, 14, 16.
- static bool [Resize](#) (const [SxPicBuf](#) &srcPic, [SxPicBuf](#) &dstPic)
Resizes the srcPic to the size given by dstPic. Both images must have the same image format. Resizing done using bilinear algorithm.
- static bool [ResizeRect](#) (const [SxPicBuf](#) &srcPic, quint32 sx, quint32 sy, quint32 sw, quint32 sh, [SxPicBuf](#) &dstPic, quint32 dx, quint32 dy, quint32 dw, quint32 dh)
Resizes the rectangular section within srcPic to the rectangle within dstPic. Both images must have the same image format. Resizing done using bilinear algorithm.
- static bool [ShiftBits](#) ([SxPicBuf](#) &srcPic, int iShift)
Shift pixel bits right (iShift > 0) or left (iShift < 0).
- static bool [Subtract](#) ([SxPicBuf](#) &pic1, const [SxPicBuf](#) &pic2)
Subtracts values of pic2 from pic1, leaving the result in pic1.
- static bool [Add](#) ([SxPicBuf](#) &pic1, const [SxPicBuf](#) &picToAdd)
Adds values of picToAdd to pic1, leaving the result in pic1. Pic1 can be Uint32.
- static bool [DivideByValue](#) ([SxPicBuf](#) &picSrc, int iDivideBy)
Divides value in src picture by constant. picSrc can be Uint32.

- static bool [DivideByValue](#) (const [SxPicBuf](#) &picSrc, int iDivideBy, [SxPicBuf](#) &picDst)
Divides value in src picture by constant, storing result to dst. picSrc can be Uint32.
- static bool [ExtractChannel](#) (const [SxPicBuf](#) &picSrc, [SxPicBuf](#) &picDst, const quint32 uiComp)
Extracts one channel from source image. Dst image must be mono with same size as src.
- static bool [FilterMinOnMosaicPicBuf](#) (const [SxPicBuf](#) &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight, [SxPicBuf](#) &picDst)
Applies the "min" filter on the mosaic image (RAW or HSI).
- static bool [FilterMaxOnMosaicPicBuf](#) (const [SxPicBuf](#) &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight, [SxPicBuf](#) &picDst)
Applies the "max" filter on the mosaic image (RAW or HSI).
- static bool [FilterAvgOnMosaicPicBuf](#) (const [SxPicBuf](#) &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight, [SxPicBuf](#) &picDst)
Applies the "average" filter on the mosaic image (RAW or HSI).
- static bool [FilterNormalizeOnMosaicPicBuf](#) (const [SxPicBuf](#) &picSrc, const QRect &rcSrc, int iMosaicWidth, int iMosaicHeight, [SxPicBuf](#) &picDst)
Applies the "normalization" filter on the mosaic image (RAW or HSI). The output image must have the same size as input image.
- static bool [Median](#) (const [SxPicBuf](#) &picSrc, [SxPicBuf](#) &picDst, int iKernelRadius, void *pBuffer=NULL)
Performs a median filtering on the picSrc with given kernelRadius.
- static bool [MedianAllocBuffer](#) (const [SxPicBufInfo](#) &picSrc, int iKernelRadius, void **pBuffer)
Allocates a temporary buffer used by [Median\(\)](#) function.
- static void [MedianFreeBuffer](#) (void *pBuffer)
Frees the buffer allocated by [MedianAllocBuffer\(\)](#).

4.38.1 Member Function Documentation

- 4.38.1.1 static bool CxPicBufAPI::AllocPicBufFromHeap ([SxPicBuf](#) & pic, quint32 uiWidth, quint32 uiHeight, quint32 uiChannels, ExImageDataType eDataType, quint32 uiStride = 0, quint32 uiColorBitDepth = 0, ExDataStorageFormat eFormat = exdataStoragePixel, ExColorFilterArray eCfa = excfaNone, bool bResetData = false, uint uiGenTL = 0) [static]

Allocates and initializes the [SxPicBuf](#) structure.

Remarks

If the uiStride is zero then it will be calculated. If the bResetData is true then the allocated data will be set to zero. The allocation is preformed using the malloc() function from heap.

- 4.38.1.2 static bool CxPicBufAPI::AllocPicBufFromPool (TxMemPoolHandle hPool, QObject * pMemPoolClient, [SxPicBuf](#) & aPic, const [SxPicBufInfo](#) & aInfo, bool bResetData = false) [static]

Allocates the aPic from the memory pool given by hPool handle using the initialized aInfo structure.

Remarks

If the hPool is NULL, then the image will be allocated using [AllocPicBufFromHeap\(\)](#)

- 4.38.1.3 static bool CxPicBufAPI::ConvertFFToRGB (const [SxPicBuf](#) & picRe, const [SxPicBuf](#) & picIm, [SxPicBuf](#) & picRGB, [SxPicBuf](#) * pTmp = NULL) [static]

Converts the source FFT images to a picbuf in the format of RGB.

Remarks

In fact the function calculates an amplitued of the FFT and converts it to RGB. Internally it calls [ConvertFFToRGB32\(\)](#) or [ConvertFFToRGB24\(\)](#) or [ConvertFFToRGB565\(\)](#).

4.38.1.4 `static bool CxPicBufAPI::ConvertFFTTToRGB24 ( const SxPicBuf & picRe, const SxPicBuf & picIm, SxPicBuf & picRGB24, SxPicBuf * pTmp = NULL ) [static]`

Converts the source FFT images to a picbuf in the format of RGB24.

Remarks

In fact the function calculates an amplitued of the FFT and converts it to RGB24.

4.38.1.5 `static bool CxPicBufAPI::ConvertFFTTToRGB32 ( const SxPicBuf & picRe, const SxPicBuf & picIm, SxPicBuf & picRGB32, SxPicBuf * pTmp = NULL ) [static]`

Converts the source FFT images to a picbuf in the format of RGB32.

Remarks

In fact the function calculates an amplitued of the FFT and converts it to RGB32.

4.38.1.6 `static bool CxPicBufAPI::ConvertFFTTToRGB565 ( const SxPicBuf & picRe, const SxPicBuf & picIm, SxPicBuf & picRGB565, SxPicBuf * pTmp = NULL ) [static]`

Converts the source FFT images to a picbuf in the format of RGB565.

Remarks

In fact the function calculates an amplitued of the FFT and converts it to RGB565.

4.38.1.7 `static bool CxPicBufAPI::ConvertPicToDifBpc ( const SxPicBuf & picSrc, SxPicBuf & picDst, bool bStretchValues = false ) [static]`

Converts the source picbuf to picbuf with uint32 or uin16 type and same number of components.

Remarks

Suitable for averaging, see [Add](#) and [DivideByValue](#).

4.38.1.8 `static bool CxPicBufAPI::DetectColorBitDepth ( const SxPicBuf & picBuf, uint & uiBitDepth ) [static]`

The function detect the color bit depth of the image (per channel). Possible values that can be returned: 1, 8, 10, 12, 14, 16.

Remarks

The function find maximum value in each channel and from the total maximum will guess the color bit depth. The color bit depth will be detected on images with 1 or 2 bytes per channel only.

4.38.1.9 `static bool CxPicBufAPI::FilterAvgOnMosaicPicBuf ( const SxPicBuf & picSrc, const QRect & rcSrc, int iMosaicWidth, int iMosaicHeight, SxPicBuf & picDst ) [static]`

Applies the "average" filter on the mosaic image (RAW or HSI).

Remarks

The output image can have either the size of:

- [picSrc.m_uiWidth, picSrc.m_uiHeight] for interpolated image,
- [rcSrc.width()/iMosaicWidth, rcSrc.height()/iMosaicHeight]. For interpolated image the inactive area of the input image is copied to the output image.

4.38.1.10 `static bool CxPicBufAPI::FilterMaxOnMosaicPicBuf ( const SxPicBuf & picSrc, const QRect & rcSrc, int iMosaicWidth, int iMosaicHeight, SxPicBuf & picDst ) [static]`

Applies the "max" filter on the mosaic image (RAW or HSI).

Remarks

The output image can have either the size of:

- [picSrc.m_uiWidth, picSrc.m_uiHeight] for interpolated image,
- [rcSrc.width()/iMosaicWidth, rcSrc.height()/iMosaicHeight]. For interpolated image the inactive area of the input image is copied to the output image.

4.38.1.11 `static bool CxPicBufAPI::FilterMinOnMosaicPicBuf ( const SxPicBuf & picSrc, const QRect & rcSrc, int iMosaicWidth, int iMosaicHeight, SxPicBuf & picDst ) [static]`

Applies the "min" filter on the mosaic image (RAW or HSI).

Remarks

The output image can have either the size of:

- [picSrc.m_uiWidth, picSrc.m_uiHeight] for interpolated image,
- [rcSrc.width()/iMosaicWidth, rcSrc.height()/iMosaicHeight]. For interpolated image the inactive area of the input image is copied to the output image.

4.38.1.12 `static bool CxPicBufAPI::FilterNormalizeOnMosaicPicBuf ( const SxPicBuf & picSrc, const QRect & rcSrc, int iMosaicWidth, int iMosaicHeight, SxPicBuf & picDst ) [static]`

Applies the "normalization" filter on the mosaic image (RAW or HSI). The output image must have the same size as input image.

Remarks

The normalization is performed per hyperpixel i.e. the hyperpixel's intensities are normalized to $<0; \max \leftarrow \text{Intensity} >$. An eventual inactive area of the input image is copied to the output image.

4.38.1.13 `static void CxPicBufAPI::FreePicBuf ( SxPicBuf & pic, bool bReset = true ) [static]`

Deallocates the given image depending on SxPicBuf::m_bDataOwner.

Remarks

If the SxPicBuf::m_eDataOwner member is not [SxPicBuf::edoThisPicBuf](#) then the function tries to release the buffer calling [CxMemoryManager::freePicBuf\(\)](#). The member [SxPicBuf::m_pData](#) is always set to NULL.

4.38.1.14 `static bool CxPicBufAPI::Median ( const SxPicBuf & picSrc, SxPicBuf & picDst, int iKernelRadius, void * pBuffer = NULL ) [static]`

Performs a median filtering on the `picSrc` with given `kernelRadius`.

Remarks

The kernel size used by the algorithm is $(\text{kernelRadius} * 2 + 1)$. The function operates with 8,10,12,14 and 16bit RGB or MONO images. The mosaic images (HSI, RAW) are not supported.

Temporary buffer `pBuffer`

If the `pBuffer` is `NULL` then the function allocates and frees it internally. For batch processing it is recommended to prealloc it before the batch starts. For allocation of the buffer use [MedianAllocBuffer\(\)](#)

4.38.1.15 `static bool CxPicBufAPI::ReverseColOrder ( SxPicBuf & pic ) [static]`

The function reverses the order of columns of the image.

Warning

The function does not change the eventual [SxPicBuf::m_eColorFilterArray](#) for RAW images!!!

4.38.1.16 `static bool CxPicBufAPI::ReverseRowOrder ( SxPicBuf & pic ) [static]`

The function reverses the row order of the image.

Warning

The function does not change the eventual [SxPicBuf::m_eColorFilterArray](#) for RAW images!!!

The documentation for this class was generated from the following file:

- `PicBufAPI.h`

4.39 CxPicBufMeas Class Reference

Static Public Member Functions

- static bool **measureHistogram** ([CxHistoData](#) &rData, const [CxImageData](#) *pImageData, const QRect &rc, `CxHistoData::ExHistoDataFlags` aFlags)
- static bool **measureLineProfile** ([CxProfileData](#) &rData, const [CxImageData](#) *pImageData, const QLine &line, int iBandSize)

Static Private Member Functions

- static void **fillMeasChannels** (const [SxPicBuf](#) &rPic, `QVector< QColor >` &rvecChannelColor, `QVector< QString >` &rvecChannelName, bool bTreatRawAsMono)

The documentation for this class was generated from the following file:

- `PicBufMeas.h`

4.40 CxHistoData::CxPixelStats Class Reference

Class for measurement of data in the ROI/image per channel.

Public Attributes

- quint32 **m_uiMin**
- quint32 **m_uiMax**
min, max
- double **m_dMean**
- double **m_dStDev**
mean, standard deviation
- double **m_dRms**
*root mean square: $\sqrt{\text{sum}(x*x)/n}$*

4.40.1 Detailed Description

Class for measurement of data in the ROI/image per channel.


The documentation for this class was generated from the following file:

- PicBufMeas.h

4.41 CxPlugin Class Reference

Base class of all plugins. For more informations about plugin development see the plugins.

Inheritance diagram for CxPlugin:



Public Member Functions

- virtual bool **init** ()
Called on application start.
- virtual void **finish** ()
Called on application exit.
- virtual void **loadSettings** (QSettings *pSettings)
*Called right after **init**(), and then while loading program state.*
- virtual void **saveSettings** (QSettings *pSettings)
*Called before **finish**(), and when program state is saved.*
- virtual QString **copyrightNotice** () const
In case plugin needs to write some legal info to Help - About - Legal...
- virtual QString **name** () const
Plugin name displayed in the Plugin Manager view.
- virtual QString **author** () const
Name of the author (company) of this plugin.

- virtual QString [description](#) () const
Short description of the plugin, notes, warnings etc.
- virtual void [version](#) (int &major, int &minor, int &build) const
Version of the plugin.
- virtual QString [version](#) () const
Returns a version string. The default implementation returns a string formatted as M.m.b (i.e. Major.minor.build).
- virtual QString [website](#) () const
The author's website adress.
- virtual void [activeViewChanged](#) (int iViewId)
When new MDI view activated (or created)
- virtual void [viewAboutToClose](#) (int iViewId)
The view will closed and deleted immediatelly.
- virtual void [newCameraConnected](#) (XICORE_HANDLE hCamera)
After new camera connected after app start.
- virtual void [cameraDisconnected](#) (XICORE_HANDLE hCamera)
After new camera disconnected during app session (not when exiting)
- virtual void [viewGotNewImage](#) (int iViewId, const [CxImageData](#) *pImageData)
New image received when acqusition running.
- virtual bool [viewMouseButtonPressEvent](#) (int iViewId, QGraphicsSceneMouseEvent *event)
Return true to prevent further processing.
- virtual bool [viewMouseMoveEvent](#) (int iViewId, QGraphicsSceneMouseEvent *event)
Return true to prevent further processing.
- virtual bool [viewMouseButtonReleaseEvent](#) (int iViewId, QGraphicsSceneMouseEvent *event)
Return true to prevent further processing.
- virtual bool [viewMouseDoubleClickEvent](#) (int iViewId, QGraphicsSceneMouseEvent *event)
Return true to prevent further processing.
- virtual bool [customCommand](#) (const QString &sCmdName, int iParam1, int iParam2, void *iParam1, void *iParam2)
for internal use.

4.41.1 Detailed Description

Base class of all plugins. For more informations about plugin development see the plugins.

The documentation for this class was generated from the following file:

- [Plugin.h](#)

4.42 CxProfileData Class Reference

Class holding measured line profile data.

Public Member Functions

- void [freeChannelData](#) ()

Public Attributes

- quint32 [m_uiMaxVal](#)
max value according to bpc (i.e. 1024 for 10 bit images)
- int [m_iLength](#)
size of arrays in m_vecChannelData
- QVector< float * > [m_vecChannelData](#)
measured profile per channel (float because of averaging)
- QVector< QColor > [m_vecChannelColor](#)
suggested channel colors
- QVector< QString > [m_vecChannelName](#)
suggested channel names

4.42.1 Detailed Description

Class holding measured line profile data.

The documentation for this class was generated from the following file:

- PicBufMeas.h

4.43 CxLut::CxPseudoColorTable Class Reference

Public Member Functions

- bool **loadFromFile** (const QString &sFile)
- bool **saveToFile** (const QString &sFile)

Public Attributes

- QString **m_sName**
- quint32 **m_arrValues** [256]


The documentation for this class was generated from the following file:

- Lut.h

4.44 CxRgb32Convertor Class Reference

The class converts all incoming image data to RGB32 format."

Inheritance diagram for CxRgb32Convertor:



Public Types

- enum `ExConversionOutputFormat` { `ecofRGB32`, `ecofRGB24`, `ecofRGB565` }

Public Member Functions

- virtual `QString title ()` const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual bool `acceptsDataFrom (CxChainable *pPredecessor)`
Query if this object can accept data from the given chainable object.
- virtual int `buffersCountInMemoryPool ()` const
Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).
- virtual `CxChainable * clone ()`
Creates a deep copy of this object.
- virtual bool `queryOutputImageInfo (const SxPicBufInfo &picInfoInput, SxPicBufInfo &picInfoOutput, const CxImageMetadata *pMetadataInput, CxImageMetadata *pMetadataOutput=NULL)`
Query the output image info (format) for given input image format.
- void `setConversionOutputFormat (ExConversionOutputFormat eFmt)`
- `ExConversionOutputFormat conversionOutputFormat ()`
- void `setLut (CxLut *pLut)`
set conversion LUTs. Objects takes ownership of the pointer.
- `CxLut * lut ()`

Protected Member Functions

- virtual `lxChainData * processData (lxChainData *pReceivedData)`
The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.
- virtual void `setupStop ()`
Sets up the object after the object is stopped.
- bool `createConvertedPicBuf (const CxImageData *pSrcData, SxPicBuf *pDstPic, CxImageMetadata *pMetadataOutput)`

Protected Attributes

- ExConversionOutputFormat **m_eConvOutputFormat**
- [CxLut](#) * **m_pLut**
- QMutex **m_lockLuts**
- [SxPicBuf](#) **m_picFFTtmp**

Preallocated temporary buffer for FFT conversion to RGB.

Additional Inherited Members

4.44.1 Detailed Description

The class converts all incoming image data to RGB32 format."

4.44.2 Member Function Documentation

4.44.2.1 `virtual IxChainData* CxRgb32Converter::processData ( IxChainData * pReceivedData ) [protected], [virtual]`

The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Remarks

If the implementation creates (allocates) a new data here then the original (received) data should be deleted using the delete operator. The default implementation of this method returns a pointer to the received data without any processing.

Warning

If the method returns NULL, then the `pReceivedData` will be deleted automatically!

Reimplemented from [CxChainable](#).

Reimplemented in [CxViewImageReceiver](#).

4.44.2.2 `virtual void CxRgb32Converter::setupStop ( ) [protected], [virtual]`

Sets up the object after the object is stopped.

Remarks

The method is internally used by the [stop\(\)](#) method immediately after the [run\(\)](#) method is finished. For example, the [CxCameraSource](#) stops camera acquisition here. The default implementation does nothing.

Reimplemented from [CxChainable](#).

4.44.3 Member Data Documentation

4.44.3.1 `QMutex CxRgb32Converter::m_lockLuts [protected]`

Lock when working with LUTs.


The documentation for this class was generated from the following file:

- [Rgb32Converter.h](#)

4.45 CxRTTI Class Reference

Our RTTI builds upon Qt meta object system, and just adds class discovery.

Inheritance diagram for CxRTTI:



Signals

- void [itemEnableStateChanged](#) (TxRttiItem hItem, bool bEnabled)
The signal is emitted from [setItemEnabled\(\)](#) method every time the item's enable state changes.

Public Member Functions

- void [registerClass](#) (CxPlugin *pPlugin, const QMetaObject &aClassMetaObj, const QString &sInterface)
This method is called from xiRTTI_REGISTER macro, do not use directly.
- QStringList [classesForInterface](#) (const QString &sInterface, bool bEnabledOnly=true)
Queries available classes - get all classes that implement given interface.
- QString [interfaceForClass](#) (const QString &sClass)
Returns an interface for the given registered class name. It returns empty string when the class is not registered.
- QString [customNameForClass](#) (const QString &sClass)
Title that should be presented to user when enumerating the class.
- QString [uniqueIdForClass](#) (const QString &sClass)
UniqueId of the class (it may not exist and the method may return empty string in that case).
- QString [classForUniqueId](#) (const QString &sUniqueId)
Returns a classname for given UniqueID of chainable (if it exists).
- template<class T >
T [createInstanceAndCast](#) (const QString &sClass)
- bool [itemEnabled](#) (TxRttiItem hItem)
Checks whether the corresponding Rtti item is enabled (used by plugin manager).
- void [setItemEnabled](#) (TxRttiItem hItem, bool bEnable)
Sets the corresponding Rtti item enabled/disabled.
- QString [itemName](#) (TxRttiItem hItem)
Gets the "CustomName" of the corresponding Rtti item from metadata of the object.
- QMap< QString, QList< TxRttiItem > > [itemsByInterface](#) ()
Returns Rtti items grouped by interface.
- const CxPlugin * [pluginForItem](#) (TxRttiItem hItem) const
Returns a plugin from which this item was registered.
- QList< TxRttiItem > [itemsForPlugin](#) (const CxPlugin *pPlugin) const
Returns a list of rtti items registered by the plugin.
- QString [interfaceForItem](#) (TxRttiItem hItem)
Returns the interface which this item implements.
- void [saveSettings](#) (QSettings *pSettings)
Saves the current state of the Rtti items (their enable state etc.)
- void [loadSettings](#) (const QSettings *pSettings)
Loads the state of Rtti items.

Static Public Member Functions

- static [CxRTTI](#) * **instance** ()

Protected Member Functions

- [QObject](#) * **createInstanceOfClass** (const [QString](#) &sClass)
Helper function. Use createInstanceAndCast instead.
- [QString](#) **defaultNameForInterfaceClass** (const [QString](#) &interfaceClass)
If the Rtti item class has not the CustomName then this function generates a default name.
- [QString](#) **nameFromMetadataForClass** (const [QString](#) &sClass, const [QString](#) &sMetadataTag)
Returns the name from metadata for given class and metadata tag.

4.45.1 Detailed Description

Our RTTI builds upon Qt meta object system, and just adds class discovery.

Registered interfaces are: [CxChainable](#), [IxImageDataLoader](#), [IxImageDataSaver](#), [IxVideoEncoder](#), [IxVideoDecoder](#)

4.45.2 Member Function Documentation

4.45.2.1 `template<class T > T CxRTTI::createInstanceAndCast ( const QString & sClass ) [inline]`

Creates instance of the class and casts for requested interface. Important: the default constructor of that class must be declared as `Q_INVOKABLE`!


The documentation for this class was generated from the following file:

- `xiRtti.h`

4.46 CxSegmentedControl Class Reference

Set of exclusively checkable buttons next to each other, inspired by OS X.

Inheritance diagram for CxSegmentedControl:



Signals

- void **selectedSegmentChanged** ()

Public Member Functions

- **CxSegmentedControl** ([QWidget](#) *parent=0)
- void **setSegmentCount** (int nCount)
- void **setSelectedSegment** (int idx)
- int **selectedIndex** ()

- void **setSegmentText** (int ildx, const QString &sText)
- void **setSegmentIcon** (int ildx, const QIcon &rlcon, int iWidth=-1, int iHeight=-1)

Protected Member Functions

- QToolButton * **segmentBtn** (int ildx)
- virtual bool **eventFilter** (QObject *pObject, QEvent *pEvent)

Protected Attributes

- int **m_nSegmentCount**

Private Slots

- void **buttonPressed** ()

4.46.1 Detailed Description

Set of exclusively checkable buttons next to each other, inspired by OS X.

The documentation for this class was generated from the following file:

- SegmentedControl.h

4.47 CxShadingPreset Class Reference

Classes

- class [CxItem](#)

Public Member Functions

- void [saveCurrentValuesFrom](#) (CxCameraParameters *pPars)
Save values of defining camera parameters.
- void [saveShadingImages](#) (CxCameraParameters *pPars, SxPicBuf *pPicDark, SxPicBuf *pPicMid, const QString &sStorageFolder)
Save shading images.
- void [removeShadingImages](#) ()
Delete stored images.
- bool [hasSameValuesAs](#) (CxCameraParameters *pPars) const
Are camera parameters same?
- void [activateValuesTo](#) (CxCameraParameters *pPars) const
Activate shading correction.

Static Public Member Functions

- static QStringList [definingParamList](#) ()
List of defining camera parameters (exposure, downsampling, etc...)
- static bool [isShadingActive](#) (CxCameraParameters *pPars)
Returns true if shading is switched on on the camera (TODO: remove after parameter exists in xiApi)
- static void [disableShading](#) (CxCameraParameters *pPars)
Switches off shading (TODO: remove after parameter exists in xiApi)

Public Attributes

- [QString m_sName](#)
Name of the preset.
- [QList< CxItem > m_1stParamValues](#)
Important camera parameters (defined in [definingParamList\(\)](#))
- [QRect m_rcRoi](#)
ROI is always stored.
- [QString m_sDarkFilename](#)
File with dark frame for shaindg correction.
- [QString m_sMidFilename](#)
File with light frame for shaindg correction.


The documentation for this class was generated from the following file:

- CameraParameters.h

4.48 CxStaticImageSource Class Reference

The image source of static image.

Inheritance diagram for CxStaticImageSource:



Public Member Functions

- [CxStaticImageSource](#) (const [QString](#) &sPath)
Creates the instance from the existing file name.
- [CxStaticImageSource](#) (const [SxPicBuf](#) &picBuf, const [QString](#) &sName)
Creates the instance from picbuf.
- [CxStaticImageSource](#) ([CxImageData](#) *pImageData)
Creates the instance from image data.
- [~CxStaticImageSource](#) ()
Destructor.
- virtual [QString title](#) () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual int [buffersCountInMemoryPool](#) () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*

- virtual [CxChainable](#) * [clone](#) ()
Creates a deep copy of this object.
- virtual [IxChainData](#) * [data](#) (QObject *pReceiver, qint32 iFrameNo)
Gets the pointer to data processed by this object.
- virtual [IxChainData](#) * [latestData](#) (QObject *pReceiver)
Gets the latest data from the output queue. It is called automatically if `ExChainableFlag::eTakeLatestImage` is set.
- virtual bool [saveSettings](#) (QDomDocument &aDoc, QDomElement &aElm) const
The method saves the settings of this object to the DOM element.
- virtual bool [loadSettings](#) (const QDomElement &aDom)
Loads settings from DOM.
- virtual bool [currentOutputImageInfo](#) (SxPicBufInfo &picInfo, [CxImageMetadata](#) *pMetadata=NULL)
Gets the current output image info (format) provided by this object.
- void [run](#) ()
This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the `CHAINABLE_RUN_WHILE_TRUE` and `CHAINABLE_RUN_END_WHILE` macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).
- bool [isLoading](#) () const
Whether the static image is loaded correctly.
- bool [reload](#) ()
Reloads the image.
- QString [path](#) ()
Get the image full path.
- void [setPath](#) (const QString &path)
Sets the image's path. The method will not reload the image when the path differs from the original path.

Protected Member Functions

- bool [loadFromFile](#) (const QString &sPath)

Private Attributes

- [CxImageData](#) * [m_pData](#)
Pointer to image data.
- QString [m_sPath](#)
Path of file.

Additional Inherited Members

4.48.1 Detailed Description

The image source of static image.

4.48.2 Member Function Documentation

- 4.48.2.1 virtual bool [CxStaticImageSource::currentOutputImageInfo](#) ([SxPicBufInfo](#) & *picInfo*, [CxImageMetadata](#) * *pMetadata* = NULL) [virtual]

Gets the current output image info (format) provided by this object.

Remarks

The method should call `getCurrentOutputImageInfo()` on its predecessors and calculate the output image info provided by this object.

Reimplemented from [CxImageProvider](#).

4.48.2.2 `virtual lxChainData* CxStaticImageSource::data ( QObject * pReceiver, qint32 iFrameNo ) [virtual]`

Gets the pointer to data processed by this object.

See also

[discardData\(\)](#)
[CxChainable](#)

Reimplemented from [CxChainable](#).

4.48.2.3 `void CxStaticImageSource::run ( ) [virtual]`

This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).

See also

[CxChainable](#)

Reimplemented from [CxChainable](#).

4.48.2.4 `virtual bool CxStaticImageSource::saveSettings ( QDomDocument & aDoc, QDomElement & aElm ) const [virtual]`

The method saves the settings of this object to the DOM element.

Remarks

The element `aElm` is already named by a caller, do not change the name (ie. do not call the `aElm.setTag←Name()`). You can add attributes and child elements.

Reimplemented from [CxChainable](#).


The documentation for this class was generated from the following file:

- `StaticImageSource.h`

4.49 CxTuningBar Class Reference

Our own implementation of slides, working with [CxValue](#) types `etInt` and `etFloat`.

Inheritance diagram for `CxTuningBar`:



Public Types

- enum **ExBarColorScale** {
 ebcsNone, **ebcsGray**, **ebcsRed**, **ebcsGreen**,
 ebcsBlue, **ebcsColorful** }
- enum **ExBarTransform** { **ebtLinear**, **ebtLogarithmic** }

Public Slots

- void **setValue** (const [CxValue](#) &vValue)
- void **setReadOnly** (bool bReadOnly)
- void **processExternalWheelEvent** (QWheelEvent *event)

Signals

- void **valueChanged** (const [CxValue](#) &vValue)

Public Member Functions

- **CxTuningBar** (QWidget *parent=0)
- void **setRange** (const [CxValue](#) &vMin, const [CxValue](#) &vMax, const [CxValue](#) &vIncrement)
- const [CxValue](#) & **rangeMin** () const
- const [CxValue](#) & **rangeMax** () const
- const [CxValue](#) & **increment** () const
- const [CxValue](#) & **value** () const

Public Attributes

- bool **m_bShowTickMarks**
- bool **m_bSnapToTickMarks**
- ExBarColorScale **m_eColorScale**
- ExBarTransform **m_eTransform**

Protected Member Functions

- int **calcPaintXPos** (const [CxValue](#) &vValue)
- [CxValue](#) **calcValueFromPos** (int nXPos)
- bool **setValueAndEmitSignal** ([CxValue](#) vNewVal)
- [CxValue](#) **smallStep** (int iRangeFraction, bool bUp)
- virtual QSize **sizeHint** () const
- virtual void **paintEvent** (QPaintEvent *)
- virtual void **mousePressEvent** (QMouseEvent *event)
- virtual void **mouseMoveEvent** (QMouseEvent *event)
- virtual void **mouseReleaseEvent** (QMouseEvent *event)
- virtual void **wheelEvent** (QWheelEvent *event)
- virtual void **keyPressEvent** (QKeyEvent *event)

Protected Attributes

- [CxValue](#) **m_vRangeMin**
- [CxValue](#) **m_vRangeMax**
- [CxValue](#) **m_vIncrement**
- [CxValue](#) **m_vValue**
- **bool m_bReadOnly**
- **bool m_bMovingArrow**
- **int m_iBorder_X**
- **int m_iSnapDist**

Static Protected Attributes

- static QPixmap * **m_pPicArrow**
- static QPixmap * **m_pPicArrowDis**
- static QPixmap * **m_pPicBarColorful**

4.49.1 Detailed Description

Our own implementation of slides, working with [CxValue](#) types `etInt` and `etFloat`.

The documentation for this class was generated from the following file:

- `TuningBar.h`

4.50 CxUnits Class Reference

Class the defines all known units and measures, and includes conversions between them.

Public Types

- enum [ExUnitGroup](#) {
 eugInvalid, **eugTime**, **eugLength**, **eugTemperature**,
 eugRatio }
 All measures defined here.
- enum **ExTime** {
 etMicrosecond = 0, **etMillisecond** = 1, **etSecond** = 2, **etMinute** = 3,
 etHour = 4 }
- enum **ExLength** { **etMicrometre** = 0, **etMillimetre** = 1, **etCentimetre** = 2, **etMetre** = 3 }
- enum **ExTemperature** { **etCelsius** = 0, **etFahrenheit** = 1, **etKelvin** = 2 }
- enum [ExRatio](#) { **erFraction** = 0, **erPercent** = 1 }

Static Public Member Functions

- static int **countInGroup** ([ExUnitGroup](#) eUnitGroup)
 Number of units for given group.
- static QString **text** ([ExUnitGroup](#) eUnitGroup, int iUnit)
 Text for the unit.
- static int **unitFromText** ([ExUnitGroup](#) eUnitGroup, const QString &sText)
 Unit from text, returns -1 when not found.
- static double **convertValue** (double dValue, [ExUnitGroup](#) eUnitGroup, int iUnitFrom, int iUnitTo)

4.50.1 Detailed Description

Class the defines all known units and measures, and includes conversions between them.

4.50.2 Member Enumeration Documentation

4.50.2.1 enum CxUnits::ExRatio

Enumerator

- erFraction** Ratio as decimal value between 0 and 1.
- erPercent** Ratio as percent value between 0 and 100.

The documentation for this class was generated from the following file:

- Units.h

4.51 CxValue Class Reference

Class for storing values similar to variant. Used in camera parameters.

Public Types

- enum [Type](#) {
[etInvalid](#), [etInt](#), [etFloat](#), [etBool](#),
[etString](#) }

Public Member Functions

- **CxValue** (int nVal)
- **CxValue** (float fVal)
- **CxValue** (bool bVal)
- **CxValue** (const QString &sVal)
- bool **operator==** (const [CxValue](#) &rOther) const
- bool **operator!=** (const [CxValue](#) &rOther) const
- bool **operator<** (const [CxValue](#) &rOther) const
- bool **operator>** (const [CxValue](#) &rOther) const
- [CxValue](#) **operator+** (const [CxValue](#) &rOther) const
- [CxValue](#) **operator-** (const [CxValue](#) &rOther) const
- [CxValue](#) **operator*** (const [CxValue](#) &rOther) const
- [CxValue](#) **operator/** (const [CxValue](#) &rOther) const
- int **toInt** () const
Returns value as int (not implemented for etString)
- float **toFloat** () const
Returns value as float (for etInt and etFloat)
- bool **toBool** () const
Returns value as bool (value is 1)
- bool **fromString** (const QString &sString)
Loading from string using the type already in m_eType.
- QString **toString** () const
Save value to string.

Public Attributes

- enum [CxValue::Type m_eType](#)
Value type.
- union {
 int **m_intVal**
 float **m_floatVal**
 bool **m_boolVal**
};
- QString **m_stringVal**

4.51.1 Detailed Description

Class for storing values similar to variant. Used in camera parameters.

4.51.2 Member Enumeration Documentation

4.51.2.1 enum CxValue::Type

Enumerator

- etInvalid** Marks an error, or undefined value.
- etInt** Use m_intVal.
- etFloat** Use m_floatVal.
- etBool** Use m_boolVal.
- etString** Use m_stringVal.

The documentation for this class was generated from the following file:

- ValueTypes.h

4.52 CxValueRegInvalidator Class Reference

Class describing invalidator connection.

Public Member Functions

- **CxValueRegInvalidator** (const QString &sXiApiName, bool bReSet=false)
- **CxValueRegInvalidator** (const char *szXiApiName)

Public Attributes

- QString [m_sXiApiName](#)
name of the parameter that invalidates this register
- bool [m_bReSet](#)
when set to true, we should set the read value back to API again

4.52.1 Detailed Description

Class describing invalidator connection.

The documentation for this class was generated from the following file:

- CameraParameters.h

4.53 CxValueRegister Class Reference

This class is responsible for the communitaion with the camera.

Public Member Functions

- bool **readValue** (XICORE_HANDLE hCamera, [CxValue](#) &vValue) const
Note: vValue.m_eType must be correctly filler.
- bool **storeValue** (XICORE_HANDLE hCamera, const [CxValue](#) &vValue) const
- bool **isInvalidatedByParam** (const QString &sInvalidator) const
- bool **readIfInvalidated** (XICORE_HANDLE hCamera, [CxValue](#) &vValue, const QString &sInvalidator) const

Static Public Member Functions

- static bool **readValue** (XICORE_HANDLE hCamera, const QString &sXiApiParamName, [CxValue](#) &vValue)
- static bool **storeValue** (XICORE_HANDLE hCamera, const QString &sXiApiParamName, const [CxValue](#) &vValue, int *piRetVal=NULL)

Public Attributes

- QString **m_sXiApiName**
which parameter should we call (read from XML, e.g. "exposure" or "exposure:min". If empty, no update possible
- TxValueRegInvalidatorList **m_lstInvalidators**
paramaters that invalidate this register when they change

4.53.1 Detailed Description

This class is responsible for the communitaion with the camera.

The documentation for this class was generated from the following file:

- CameraParameters.h

4.54 CxVecObjContainer Class Reference

Our special class above QGraphicsScene maintains list of our vector objects and assigns IDs to them.

Public Member Functions

- **CxVecObjContainer** (QGraphicsScene *pScene)
- QRectF **sceneRect** () const
- TxVecObjID **addObject** ([CxVectorObject](#) *pObj)

- *Add object to scene, and assign it a new ID.*
- void **deleteObject** (TxVecObjID idObj)
Removes the object from scene and deletes the object.
- void **deleteAllObjects** ()
Removes all our vector object from scene.
- QVector< TxVecObjID > **listAllObjects** ()
Returns the list of all objects in the scene.
- CxVectorObject * **objectWithId** (TxVecObjID idObj)
Returns pointer from ID.
- void **scaleObjects** (double sx, double sy, bool bReflectConstraints=true)
- void **moveObjects** (double dx, double dy, bool bReflectConstraints=true)
- void **saveObjVisibilityState** ()
Save the per-object vislity for later.
- void **restoreObjVisibilityState** ()
Restores the state as saved in saveObjVisibilityState.
- void **hideAllObjects** ()
Sets all vector objects as invisible.

Protected Attributes

- QGraphicsScene * **m_pScene**
Each container has to be connected to a scene.
- TxVecObjID **m_idMax**
Object ID generator.
- TxVectorObjectMap **m_mapObjs**
Mapping from object id to the pointer.
- QMap< TxVecObjID, bool > **m_mapObjVisibilityState**
Saved visibility state using saveObjVisibilityState.

4.54.1 Detailed Description

Our special class above QGraphicsScene maintains list of our vector objects and assigns IDs to them.


The documentation for this class was generated from the following file:

- VectorContainer.h

4.55 CxVectorObject Class Reference

Base class for all our vector objects in QGraphicsScene. Needed to send signals away.

Inheritance diagram for CxVectorObject:



Public Types

- enum [ExCameraFormatChangeBehavior](#) { [excfcAttachSensor](#), [excfcAttachScene](#) }

Signals

- void [customContextMenuRequested](#) (const QPoint &pos)
Sent after right mouse click on the object.
- void **doubleClicked** ()
- void **clicked** (int idObject)

Public Member Functions

- **CxVectorObject** (QGraphicsItem *parent=0)
- TxVecObjID **id** ()
- void **setId** (TxVecObjID idObj)
- virtual void **setPenColor** (const QColor &color)
- virtual void **setPenWidth** (int iWidth)
- virtual void [setSelectedPenWidth](#) (int iWidth)
Set pen width of selected object.
- virtual void **setBrushColor** (const QColor &color)
- virtual void [setCoordIncrement](#) (int nXInc, int nYInc)
0 for no limit, 1 for rounding to nearest int, 4 for coords in multiples of 4
- virtual void [setResizable](#) (bool bEnable)
Enable object resizing with mouse.
- virtual void [setKeepInScene](#) (bool bKeep)
Forces the object to stay within image bounds.
- virtual void [scaleObject](#) (double sx, double sy, bool bReflectConstraints=true)
Object scaling.
- virtual void [moveObject](#) (double dx, double dy, bool bReflectConstraints=true)
Object translation.
- QColor **penColor** ()
- QColor **brushColor** ()
- int **penWidth** ()
- int [selectedPenWidth](#) ()
Width of a pen of selected object.
- [ExCameraFormatChangeBehavior](#) **cameraFormatChangeBehavior** ()
See [ExCameraFormatChangeBehavior](#).
- virtual void [setCameraFormatChangeBehavior](#) ([ExCameraFormatChangeBehavior](#) eMode)
See [ExCameraFormatChangeBehavior](#).

Protected Member Functions

- QPen **defaultPen** () const
- QBrush **defaultBrush** () const
- bool **reflectCoordIncrement** (QPointF &pt) const
- qreal **viewScaleAtEvent** (QGraphicsSceneEvent *pEvent)
- virtual void **mouseDoubleClickEvent** (QGraphicsSceneMouseEvent *event)
- virtual void **mousePressEvent** (QGraphicsSceneMouseEvent *event)

Protected Attributes

- TxVecObjID [m_id](#)
ID as set by [CxVecObjContainer](#).
- QColor [m_clPen](#)
Color of the pen used when painting the object.
- QColor [m_clBrush](#)
Fill color.
- int [m_iPenWidth](#)
- int [m_iSelectedPenWidth](#)
Pen width of selected object.
- int [m_nCoordXIncrement](#)
0 for no limit, 1 for rounding to nearest int, 4 for coords in multiples of 4
- int [m_nCoordYIncrement](#)
- bool [m_bIsResizable](#)
It is possible to resize the object (rectangle, line end points, ...)
- bool [m_bKeepInScene](#)
Do not allow to move the object outside the image boundaries.
- [ExCameraFormatChangeBehavior](#) [m_eCameraFormatChangeBehavior](#)
Defaults to [excfcbAttachSensor](#).

4.55.1 Detailed Description

Base class for all our vector objects in QGraphicsScene. Needed to send signals away.

4.55.2 Member Enumeration Documentation

4.55.2.1 enum CxVectorObject::ExCameraFormatChangeBehavior

Enumerator

- [excfcbAttachSensor](#)** Keep same position on chip relative to ROI and binning (default)
[excfcbAttachScene](#) Stay visually on same place in the image view.


The documentation for this class was generated from the following file:

- VectorObject.h

4.56 CxVectorObjectEllipse Class Reference

Object representing the ellipse, using the functionality implemented in rectangle.

Inheritance diagram for CxVectorObjectEllipse:



Public Member Functions

- **CxVectorObjectEllipse** ([CxVectorObject](#) *parent=0)
- virtual void **paint** (QPainter *painter, const QStyleOptionGraphicsItem *, QWidget *)

Protected Member Functions

- virtual QRectHitTest **hitTest** (qreal x, qreal y, qreal dViewScale)

Additional Inherited Members

4.56.1 Detailed Description

Object representing the ellipse, using the functionality implemented in rectangle.


The documentation for this class was generated from the following file:

- VectorObject.h

4.57 CxVectorObjectLine Class Reference

Object representing the single line, from start point to end point.

Inheritance diagram for CxVectorObjectLine:



Public Types

- enum [ExLineStyle](#) { **elsFree**, **elsHorizontal**, **elsVertical** }

Signals

- void [lineChanging](#) (const QLineF &[line](#))
Object is being changed (sent during mouse moves)
- void [lineChanged](#) (const QLineF &[line](#))
Object has new geometry (sent after mouse up)

Public Member Functions

- **CxVectorObjectLine** ([CxVectorObject](#) *parent=0)
- QLineF [line](#) ()
Returns line geometry.
- void [setLine](#) (const QLineF &[line](#))
Sets line geometry.

- virtual void [setCoordIncrement](#) (int nXInc, int nYInc)
Reimplemented to check both end points.
- virtual void [scaleObject](#) (double sx, double sy, bool bReflectConstraints=true)
Object scaling.
- virtual void [moveObject](#) (double dx, double dy, bool bReflectConstraints=true)
Object translation.
- [ExLineStyle](#) [lineStyle](#) ()
- void [setLineStyle](#) ([ExLineStyle](#) eStyle)
- virtual void [paint](#) (QPainter *painter, const QStyleOptionGraphicsItem *, QWidget *)

Protected Types

- enum [ExLineHitTest](#) { [ehtNone](#), [ehtEdge](#), [ehtStartPoint](#), [ehtEndPoint](#) }

Protected Member Functions

- double [distanceFromLine](#) (const QPointF &pt) const
- virtual [ExLineHitTest](#) [hitTest](#) (qreal x, qreal y, qreal dViewScale)
- virtual QRectF [boundingRect](#) () const
- virtual void [hoverMoveEvent](#) (QGraphicsSceneHoverEvent *event)
- virtual void [mousePressEvent](#) (QGraphicsSceneMouseEvent *event)
- virtual void [mouseMoveEvent](#) (QGraphicsSceneMouseEvent *event)
- virtual void [mouseReleaseEvent](#) (QGraphicsSceneMouseEvent *event)

Protected Attributes

- QPointF [m_ptEnd](#)
- [ExLineStyle](#) [m_eLineStyle](#)
- enum [CxVectorObjectLine::ExLineHitTest](#) [m_eLastHitTest](#)
- bool [m_bResizing](#)
- QPointF [m_ptStartDragMousePos](#)
- QPointF [m_ptStartPos](#)

4.57.1 Detailed Description

Object representing the single line, from start point to end point.

4.57.2 Member Enumeration Documentation

4.57.2.1 enum [CxVectorObjectLine::ExLineStyle](#)

Enumerator

- **[elsHorizontal](#)** Force line to be horizontal.
- **[elsVertical](#)** Force line to be vertical.


The documentation for this class was generated from the following file:

- [VectorObject.h](#)

4.58 CxVectorObjectRect Class Reference

Object representing the rectangular object. Implements mouse resizing and caption.

Inheritance diagram for CxVectorObjectRect:



Public Types

- enum **ExRectShape** { **ersSharpCorners**, **ersRoundCorners** }

Signals

- void **rectChanging** (const QRectF &rc)
Object is being changed (sent during mouse moves)
- void **rectChanged** (const QRectF &rc)
Object has new geometry (sent after mouse up)

Public Member Functions

- **CxVectorObjectRect** (**CxVectorObject** *parent=0)
- QRectF **rect** ()
- void **setRect** (const QRectF &rect)
- void **setRect** (qreal x, qreal y, qreal w, qreal h)
- qreal **width** () const
- qreal **height** () const
- void **setWidth** (qreal dWidth)
- void **setHeight** (qreal dHeight)
- void **setSizeIncrement** (int nWidthInc, int nHeightInc)
0 for no limit, 1 for rounding to nearest int, 4 for coords in multiples of 4
- void **setMinSize** (qreal dMinWidth, qreal dMinHeight)
Minimal size of the rectangle.
- QString **caption** ()
Caption is text drawn in the center of rectangle.
- void **setCaption** (const QString &sCaption)
- QColor **captionColor** ()
- void **setCaptionColor** (const QColor &color)
- virtual void **scaleObject** (double sx, double sy, bool bReflectConstraints=true)
Object scaling.
- virtual void **moveObject** (double dx, double dy, bool bReflectConstraints=true)
Object translation.
- virtual void **paint** (QPainter *painter, const QStyleOptionGraphicsItem *, QWidget *)
- virtual QRectF **boundingRect** () const

- ExRectShape **rectShape** () const
- void **setRectShape** (ExRectShape eShape)
- double **cornerRadiusX** () const
- double **cornerRadiusY** () const
- Qt::SizeMode **cornerSizeMode** () const
- void **setCornersRadius** (double dXradius, double dYradius, Qt::SizeMode mode=Qt::AbsoluteSize)
- void **setCornersRadius** (double dRadius, Qt::SizeMode mode=Qt::AbsoluteSize)

Protected Types

- enum **ExRectHitTest** {
ehtNone, **ehtMove**, **ehtTopEdge**, **ehtLeftEdge**,
ehtRightEdge, **ehtBottomEdge**, **ehtTLCorner**, **ehtTRCorner**,
ehtBLCorner, **ehtBRCorner** }

Protected Member Functions

- virtual ExRectHitTest **hitTest** (qreal x, qreal y, qreal dViewScale)
- bool **reflectSizeIncrement** (qreal &w, qreal &h) const
- virtual void **hoverMoveEvent** (QGraphicsSceneHoverEvent *event)
- virtual void **mousePressEvent** (QGraphicsSceneMouseEvent *event)
- virtual void **mouseMoveEvent** (QGraphicsSceneMouseEvent *event)
- virtual void **mouseReleaseEvent** (QGraphicsSceneMouseEvent *event)

Protected Attributes

- ExRectShape **m_eShape**
- double **m_dXradius**
- double **m_dYradius**
- Qt::SizeMode **m_eCornerSizeMode**
- qreal **m_dWidth**
- qreal **m_dHeight**
- enum CxVectorObjectRect::ExRectHitTest **m_eLastHitTest**
- int **m_nWidthIncrement**
- int **m_nHeightIncrement**
- QSizeF **m_sizeMin**
- QString **m_sCaption**
- QColor **m_clCaption**
- bool **m_bResizing**
- QPointF **m_ptStartDragMousePos**
- QPointF **m_ptStartRectPos**

4.58.1 Detailed Description

Object representing the rectangular object. Implements mouse resizing and caption.


The documentation for this class was generated from the following file:

- VectorObject.h

4.59 CxVideoCodecManager Class Reference

The class that manages video codecs used by video loaders and savers.

Inheritance diagram for CxVideoCodecManager:



Classes

- struct [SxCodecItem](#)

Public Member Functions

- [IxVideoEncoder](#) * **encoderForFourCC** (TxFourCC fcc)
- [IxVideoDecoder](#) * **decoderForFourCC** (TxFourCC fcc)
- QList< QPair< TxFourCC, QString > > **encoders** ()
Returns a list of encoders with their FourCCs and names.
- QList< QPair< TxFourCC, QString > > **decoders** ()
Returns a list of decoders with their FourCCs and names.
- void **setCodecOptions** (TxFourCC fccCodec, const TxImageFormatOptions &aOpt)
set defaults codec
- void **saveSettings** (QSettings *pSettings)
store format options to app settings
- void **loadSettings** (QSettings *pSettings)
load format options from app settings

Static Public Member Functions

- static [CxVideoCodecManager](#) * **instance** ()
Singleton class, this is the way to obtain the pointer.

Private Slots

- void **onRttiltemEnableStateChanged** (TxRttiltem hItem, bool bEnabled)
Handle the [CxRTTI](#) changes.

Private Member Functions

- void **resetEncoders** ()
- void **resetDecoders** ()
- void **loadAvailableEncoders** ()
- void **loadAvailableDecoders** ()

Static Private Member Functions

- static QString **findClassForFourCC** (const QList< [SxCodecItem](#) > &lstCodecs, TxFourCC fourCC)

Private Attributes

- bool **m_bEncodersLoaded**
- bool **m_bDecodersLoaded**
- QList< [SxCodecItem](#) > **m_lstEncoders**
- QList< [SxCodecItem](#) > **m_lstDecoders**
- QMap< TxFourCC, TxImageFormatOptions > **m_mapCodecOptions**
maps FourCC of a codec to its options

4.59.1 Detailed Description

The class that manages video codecs used by video loaders and savers.


The documentation for this class was generated from the following file:

- VideoCodecManager.h

4.60 CxVideoSource Class Reference

Provides image sequences.

Inheritance diagram for CxVideoSource:



Public Member Functions

- virtual QString **title** () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual int **buffersCountInMemoryPool** () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*
- virtual [CxChainable](#) * **clone** ()
Creates a deep copy of this object.
- virtual void **run** ()
This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).
- virtual bool **saveSettings** (QDomDocument &aDoc, QDomElement &aElm) const
The method saves the settings of this object to the DOM element.

- virtual bool [loadSettings](#) (const QDomElement &aDom)
Loads settings from DOM.
- virtual bool [currentOutputImageInfo](#) (SxPicBufInfo &picInfo, CxImageMetadata *pMetadata=NULL)
Gets the current output image info (format) provided by this object.
- virtual int **framesCount** ()
- virtual void **requestFrame** (int iFrameNo)
- virtual double [frameRate](#) ()
Fps, it defines the duration of the video.
- virtual void **setDataLoader** (IxImageDataLoader *pLoader)

Protected Attributes

- [IxImageDataLoader](#) * **m_pLoader**
Our source of the images.
- [SxPicBufInfo](#) **m_almgFormat**
- qint32 **m_iFrameCount**
- double **m_dFps**

Additional Inherited Members

4.60.1 Detailed Description

Provides image sequences.

4.60.2 Member Function Documentation

4.60.2.1 virtual bool CxVideoSource::currentOutputImageInfo (SxPicBufInfo & *picInfo*, CxImageMetadata * *pMetadata* = NULL) [virtual]

Gets the current output image info (format) provided by this object.

Remarks

The method should call [getCurrentOutputImageInfo\(\)](#) on its predecessors and calculate the output image info provided by this object.

Reimplemented from [CxImageProvider](#).

4.60.2.2 virtual void CxVideoSource::run () [virtual]

This method is invoked automatically by the [start\(\)](#) method. If the method is implemented as an infinite loop then it must use the [CHAINABLE_RUN_WHILE_TRUE](#) and [CHAINABLE_RUN_END_WHILE](#) macros to start and finish the infinite loop. For most of the chainables object it may not be implemented. It should be implemented for image sources ([CxImageSource](#)).

See also

[CxChainable](#)

Reimplemented from [CxChainable](#).

4.60.2.3 `virtual bool CxVideoSource::saveSettings ( QDomDocument & aDoc, QDomElement & aElm ) const [virtual]`

The method saves the settings of this object to the DOM element.

Remarks

The element `aElm` is already named by a caller, do not change the name (ie. do not call the `aElm.setTag←Name()`). You can add attributes and child elements.

Reimplemented from [CxChainable](#).


The documentation for this class was generated from the following file:

- VideoSource.h

4.61 CxVideoToFile Class Reference

The chainable class that allows to save images coming from camera to a video file.

Inheritance diagram for CxVideoToFile:



Signals

- void **frameSaved** (int iSeqIndex, quint64 uiTimeStamp)

Public Member Functions

- **CxVideoToFile** ([IxImageDataSaver](#) *pVideoSaver)
- virtual QString **title** () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual bool **acceptsDataFrom** ([CxChainable](#) *pPredecessor)
Query if this object can accept data from the given chainable object.
- virtual int **buffersCountInMemoryPool** () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*
- virtual [CxChainable](#) * **clone** ()
Creates a deep copy of this object.

Protected Slots

- virtual void **onDataAvailable** ([CxChainable](#) *pSender, int iFrameNo)

Protected Member Functions

- virtual bool [init](#) (QString *pErrMsg=NULL)
Initializes the object. When calling this method, the object MUST be stopped.
- virtual bool [initialized](#) ()

Private Attributes

- [IxImageDataSaver](#) * [m_pVideoSaver](#)
- int [m_iSeqIndex](#)

Additional Inherited Members

4.61.1 Detailed Description

The chainable class that allows to save images coming from camera to a video file.

4.61.2 Member Function Documentation

4.61.2.1 virtual bool CxVideoToFile::init (QString * *pErrMsg* = NULL) [protected],[virtual]

Initializes the object. When calling this method, the object MUST be stopped.

Remarks

Methods [init\(\)](#) and [initialized\(\)](#) are internally used by the [start\(\)](#) method. It should initialize the chainable object so that it is ready to start ([start\(\)](#)) . If this object is the [CxImageProvider](#) then the method initializes internal memory pool.

Implements [CxChainable](#).


The documentation for this class was generated from the following file:

- VideoToFile.h

4.62 CxViewImageData Class Reference

Implementation of [CxImageData](#) intended for displaying in views.

Inheritance diagram for CxViewImageData:



Public Member Functions

- [CxViewImageData](#) ([CxImageData](#) *pOrigData, [SxPicBuf](#) *pViewPicBufRgb32)
Constructor.
- virtual [IxChainData](#) * [clone](#) (TxMemPoolHandle hMemPool, QObject *pMemPoolClient) const
Creates the deep copy of this object. If the hMemPool is not NULL and it is a valid memory pool, then the data will be allocated from the memory pool. If the hMemPool is NULL then the data will be allocated from heap (operator new).
- virtual bool [isFromMemoryPool](#) (TxMemPoolHandle hMemPool) const
Returns true if the original picbuf ([origImageData\(\)](#)->[picBuf\(\)](#)) or the current picbuf is allocated from the given memory pool.
- const [CxImageData](#) * [origImageData](#) () const
Returns the original image data.

Private Attributes

- [CxImageData](#) * [m_pOrigData](#)

Additional Inherited Members

4.62.1 Detailed Description

Implementation of [CxImageData](#) intended for displaying in views.

Remarks

The class keeps pointer to the original data and to the data converted to RGB32 which can be displayed by QGraphicsView.

4.62.2 Constructor & Destructor Documentation

4.62.2.1 CxViewImageData::CxViewImageData (CxImageData * pOrigData, SxPicBuf * pViewPicBufRgb32)

Constructor.

Remarks

The object takes ownership of the pOrigData and pViewPicBuf. If the pViewPicBuf is NULL then the object internally uses pOrigData->[picBuf\(\)](#) as a view picbuf.

4.62.3 Member Data Documentation

4.62.3.1 CxImageData* CxViewImageData::m_pOrigData [private]

Original image data.


The documentation for this class was generated from the following file:

- ImageData.h

4.63 CxViewImageReceiver Class Reference

This class is primarily used by images views to receive a image data from a chain.

Inheritance diagram for CxViewImageReceiver:



Public Member Functions

- void [setRenderingEnabled](#) (bool bEnable)
Enable rendering in view. Can be disabled when wanting to get highest possible FPS on camera.
- bool [isRenderingEnabled](#) ()
Returns true when rendering is enabled.
- void [setForceFpsLimit](#) (bool bForce, int iFpsLimit=-1)
User-defined fps limit. Limit should be above zero.
- bool [isFpsLimitForced](#) (int *piFpsLimit=NULL)
- virtual QString [title](#) () const
Gets a title of this chainable object. The title is usually used to display a window title in the GUI etc. It should be a localized string.
- virtual [CxChainable](#) * [clone](#) ()
Creates a deep copy of this object.
- virtual int [buffersCountInMemoryPool](#) () const
*Gets a number of buffers which **this** object needs to have reserved in memory pool. If this object allocates no buffers then the method should return 0. See the [CxChainable](#).*
- virtual int [queryImagesCountInMemoryPool](#) (const [SxPicBuInfo](#) &picInfo, TxMemPoolHandle hPool)
Gets the required number of images reserved in the memory pool by its successors.
- virtual QWidget * [configurationWidget](#) ()
Displays the object's configuration dialog.
- virtual bool [hasConfigurationWidget](#) () const
Returns true when the the object has configuration widget (ie. the [configurationWidget\(\)](#) returns not NULL).
- virtual void [setMaxFps](#) (double dFps)
Sets max. FPS (when zero or negative, then the FPS will not be regulated)

Protected Member Functions

- virtual [IxChainData](#) * [processData](#) ([IxChainData](#) *pReceivedData)
The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Protected Attributes

- bool [m_bRenderingEnabled](#)
Enable rendering in view.
- bool [m_bForceFpsLimit](#)
User-set fps limit (when false, our automatics code is used)
- int [m_iFpsLimitForced](#)
User-set fps limit.
- double [m_dFpsLimitAutomatic](#)
Fps limit as defined by CamTool (from monitor frequency)

Additional Inherited Members

4.63.1 Detailed Description

This class is primarily used by images views to receive a image data from a chain.

Remarks

It converts the incoming data to RGB32 and creates an instance of [CxViewImageData](#) for the view.

Warning

It is not recommended to use this class in plugins or chainables because the usage of this classe somewhere in the middle of the chain can have a lot of unwanted side effects! Just forget this class exists and always work with [CxImageData](#).

4.63.2 Member Function Documentation

4.63.2.1 `virtual QWidget* CxViewImageReceiver::configurationWidget ( ) [virtual]`

Displays the object's configuration dialog.

Remarks

The configuration widget should not contain the Ok and Cancel button as they will be added automatically by the application.

Warning

When this chainable object changes the image format, don't forget to stop the current chain before applying new value from the dialog!
This method MUST be called from GUI (main) thread only!!

Returns

The default value is NULL. When overridden then it should return the pointer to the configuration widget.

Reimplemented from [CxChainable](#).

4.63.2.2 `virtual IxChainData* CxViewImageReceiver::processData ( IxChainData * pReceivedData ) [protected], [virtual]`

The method processes the received data. It is called automatically in the default implementation of the [onDataAvailable\(\)](#) slot.

Remarks

If the implementation creates (allocates) a new data here then the original (received) data should be deleted using the delete operator. The default implementation of this method returns a pointer to the received data without any processing.

Warning

If the method returns NULL, then the `pReceivedData` will be deleted automatically!

Reimplemented from [CxRgb32Convertor](#).

4.63.2.3 `virtual int CxViewImageReceiver::queryImagesCountInMemoryPool ( const SxPicBufInfo & picInfo, TxMemPoolHandle hPool ) [virtual]`

Gets the required number of images reserved in the memory pool by its successors.

Remarks

If this object uses its own memory pool then the method allways returns zero!


Reimplemented from [CxImageProvider](#).

The documentation for this class was generated from the following file:

- ViewImageReceiver.h

4.64 CxViewReceiverConfWidget Class Reference

Inheritance diagram for CxViewReceiverConfWidget:



Public Member Functions

- **CxViewReceiverConfWidget** ([CxViewImageReceiver](#) *pReceiver, QWidget *parent=0)

Private Slots

- void **on_chkRendering_clicked** (bool bChecked)
- void **on_rbFpsAuto_clicked** (bool bChecked)
- void **on_rbFpsForced_clicked** (bool bChecked)
- void **on_spinFps_edited** (int iValue)

Private Attributes

- [CxViewImageReceiver](#) * **m_pReceiver**

The documentation for this class was generated from the following file:

- ViewImageReceiver.h

4.65 CxHsiCameraParams::CxVirtualBand Class Reference

< Band record inside <virtual_bands>

Public Member Functions

- bool **isValid** () const
- QString **name** () const
Returns virtual band name as "123 nm".

Public Attributes

- double **m_dWaveLength**
Central wavelength of the virtual band in nm.
- double **m_dFwhm**
Full width of the virtual band at half the band's maximum, in nm.
- QVector< double > **m_vecCorrectionData**
Correction data, size is equal to vecBands.count.

4.65.1 Detailed Description

< Band record inside <virtual_bands>

The documentation for this class was generated from the following file:

- HsiCameraParams.h

4.66 IxAppDelegate Class Reference

Class for plugins to communicate with the application. Use [IxAppDelegate::instance\(\)](#) to get the pointer.

Public Types

- enum [ExSound](#) { [eSoundCapture](#), [eSoundRecord](#) }

Public Member Functions

- virtual QMainWindow * **mainWindow** ()=0
For adding something to main menu, or anything.
- virtual QToolBar * **mainToolBar** ()=0
For adding toolbar buttons.

- virtual QMenu * [mainMenu](#) (ExMainMenu eMenu)=0
For adding actions to main menu.
- virtual QSettings * [settings](#) ()=0
Use for storing plugin settings.
- virtual QString [appDataDirectory](#) ()=0
Store there your internal data, preserved for future sessions.
- virtual QString [lastOpenedDir](#) ()=0
Directory use last used.
- virtual void [setLastOpenedDir](#) (const QString &sDir)=0
Set directory use last used.
- virtual void [addDockedWidget](#) (QWidget *pWidget, Qt::DockWidgetArea eDockArea, Qt::DockWidgetAreas aAllowedAreas, bool bShow)=0
Add new widget to main window. QDockWidget is created around it and docked in correct place.
- virtual void [showNotification](#) (const QString &sMessage)=0
Popup notification in top right corner.
- virtual CxChain * [activeChain](#) ()=0
Chain of the active view.
- virtual int [activeViewId](#) ()=0
iViewId of the current view
- virtual void [activateView](#) (int iViewId)=0
Let the view become selected.
- virtual void [enableGUI](#) (bool bEnable)=0
Disables or enables application GUI.
- virtual CxChain * [loadChainFromXML](#) (const QString &sXml, bool bStart=true, CxImageSource *pImage←Source=NULL, QString *pErrMsg=NULL)=0
Loads chain from XML, attaches all views and starts it. If pImageSource is not NULL then the chain will be appended to it and the image source from the XML will be omitted.
- virtual void [showChainableConfiguration](#) (CxChainable *pChainable, QWidget *pParent)=0
Displays a modal configuration dialog for given widget if it exists.
- virtual CxMdiView * [viewPtr](#) (int iViewId)=0
Pointer to view by its id.
- virtual CxChain * [viewChain](#) (int iViewId)
Pointer to chain the view is belonging to.
- virtual bool [captureImageInView](#) (int iViewId=-1)=0
Duplicates the image in active view (as in menu File - Capture)
- virtual bool [runAcquisition](#) (bool bStart, XICORE_HANDLE hCamera=NULL)=0
Start or stop acquisition. Camera can be NULL in case of one camera connected.
- virtual XICORE_HANDLE [cameraHandleFromView](#) (int iViewId)=0
Camera used in the given view (returns NULL when no camera involved)
- virtual QList< int > [viewIdsForCameraHandle](#) (XICORE_HANDLE hCamera)=0
All views associated with given camera handle.
- virtual SxPicBuf * [capturePictureWithGUI](#) (QWidget *pParent, XICORE_HANDLE hCamera, const QString &sWndTitle, const QString &sInstructions, bool bEnableAveraging, int &iAverageCount, CxImageMetadata *pDstMetadata=NULL)=0
- virtual bool [openFile](#) (const QString &sFilename, bool bAddToMRU=true)=0
Open file as new view. Optionally add the file name to recently used list.
- virtual bool [createViewWithImage](#) (CxImageData *pImgData, const QString &sTitle)=0
Creates new view using the image. View takes ownership of the data and will free it when needed.
- virtual const CxImageData * [viewCurrentImage](#) (int iViewId)=0
Image currently displayed in view.
- virtual TxVecObjID [addVectorObject](#) (int iViewId, CxVectorObject *pObj)=0

- Adds vector object to view. You should store the returned ID instead of object pointer.*
- virtual [CxVectorObject](#) * [vectorObjectPtr](#) (int iViewId, TxVecObjID idObj)=0
In case you want to manipulate the object directly.
- virtual void [removeVectorObject](#) (int iViewId, TxVecObjID idObj)=0
Deletes the object.
- virtual QVector< TxVecObjID > [listAllVectorObjects](#) (int iViewId)=0
Returns the list of all objects in the scene.
- virtual void [overrideViewCursor](#) (int iViewId, QCursor *pCursor)=0
Pass NULL to stop overriding.
- virtual bool [mapPositionToImage](#) (int iViewId, const QPoint &ptScreenPos, QPointF *pPtInImage)=0
Recalculate to position inside the image. Returns true when inside.
- virtual TxBinaryLayerID [addBinaryLayer](#) (int iViewId, [CxBinaryLayer](#) *pLayer, bool bAutoColor=true)=0
Add new binary layer.
- virtual [CxBinaryLayer](#) * [binaryLayerPtr](#) (int iViewId, TxBinaryLayerID idLayer)=0
Returns a pointer to the binary layer with given id.
- virtual void [removeBinaryLayer](#) (int iViewId, TxBinaryLayerID idLayer)=0
Deletes the binary layers.
- virtual QList< TxBinaryLayerID > [listAllBinaryLayers](#) (int iViewId)=0
Returns a list of all binary layers in the view.
- virtual QList< TxBinaryLayerID > [listVisibleBinaryLayers](#) (int iViewId)=0
Returns a list of visible binary layers in the view.
- virtual void [setCustomPlotColors](#) (QCustomPlot *pPlot)=0
Sets app-style colors to graph grawing component.
- virtual void [playSound](#) ([ExSound](#) eSound)=0
Plays the application predefined sound.
- virtual void [playSound](#) (const QString &sFilename)=0
Plays a sound in WAV format.

Static Public Member Functions

- static [IxAppDelegate](#) * [instance](#) ()
Returns the pointer to actual application delegate.
- static QByteArray [fileSystemRepresentation](#) (const QString &sFilename)
- static QString [stringFromFileSystemRepresentation](#) (const QByteArray &sEncodedFilename)
- static void [setInstance](#) ([IxAppDelegate](#) *pInstance)
Reserved for internal use.

4.66.1 Detailed Description

Class for plugins to communicate with the application. Use [IxAppDelegate::instance\(\)](#) to get the pointer.

Warning

This class should always be used in the main (GUI) thread. You should never use it in [CxChainable::processData\(\)](#). In this case the application behaviour is undefined.

4.66.2 Member Enumeration Documentation

4.66.2.1 enum IxAppDelegate::ExSound

Enumerator

- eSoundCapture** When image is captured or autosaved (instant save).
- eSoundRecord** When instant video recordint is started or stopped.

4.66.3 Member Function Documentation

4.66.3.1 `virtual SxPicBuf* IxAppDelegate::capturePictureWithGUI ( QWidget * pParent, XICORE_HANDLE hCamera, const QString & sWndTitle, const QString & sInstructions, bool bEnableAveraging, int & iAverageCount, CxImageMetadata * pDstMetadata = NULL ) [pure virtual]`

Show dialog that lets user to capture a image from camera, returning the picture. Caller must free and delete the picture). If you leave pParent NULL, mainWnd is used as parent.


The documentation for this class was generated from the following file:

- AppDelegate.h

4.67 IxChainData Class Reference

The base class for all data exchanged by [CxChainable](#) objects.

Inheritance diagram for IxChainData:



Public Member Functions

- `virtual quint32 frameNo () const =0`
Gets the frame number of the data. For "static" data (e.g. static image) it is always zero.
- `virtual IxChainData * clone (TxMemPoolHandle hMemPool, QObject *pMemPoolClient) const =0`
Creates the deep copy of this object. If the hMemPool is not NULL and it is a valid memory pool, then the data will be allocated from the memory pool. If the hMemPool is NULL then the data will be allocated from heap (operator new).
- `virtual IxChainData * clone (QObject *pMemPoolClient) const =0`
Creates the deep copy of this object from the same memory source as the original data.
- `virtual bool isFromMemoryPool (TxMemPoolHandle hMemPool) const =0`
Checks if the data is from given memory pool.

4.67.1 Detailed Description

The base class for all data exchanged by [CxChainable](#) objects.


The documentation for this class was generated from the following file:

- Data.h

4.68 IxChainDataLoader Class Reference

The base class for any kind of chain data loading.

Inheritance diagram for IxChainDataLoader:



Public Member Functions

- virtual bool **load** ()=0
Loads data from some source.
- virtual quint32 **sequenceSize** ()=0
In case the source has multiple data items, their count is returns here.
- virtual **IxChainData** * **loadedData** (quint32 iSeqIdx=0)=0
*Gets the pointer to the current data (after the **load()** method was called).*
- virtual QString **inputFile** ()=0
If the loader load data from file then this method returns the the file's name.
- virtual void **setInputFile** (const QString &sFileName)=0

4.68.1 Detailed Description

The base class for any kind of chain data loading.

4.68.2 Member Function Documentation

4.68.2.1 virtual IxChainData* IxChainDataLoader::loadedData (quint32 iSeqIdx = 0) [pure virtual]

Gets the pointer to the current data (after the **load()** method was called).

Remarks

The caller of this method is the owner of the returned data so that it has to take care about its deletion using the **delete** operator.

Implemented in **CxCircularBuffer**.


The documentation for this class was generated from the following file:

- DataLoader.h

4.69 IxChainDataSaver Class Reference

The base class for any kind of class for saving **IxChainData**.

Inheritance diagram for IxChainDataSaver:



Public Member Functions

- virtual void [setSequenceSize](#) (qint32 iSeqSize)
In case the source has multiple data items, we let saver know the number of items.
- virtual bool [save](#) (const [IxChainData](#) *pData, qint32 iSeqIdx=0)=0
Saves the current data. May be part of sequence.

4.69.1 Detailed Description

The base class for any kind of class for saving [IxChainData](#).


The documentation for this class was generated from the following file:

- [DataSaver.h](#)

4.70 IxImageDataLoader Class Reference

The base class for image data loading.

Inheritance diagram for IxImageDataLoader:



Public Member Functions

- virtual double [loadedFps](#) ()=0
In case the source is a sequence, this should get its fps. Return 0 if uncertain.
- virtual [TxImageFormatList](#) [supportedImageFormats](#) ()=0
Image formats capable of loading.

4.70.1 Detailed Description

The base class for image data loading.

Remarks

The implementations of this class should be registered by CxRtti so that the application could use them.


The documentation for this class was generated from the following file:

- [DataLoader.h](#)

4.71 IxImageDataSaver Class Reference

The base class for any kind of image data saver.

Inheritance diagram for IxImageDataSaver:



Public Member Functions

- virtual TxImageFormatList [supportedImageFormats](#) ()=0
Image formats capable of saving.
- virtual void **setOutputFile** (const QString &sFileName)=0
- virtual bool **querySavedImageFormat** (const [SxPicBufInfo](#) &aInput, [SxPicBufInfo](#) &aOutput)
- virtual void **setSaverOptions** (const TxImageFormatOptions &aOpt)
- virtual QWidget * **widgetForSaveDlg** (QWidget *pParent)
- virtual TxImageFormatOptions **optionsFromWidgetInSaveDlg** (QWidget *pWidget)
- virtual TxImageFormatOptions **codecOptionsFromWidgetInSaveDlg** (QWidget *pWidget, TxFourCC &fcc↔ Codec)

4.71.1 Detailed Description

The base class for any kind of image data saver.

Remarks

The implementations of this class should be registered by CxRtti so that the application could use them.

The documentation for this class was generated from the following file:

- DataSaver.h

4.72 IxVideoAccessor Class Reference

The class for accessing the raw video data buffers.

Public Member Functions

- virtual qint32 [framesCount](#) ()=0
Number of frames of the video sequence.
- virtual size_t [frameDataSize](#) (qint32 iSeqIndex)=0
Size of the frame with the given index in bytes.
- virtual bool [frameData](#) (qint32 iSeqIndex, void *pData, size_t &ruiSize)=0
Copies the raw frame data with given index to the preallocated pData.
- virtual qint32 [lastProcessedFrame](#) ()=0
Returns the index of last processed (decoded) frame.
- virtual void [keyFrames](#) (QList< qint32 > &rKeyFrames)=0
Returns a list of keyframes of the video sequence.

4.72.1 Detailed Description

The class for accessing the raw video data buffers.

Remarks

This class should implemented by video loaders. It is passed to [IxVideoDecoder](#) instances to access the video data.

The documentation for this class was generated from the following file:

- VideoAccessor.h

4.73 IxVideoDecoder Class Reference

The base class of all video decoders.

Public Member Functions

- virtual QString [name](#) ()=0
Localized name of the decoder visible in GUI.
- virtual TxFourCC [fourCC](#) ()=0
- virtual bool [decode](#) (qint32 iFrameIndex, [IxVideoAccessor](#) *pVideo, [SxPicBuf](#) &rDecodedPic)=0
The FourCC identifier of the codec.

4.73.1 Detailed Description

The base class of all video decoders.

Remarks

The implementations of this class should be registered by CxRtti so that the application could use them.

4.73.2 Member Function Documentation

4.73.2.1 virtual bool IxVideoDecoder::decode (qint32 iFrameIndex, IxVideoAccessor * pVideo, SxPicBuf & rDecodedPic) [pure virtual]

The FourCC identifier of the codec.

The method decodes a video frame given by iFrameIndex.

Remarks

If the rDecodedPicBuf.m_pData==NULL, the decoder will allocate it but the image must be initialized correctly (width, height, channels, datatype...)

The documentation for this class was generated from the following file:

- VideoDecoder.h

4.74 IxVideoEncoder Class Reference

The base class of all video encoders.

Public Member Functions

- virtual QString **name** ()=0
Localized name name visible in GUI.
- virtual TxFourCC **fourCC** ()=0
The FourCC identifier of the codec.
- virtual bool **encode** (const SxPicBuf &rPic, void *pEncoded, size_t &ruiSize, bool *pbKeyFrame=NULL)=0
Encodes next image (rPic) in the stream, output is stored to pEncoded buffer with size ruiSize.
- virtual size_t **maxEncodedSize** (const SxPicBuflInfo &rPicInfo)=0
- virtual void **queryOuputFormat** (const SxPicBuflInfo &rInput, SxPicBuflInfo &rOutput)=0
For the given input format it sets the output format.
- virtual void **setEncoderOptions** (const TxImageFormatOptions &aOpt)
- virtual QWidget * **widgetForSaveDlg** (QWidget *pParent)
- virtual TxImageFormatOptions **optionsFromWidgetInSaveDlg** (QWidget *pWidget)

4.74.1 Detailed Description

The base class of all video encoders.

Remarks

The implementations of this class should be registered by CxRtti so that the application could use them.

4.74.2 Member Function Documentation

- 4.74.2.1 virtual bool lxVideoEncoder::encode (const SxPicBuf & rPic, void * pEncoded, size_t & ruiSize, bool * pbKeyFrame = NULL) [pure virtual]

Encodes next image (rPic) in the stream, output is stored to pEncoded buffer with size ruiSize.

Parameters

<i>pbKeyFrame</i>	is optional. When not NULL and set to true on input, then the encoder should force this image to be a key frame. On the output it gives the information whether image was encoded as the key frame.
-------------------	---

The documentation for this class was generated from the following file:

- VideoEncoder.h

4.75 CxHsiCameraParams::SxBandIdx Struct Reference

Public Attributes

- int **zone**
- int **band**

The documentation for this struct was generated from the following file:

- HsiCameraParams.h

4.76 SxCameraInfo Struct Reference

Information about camera device.

Public Attributes

- [QString m_sName](#)
Camera model name (XI_PRM_DEVICE_NAME)
- [QString m_sType](#)
Device type (XI_PRM_DEVICE_TYPE)
- [int m_iDeviceModelId](#)
Camera model ID (XI_PRM_DEVICE_MODEL_ID)
- [QString m_sSerial](#)
Serial number. We expect it to be unique. (XI_PRM_DEVICE_SN)
- [QString m_sSensorSerial](#)
Serial number of the sensor. We expect it to be unique. (XI_PRM_DEVICE_SENS_SN)
- [QString m_sUserId](#)
Camera user id (as set by user, and loaded later opening the camera)
- [QString m_sInstancePath](#)
Instance path of the device. (XI_PRM_DEVICE_INSTANCE_PATH)
- [ExColorFilterArray m_eColorFilterArray](#)
Color filter array type of RAW data (XI_PRM_COLOR_FILTER_ARRAY)

4.76.1 Detailed Description

Information about camera device.

The documentation for this struct was generated from the following file:

- CameraParameters.h

4.77 CxVideoCodecManager::SxCodecItem Struct Reference

Public Attributes

- [QString m_sRttiClass](#)
- [QString m_sName](#)
Localized name of the codec.
- [TxFourCC m_fcc](#)

The documentation for this struct was generated from the following file:

- VideoCodecManager.h

4.78 CxCameraParam::SxEnumItem Struct Reference

Struct holding all info for each item of enumerated parameter (ecpEnum)

Public Member Functions

- [SxEnumItem](#) (const [CxValue](#) &vVal, const [QString](#) &sCaption)

Public Attributes

- [CxValue](#) **m_vValue**
- [QString](#) **m_sCaption**

4.78.1 Detailed Description

Struct holding all info for each item of enumerated parameter (ecpEnum)

The documentation for this struct was generated from the following file:

- CameraParameters.h

4.79 CxMemoryManager::SxLimitDef Struct Reference

Public Member Functions

- **bool operator<** (const [SxLimitDef](#) &other) const

Public Attributes

- [QString](#) **m_sOS**
- [QString](#) **m_sBits**
- [QMap](#)< [QString](#), [QString](#) > **m_mapParams**
- [quint64](#) **m_uiLimitBytes**


The documentation for this struct was generated from the following file:

- MemoryManager.h

4.80 SxPicBuf Struct Reference

Picture buffer.

Inheritance diagram for SxPicBuf:



Public Types

- enum [ExDataOwner](#) { [edoNone](#), [edoThisPicBuf](#), [edoMemoryPool](#) }

Public Member Functions

- [SxPicBuf](#) ()
Default constructor, initialies all to zero.
- [size_t](#) [dataSize](#) () const

- Returns the size of pixel data ([m_pData](#)).
- void [setPicFormat](#) (const [SxPicBufInfo](#) &rInfo)
Sets the width, height, format according to the rInfo parameter.
- bool [isFromMemoryPool](#) (TxMemPoolHandle hMemPool) const
Data comes from some memory pool.
- void [setAsShallowCopyOf](#) (const [SxPicBuf](#) &rOther)
Use with caution while hacking. Should not be needed in 99% of cases.
- void [copyPicAdditionalData](#) (const [SxPicBuf](#) &rOther)
Copy frame number, time stamps and other info values.

Public Attributes

- void * [m_pData](#)
- size_t [m_uiAllocSize](#)
- enum [SxPicBuf::ExDataOwner](#) [m_eDataOwner](#)
- TxMemPoolHandle [m_hMemoryPoolOwner](#)
- QObject * [m_pMemoryPoolClient](#)
- quint32 [m_iFrameNo](#)
- quint64 [m_uiTimeStamp](#)
- quint32 [m_iXiApiFrameNo](#)

Protected Member Functions

- [SxPicBuf](#) (const [SxPicBuf](#) &rOther)
Copy constructor is forbidden to use. See [operator=](#).
- [SxPicBuf](#) & [operator=](#) (const [SxPicBuf](#) &rOther)
[operator=](#) is forbidden to use. It is dangerous to copy pointer allocated data.

4.80.1 Detailed Description

Picture buffer.

4.80.2 Member Enumeration Documentation

4.80.2.1 enum [SxPicBuf::ExDataOwner](#)

Enumerator

- [edoNone](#)** You should delete the [m_pData](#) yourself
- [edoThisPicBuf](#)** Data belong to this picbuf.
- [edoMemoryPool](#)** Data were allocated inside memory pool, [m_hMemoryPoolOwner](#)

4.80.3 Member Data Documentation

4.80.3.1 TxMemPoolHandle [SxPicBuf::m_hMemoryPoolOwner](#)

Data is part of this memory pool. Do not free.

4.80.3.2 quint32 [SxPicBuf::m_iFrameNo](#)

Frame number is set when image comes from camera or an image sequence. Else zero.

4.80.3.3 qint32 SxPicBuf::m_iXApiFrameNo

Frame number as set by xiApi. May be different to our own counter m_iFrameNo. -1 for undefined.

4.80.3.4 void* SxPicBuf::m_pData

Pointer to pixel data.

4.80.3.5 QObject* SxPicBuf::m_pMemoryPoolClient

Who allocated the buffer in the pool (may have requests).

4.80.3.6 size_t SxPicBuf::m_uiAllocSize

The allocated size of data (i.e. available size)

4.80.3.7 quint64 SxPicBuf::m_uiTimeStamp

The timestamp in micro seconds.


The documentation for this struct was generated from the following file:

- PicBuf.h

4.81 SxPicBufInfo Class Reference

Information about picture buffer.

Inheritance diagram for SxPicBufInfo:



Public Member Functions

- [SxPicBufInfo \(\)](#)
Default constructor.
- quint32 [bytesPerPixel \(\)](#) const
Returns the number of bytes per pixel.

Public Attributes

- quint32 [m_uiWidth](#)
- quint32 [m_uiHeight](#)
- quint32 [m_uiStride](#)
- quint32 [m_uiChannels](#)
- quint32 [m_uiBpc](#)
- ExImageDataType [m_eDataType](#)

- ExDataStorageFormat [m_eFormat](#)
- ExColorFilterArray [m_eColorFilterArray](#)

4.81.1 Detailed Description

Information about picture buffer.

4.81.2 Member Data Documentation

4.81.2.1 ExColorFilterArray SxPicBufInfo::m_eColorFilterArray

Bayer pattern

4.81.2.2 ExImageDataType SxPicBufInfo::m_eDataType

Data type of the image.

4.81.2.3 ExDataStorageFormat SxPicBufInfo::m_eFormat

Format of the image.

4.81.2.4 quint32 SxPicBufInfo::m_uiBpc

Bits per channel.

4.81.2.5 quint32 SxPicBufInfo::m_uiChannels

Number of color channels.

4.81.2.6 quint32 SxPicBufInfo::m_uiHeight

Height in pixels.

4.81.2.7 quint32 SxPicBufInfo::m_uiStride

Line width in bytes.

4.81.2.8 quint32 SxPicBufInfo::m_uiWidth

Width in pixels.

The documentation for this class was generated from the following file:

- PicBuf.h

Chapter 5

File Documentation

5.1 Chainable.h File Reference

Classes

- class [CxChainable](#)
A base class of all chainable objects.
- class [CxImageProvider](#)
The base class for any chainable object which provides images.
- class [CxImageSource](#)
The base class of any image provider which is an image source, i.e. it has no predecessors and it provides images.

Macros

- `#define CHAINABLE_RUN_WHILE_TRUE`
The macro must be used to start a loop in the [CxChainable::run\(\)](#) in a case that the method is implemented as an infinite loop ([CxChainable::eRunsInfiniteLoop](#)).
- `#define CHAINABLE_RUN_END_WHILE`
The macro must be used to finish a loop in the [CxChainable::run\(\)](#) in a case that the method is implemented as an infinite loop ([CxChainable::eRunsInfiniteLoop](#)).

5.1.1 Macro Definition Documentation

5.1.1.1 `#define CHAINABLE_RUN_END_WHILE`

Value:

```
} \
    if(__bStateChanged__)\
    {\
        if(__eState__ == eStopped)\
        {\
            LOG_TRACE(ERR_MSG(CxChainable::run\(\)), "Releasing the m_semStartStop (stopped).")\
        ;\
            m_semStartStop.release();\
        }\
    }
```

The macro must be used to finish a loop in the [CxChainable::run\(\)](#) in a case that the method is implemented as an infinite loop ([CxChainable::eRunsInfiniteLoop](#)).

See also

[CHAINABLE_RUN_WHILE_TRUE](#)

5.1.1.2 #define CHAINABLE_RUN_WHILE_TRUE

Value:

```
Q_ASSERT(m_eFlags.testFlag(eRunIsInfiniteLoop));\
LOG_TRACE(ERR_MSG(CxChainable::run(), "Releasing the m_semStartStop (started).")); \
m_semStartStop.release();\
bool __bStateChanged__ = false;\
CxChainable::ExChainableState __eState__ = eUndefined;\
while(true)\
{\
    {\
        QMutexLocker lock(&m_lock);\
        if(m_eReqState == eStopped)\
        {\
            m_eState = eStopped;\
            __eState__ = m_eState;\
            __bStateChanged__ = true;\
            m_eReqState = eUndefined;\
            break;\
        }\
    }\
    if(chain()->stopped())\
    {\
        continue;\
    }\
}
```

The macro must be used to start a loop in the [CxChainable::run\(\)](#) in a case that the method is implemented as an infinite loop ([CxChainable::eRunsInfiniteLoop](#)).

See also

[CHAINABLE_RUN_END_WHILE](#)